



# Safety Data Sheet

## Clinical No-Rinse Spray Cleanser

### Section 1. Identification

Product Identifier	Clinical No-Rinse Spray Cleanser		
Synonyms	MSC092204; MSD_SDS0072; MSC092208		
Manufacturer Stock Numbers	MSC092208; MSC092204		
Recommended use	N/A		
Uses advised against	N/A		
Manufacturer Contact			
Address	Medline Industries, Inc. 3 Lakes Drive Northfield, IL, 60093 USA		
	Phone	Emergency Phone	Fax
	(800) 633-5463	(800) 424-9300 CHEMTREC	(847) 643-4436
	Website		
	www.Medline.com		

### Section 2. Hazards Identification

Classification	No OSHA Hazard Classifications Applicable - Category N.A.
Signal Word	
Pictogram	
Hazard Statements	No OSHA Hazard Classifications Applicable
Precautionary Statements	
Response	N/A
Prevention	N/A
Storage	N/A

Disposal	N/A
Ingredients of unknown toxicity	0%
Hazards not Otherwise Classified	No Data Available

### Section 3. Ingredients

CAS	Ingredient Name	Weight %
-----	-----------------	----------

Occupational exposure limits, if available, are listed in Section 8.

### Section 4. First-Aid Measures

General Advice:	Take proper precautions to ensure your own health and safety before attempting rescue and providing first aid. Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.
Eye Contact:	If irritation develops flush eyes with large amounts of water for at least 15 minutes. Remove contact lenses, if worn. If irritation persists, seek medical attention.
Skin contact:	If irritation develops wash area with water. Get medical attention if irritation persists.
Inhalation:	Remove victim to fresh air and keep at rest in a position comfortable for breathing. Seek medical attention if discomfort continues or if you feel unwell.
Ingestion:	Never give anything by mouth to an unconscious person. Consult a physician if necessary.

### Section 5. Fire Fighting Measures

Suitable Extinguishing Media	In case of fire, use Water spray, Carbon Dioxide (CO <sub>2</sub> ), Foam, Dry Chemical
Unsuitable Extinguishing Media	N.D.
Special Fire Fighting Procedures:	N.D.
Unusual Fire and Explosion Hazards:	N.D.

### Section 6. Accidental Release Measures

**Non-hazardous Small Spill:** Should be cleaned up at the time of the spill. Take all necessary precautions and wear any personal protective equipment that is applicable. Dispose of per local and state regulations.

**Non-hazardous Large Spill:** Should be cleaned up at the time of the spill. May require special treatment, equipment and/or emergency assistance. Dispose of per local and state regulations.

## Section 7. Handling and Storage

**Handling Non-hazardous:** This product is considered to be an article which does not release or otherwise result in exposure to a hazardous chemical under normal use conditions.

**Storage Non-hazardous:** Store in original package. Discard unused portion. Do not reuse. No special storage conditions required.

## Section 8. Exposure Controls/Personal Protection

Occupational Exposure Limits	Ingredient Name	ACGIH TLV	OSHA PEL	STEL
Personal Protective Equipment	N/A			
Engineering Controls:	No specific measures are required provided the product is handled in accordance with the general rules of occupational hygiene and safety. Use in a well ventilated area to prevent exposure. Maintain eyewash fountain and quick-drench facilities in work areas.			
Respiratory Protection:	None required under normal use conditions.			
Skin Protection:	None required for normal use. For prolonged exposure, use appropriate goggles, protective clothing and gloves.			
Eye Protection:	None required for normal use. For prolonged exposure, use appropriate goggles, protective clothing and gloves.			
Other protective equipment:	Not expected to be necessary under normal conditions of use. Where exposure cannot be adequately controlled, use appropriate protective clothing or equipment.			
Work/hygienic practices:	Handle in accordance with good industrial hygiene and safety practices. Wash thoroughly with soap and water after handling and before eating, drinking, or using tobacco. Safety shower and eye wash should be available close to work areas.			

## Section 9. Physical and Chemical Properties

Physical State	Water-thin liquid
Color	Light yellow
Odor	Light vanilla
Odor Threshold	N.D.
Solubility	Soluble
Partition coefficient Water/n-octanol	N.D.
VOC%	N/A
Viscosity	N.D.

Specific Gravity	1.035
Density lbs/Gal	N/A
Pounds per Cubic Foot	N/A
Flash Point	N.D.
FP Method	N.D.
pH	N.D.
Melting Point	N.D.
Boiling Point	N.D.
Boiling Range	N.D.
LEL	N/A
UEL	N/A
Evaporation Rate	N.D.
Flammability	N.D.
Decomposition Temperature	N.D.
Auto-ignition Temperature	N.D.
Vapor Pressure	N.D.
Vapor Density	N.D.

## Section 10. Stability and Reactivity

Stability: Stable.

Incompatibility (Materials to avoid): N.D.

Hazardous Decomposition or Byproducts: None.

Hazardous Polymerization: N.D.

Conditions to avoid: None.

## Section 11. Toxicological Information

Acute Toxicity: Eye Contact Non-irritating based upon the available data.

Acute Toxicity: Skin Contact Non-irritating based upon the available data.

Acute Toxicity: Ingestion May be harmful if swallowed.

Carcinogenicity: None known.

Chronic Toxicity: None known.

## Section 12. Ecological Information

Toxicity: N.D.

Persistence and degradability: N.D.

Bioaccumulative potential: N.D.

Mobility in soil: N.D.

Other adverse effects: N.D.

## Section 13. Disposal

**Waste Disposal:** Waste must be disposed of in accordance with federal, state and local environmental control regulations.

## Section 14. Transport Information

UN Number N/A  
UN Proper Shipping Name Not Regulated  
DOT Classification Not Regulated  
Packing Group Not Regulated

## Section 15. Regulatory Information

SARA 311/312: Refer to Section 2 of the SDS.  
SARA 302: N.A.  
SARA 304: N.A.  
SARA 313: N.A.  
TSCA: All components are listed or exempt.  
CERCLA Hazardous Substance List: N.A.  
Clean Air Act (CAA) Section 112, 112 (r): N.A.  
  
State Regulations: N.A.

## Section 16. Other Information

Revision Date 8/23/2023  
  
Legend N.A. - Not Applicable  
N.E. - Not Established  
N.D. - Not Determined  
  
HMIS (U.S.A.): Health Hazard 0  
HMIS (U.S.A.): Flammability 0  
HMIS (U.S.A.): Reactivity 0  
HMIS (U.S.A.): Personal Protection X

National Fire Protection Association (U.S.A): Health Hazard

0

National Fire Protection Association (U.S.A): Fire Hazard

0

National Fire Protection Association (U.S.A): Instability Hazard

0

National Fire Protection Association (U.S.A): Special Hazard

0

Additional Information:

The information contained herein is furnished without warranty or legal responsibility of any kind. Employers should use this information only as a supplement to other information gathered by them and must make independent determination of suitability and completeness of information from all sources to assure proper use of these materials and the safety and health of employees.