

Surgical Sutures

McKesson surgical sutures are offered in a variety of materials, lengths and needle types. To provide options in wound and surgical closure, we offer both absorbable and non-absorbable materials. McKesson surgical sutures are sterile and individually-packaged in peel-open foil packs.

Features and Benefits

Non-Absorbable

- **Nylon:** High tensile strength with a uniform smooth surface to permit easy tissue passage; resists breakage and retains strength for long periods of time. Available as black monofilament.
- **Silk:** Excellent handling and tying characteristics. Modern braiding technique provides a uniform surface and high tensile strength. Available in black braided.
- **Polypropylene:** Versatile monofilament suture with exceptional tensile strength. Smooth and easy to handle with excellent knot security. Available as blue monofilament.
- **Polyester:** Silicone (impregnated) coating designed for precise handling and pliability. Uniform braiding process designed to provide greater resistance to breaking. Excellent knot security.



Absorbable

- **Chromic Gut:** Ensures predictable absorption and high tensile strength. Each strand has a uniform diameter for secure knotting.
- **Polyglycolic Acid:** Offers a smooth and low friction coating to reduce tissue drag. These absorbable sutures retain 50% of their tensile strength at two weeks and 20% of their strength at three weeks. Available as undyed monofilament.
- **PGA-PCL (glycolide & e-caprolactone):** Excellent absorption, pliability and handling. Popular for subcuticular closure and soft tissue approximation with 64%-76% retention of tensile strength at two week and 40% tensile strength retention at three weeks. Available as undyed monofilament.

Satisfaction Guaranteed
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them – making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

DESCRIPTION	NEEDLE TYPE	MFR #	QTY
NON-ABSORBABLE SURGICAL SUTURES			
MONOFILAMENT NYLON SUTURE			
6-0 Monofilament Nylon, 0.7 Metric, 18", Black, P-11, 3/8, 16 mm	Premium Reverse Cutting	S1665GX	12ea/bx
5-0 Monofilament Nylon, 1 Metric, 18", Black, P-12, 3/8, 19 mm	Premium Reverse Cutting	S1666G	12ea/bx
4-0 Monofilament Nylon, 1.5 Metric, 18", Black, P-12, 3/8, 19 mm	Premium Reverse Cutting	S1667GX	12ea/bx
5-0 Monofilament Nylon, 1 Metric, 18", Black, P-11, 3/8, 16 mm	Premium Reverse Cutting	S1668GX	12ea/bx
3-0 Monofilament Nylon, 2 Metric, 18", Black, P-12, 3/8, 19 mm	Premium Reverse Cutting	S1669H	36ea/bx
6-0 Monofilament Nylon, 0.7 Metric, 18", Black, P-13, 3/8, 13 mm	Premium Reverse Cutting	S1698GX	12ea/bx
5-0 Monofilament Nylon, 1 Metric, 18", Black, C-13, 3/8, 19 mm	Reverse Cutting	S661GX	12ea/bx
4-0 Monofilament Nylon, 1.5 Metric, 18", Black, C-13, 3/8, 19 mm	Reverse Cutting	S662GX	12ea/bx
3-0 Monofilament Nylon, 2 Metric, 18", Black, C-14, 3/8, 24 mm	Reverse Cutting	S663GX	12ea/bx
2-0 Monofilament Nylon, 3 Metric, 18", Black, C-15, 3/8, 26 mm	Reverse Cutting	S664GX	12ea/bx
6-0 Monofilament Nylon, 0.7 Metric, 18", Black, P-10, 3/8, 11 mm	Premium Reverse Cutting	S697GX	12ea/bx
5-0 Monofilament Nylon, 1 Metric, 18", Black, P-13, 3/8, 13 mm	Premium Reverse Cutting	S698GX	12ea/bx
4-0 Monofilament Nylon, 1.5 Metric, 18", Black, P-13, 3/8, 13 mm	Reverse Cutting	S699GX	12ea/bx
5-0 Monofilament Nylon, 1 Metric, 18", Black, C-1, 3/8, 12 mm	Reverse Cutting	S913BX	36ea/bx
4-0 Monofilament Nylon, 1.5 Metric, 18", Black, C-1, 3/8, 12 mm	Reverse Cutting	S931BX	36ea/bx
WAX COATED BRAIDED SILK SUTURE			
3-0 Wax Coated Braided Silk, 2 Metric, 30", Black, SC-2, 60 mm	Cutting	S622H	36ea/bx
4-0 Wax Coated Braided Silk, 1.5 Metric, 18", Black, C-13, 3/8, 19 mm	Reverse Cutting	S683GX	12ea/bx
MONOFILAMENT POLYPROPYLENE SUTURE			
5-0 Monofilament Polypropylene, 1 Metric, 36", Blue, CV-23, 1/2, 17 mm	Taper Point	S8556H	36ea/bx
3-0 Monofilament Polypropylene, 2 Metric, 18", Blue, C-13, 3/8, 19 mm	Reverse Cutting	S8665G	12ea/bx
4-0 Monofilament Polypropylene, 1.5 Metric, 18", Blue, P-12, 3/8, 19 mm	Premium Reverse Cutting	S8682GX	12ea/bx
4-0 Monofilament Polypropylene, 1.5 Metric, 18", Blue, C-13, 3/8, 19 mm	Reverse Cutting	S8683GX	12ea/bx
6-0 Monofilament Polypropylene, 0.7 Metric, 18", Blue, P-13, 3/8, 13 mm	Premium Reverse Cutting	S8695GX	12ea/bx
6-0 Monofilament Polypropylene, 0.7 Metric, 18", Blue, P-10, 3/8, 11 mm	Premium Reverse Cutting	S8697G	12ea/bx
5-0 Monofilament Polypropylene, 1 Metric, 18", Blue, P-13, 3/8, 13 mm	Premium Reverse Cutting	S8698GX	12ea/bx
4-0 Monofilament Polypropylene, 1.5 Metric, 18", Blue, P-13, 3/8, 13 mm	Premium Reverse Cutting	S8699GX	12ea/bx
COATED BRAIDED POLYESTER SUTURE			
5-0 Coated Braided Polyester, 1 Metric, 36", Blue, CV-331, 1/2, 16 mm	Taper Point	SX556H	36ea/bx

ORDERING INFORMATION

DESCRIPTION	NEEDLE TYPE	MFR #	QTY
ABSORBABLE SURGICAL SUTURES			
CHROMIC GUT ABSORBABLE SUTURE			
4-0 Chromic Gut, 2 Metric, 18", P-13, 3/8, 13 mm	Premium Reverse Cutting	S1654G	12ea/bx
4-0 Chromic Gut, 2 Metric, 30", C-13, 3/8, 19 mm	Reverse Cutting	S635GX	12ea/bx
3-0 Chromic Gut, 3 Metric, 30", C-13, 3/8, 19 mm	Reverse Cutting	S636GX	12ea/bx
5-0 Chromic Gut, 1.5 Metric, 18", P-13, 3/8, 13 mm	Premium Reverse Cutting	S687G	12ea/bx
3-0 Chromic Gut, 3 Metric, 30", V-20, 1/2, 26 mm	Taper Point	SG122H	36ea/bx
POLYGLYCOLIC ACID BRAIDED ABSORBABLE SUTURE			
2-0 Polyglycolic Acid Braided, 3 Metric, 30", Undyed, GS-21, 1/2, 37 mm	Taper Point	SJ259H	36ea/bx
2-0 Polyglycolic Acid Braided, 3 Metric, 30", Undyed, GS-22, 1/2, 27 mm	Taper Point	SJ269H	36ea/bx
3-0 Polyglycolic Acid Braided, 2 Metric, 30", Undyed, V-20, 1/2, 26 mm	Taper Point	SJ416H	36ea/bx
2-0 Polyglycolic Acid Braided, 3 Metric, 30", Undyed, V-20, 1/2, 26 mm	Taper Point	SJ417H	36ea/bx
5-0 Polyglycolic Acid Braided, 1 Metric, 30", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SJ421GX	36ea/bx
4-0 Polyglycolic Acid Braided, 1.5 Metric, 30", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SJ422H	36ea/bx
3-0 Polyglycolic Acid Braided, 2 Metric, 30", Undyed, C-13, 3/8, 19 mm	Reverse Cutting	SJ423H	36ea/bx
2-0 Polyglycolic Acid Braided, 3 Metric, 30", Undyed, C-14, 3/8, 24 mm	Reverse Cutting	SJ443H	36ea/bx
6-0 Polyglycolic Acid Braided, 0.7 Metric, 18", Undyed, P-13, 3/8, 13 mm	Premium Reverse Cutting	SJ492GX	12ea/bx
5-0 Polyglycolic Acid Braided, 1 Metric, 18", Undyed, P-13, 3/8, 13 mm	Premium Reverse Cutting	SJ493GX	12ea/bx
4-0 Polyglycolic Acid Braided, 1.5 Metric, 18", Undyed, P-13, 3/8, 13 mm	Premium Reverse Cutting	SJ494GX	12ea/bx
4-0 Polyglycolic Acid Braided, 1.5 Metric, 18", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SJ496GX	12ea/bx
3-0 Polyglycolic Acid Braided, 2 Metric, 18", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SJ497GX	12ea/bx
PGA-PCL MONOFILAMENT ABSORBABLE SUTURE			
4-0 PGA-PCL, 1.5 Metric, 30", Undyed, P-12g, 3/8, 19 mm	Premium Reverse Cutting	SY426H	12ea/bx
3-0 PGA-PCL, 2 Metric, 30", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SY427H	12ea/bx
5-0 PGA-PCL, 1 Metric, 18", Undyed, P-13, 3/8, 13 mm	Premium Reverse Cutting	SY493GX	12ea/bx
4-0 PGA-PCL, 1.5 Metric, 18", Undyed, P-13, 3/8, 13 mm	Premium Reverse Cutting	SY494GX	12ea/bx
4-0 PGA-PCL, 1.5 Metric, 18", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SY496GX	12ea/bx
3-0 PGA-PCL, 2 Metric, 18", Undyed, P-12, 3/8, 19 mm	Premium Reverse Cutting	SY497GX	12ea/bx
MONOFILAMENT ABSORBABLE SUTURE			
4-0 Monofilament, 1.5 Metric, 18", Clear, P-13, 3/8, 13 mm	Premium Reverse Cutting	SZ494G	36ea/bx
3-0 Monofilament, 2 Metric, 18", Clear, P-12, 3/8, 19 mm	Premium Reverse Cutting	SZ497G	36ea/bx