

Test Tubes

McKesson test tubes are molded from premium-grade polystyrene (PS) or polypropylene (PP). They are both non-graduated, translucent tubes with round bottoms.

Features and Benefits

MFR# 177-110410

- Molded from premium-grade polystyrene (PS) for clarity
- Can withstand 1500 RPM during centrifugation
- Brittleness temperature: 20°C
- Maximum temperature: 90°C
- Not made with natural rubber latex

MFR# 177-110441

- Molded from premium-grade polypropylene (PP) for impact strength and stability to common acids, alkalies and solvents
- Autoclavable, freezable and can withstand 3000 RPM during centrifugation
- Brittleness temperature: 20°C
- Maximum temperature: 90°C
- Not made with natural rubber latex



MFR# 177-110410

MFR# 177-110441

ORDERING INFORMATION

TEST TUBES			
MFR #	DESCRIPTION	SPECS	QTY
177-110410	Polystyrene (PS), Plastic	12 mm x 75 mm, 5 mL	1000/bag
177-110441	Polypropylene (PP), Plastic	12 mm x 75 mm, 5 mL	1000/bag

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2018 McKesson Medical-Surgical Inc.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908