

Lab Caps and Closures

McKesson's line of laboratory supplies includes caps and closures to recap round cuvettes, centrifuge tubes, blood drawing tubes, glass test tubes, plastic culture tubes and more. McKesson's laboratory caps are color-coded for easy identification. They protect samples from cross-contamination and help safeguard laboratory areas against infections and spills.

Features and Benefits

Flanged Plug Caps

- Provide exceptional leak resistance for 12 mm test tubes, centrifuge tubes and round cuvettes
- Feature two flexible flanges that press firmly against the tube's inner wall for enhanced leak resistance compared to standard plug caps
- Not made with natural rubber latex



Universal Snap Caps

- Ideal for recapping most 12 mm, 13 mm and 16 mm evacuated glass blood collection tubes
- Designed for easy-on and easy-off use, featuring double thumb tabs
- Flexible closures protect from aerosols of highly infectious microorganisms
- Precision-molded, low-density polyethylene (LDPE) designed to fit firmly on the tube and remain sturdy under most rigorous laboratory procedures, including centrifugation
- Color-coded for easy sample identification
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Lab Caps and Closures

Features and Benefits

Snap Caps

- Effective solution for recapping blood drawing tubes, glass test tubes or plastic culture tubes
- Precision-molded, low-density polyethylene (LDPE) to maintain a tight seal during centrifugation, agitation and refrigeration procedures
- Color-coded for easy sample identification
- Tabbed design for easy cap removal
- Single-thumb tab design allows for the tubes to be placed closely in test tube racks
- Not made with natural rubber latex



ORDERING INFORMATION

| MFR # | PRODUCT | MATERIAL | COLOR | SIZE | TAB | QTY |
|-------------|-------------------|--------------|----------|-----------------------------|------------------|----------|
| 177-118127C | Flanged Plug Caps | Polyethylene | N/A | 12 mm | N/A | 1000/bag |
| 177-118115A | Snap Caps | Polyethylene | Gray | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115B | Snap Caps | Polyethylene | Blue | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115C | Snap Caps | Polyethylene | Natural | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115G | Snap Caps | Polyethylene | Green | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115L | Snap Caps | Polyethylene | Lavender | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115R | Snap Caps | Polyethylene | Red | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-118115Y | Snap Caps | Polyethylene | Yellow | Universal, 12, 13, or 16 mm | Double Thumb Tab | 1000/bag |
| 177-113146A | Snap Caps | Polyethylene | Gray | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113146C | Snap Caps | Polyethylene | Natural | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113146L | Snap Caps | Polyethylene | Lavender | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113146R | Snap Caps | Polyethylene | Red | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113146W | Snap Caps | Polyethylene | White | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113146Y | Snap Caps | Polyethylene | Yellow | 13 mm | Single Thumb Tab | 1000/bag |
| 177-113148B | Snap Caps | Polyethylene | Blue | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148C | Snap Caps | Polyethylene | Natural | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148G | Snap Caps | Polyethylene | Green | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148L | Snap Caps | Polyethylene | Lavender | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148R | Snap Caps | Polyethylene | Red | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148W | Snap Caps | Polyethylene | White | 16 mm | Single Thumb Tab | 1000/bag |
| 177-113148Y | Snap Caps | Polyethylene | Yellow | 16 mm | Single Thumb Tab | 1000/bag |

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908