


Single Sterile Cover Drapes

McKesson single sterile cover drapes are protective barrier covers intended to be used to cover surgical equipment and provide a protective barrier for that equipment. They deliver optimal protection and performance for creating and maintaining a sterile field in surgery. The covers are individually packaged and offer single sterile use with consistent quality.

Features and Benefits

- PE film-based drapes used to provide a sterile barrier over back tables
- Resistant to tears, punctures, and fiber strains under normal conditions of use
- Flat closure
- Low-linting SMS fabric properties reduces particle generation
- Fan-folded drape designs with clear instructional stamps promote ease of use, aseptic presentation, and intuitive placement
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	ATTRIBUTES	SIZE	QTY
16-i80-13101-S	Back Table Cover, Heavy Duty	Zone-Reinforced, Sterile	44 in x 75 in	1 pk, 20 pk/cs
16-i80-13103-S	Back Table Cover, Heavy Duty	Zone-Reinforced, Sterile	50 in x 90 in	1 pk, 20 pk/cs
16-i80-13310-S	Back Table Cover	Zone-Reinforced, Sterile	44 in x 90 in	1 pk, 20 pk/cs
16-i80-13105-S	Mayo Stand Cover	SMS Reinforced, Sterile	23 in x 54 in	1 pk, 30 pk/cs
16-i80-13106-S	Mayo Stand Cover	SMS Reinforced, Sterile	29.5 in x 57 in	1 pk, 30 pk/cs

Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.