

Suction Canisters and Liners

McKesson disposable suction canisters are single-use, designed for convenience by removing the need for a secondary waste liner. The disposable canisters are available in 800 mL, 1200 mL, 2000 mL and 3300 mL sizes. Canisters include a integrated filter shut-off valve that supports best practice by preventing back flow of aerosols and fluids into the vacuum system.

The reusable canisters help manage suction waste by providing simple, safe and complete solutions. The 1000 mL canister and liner have been developed for low-volume waste fluid collection. The small, compact size helps maximize space efficiency and provides a cost-effective solution for portable suction systems.

Features and Benefits

Disposable Suction Canisters

- Lid attaches firmly to the canister, maximizing vacuum power and minimizing risk of leakage
- Easy-to-read volume marks
- Reflux protector on patient port reduces contamination exposure risk
- Shut-off valve includes a 0.3 micron bacterial eliminating filter effective in removing 99.9% of pseudomonas
- Non-Sterile, disposable
- Not made with natural rubber latex



Reusable Suction Canisters

- Made of polycarbonate material, Canisters are designed with durability and efficiency in mind
- Easy-to-read volume marks
- Large patient label strip for convenient documentation
- Canister “T” designates vacuum source connection to simplify canister set up
- Color coded by size
- Non-Sterile, reusable
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Suction Canisters and Liners

Mckesson liners are disposable and available in 1000 mL, 1500 mL, 2000 mL and 3000 mL sizes. Permanently attached lids and a multifunctional lid tube help decrease the risk of fluid exposure and improve ease-of-use. The disposable liner sizes are to be used with the same size Receptal® Reusable Canisters. The disposable liners with integrated Vac-Gard® filter shut-off valve support best practice by preventing back flow of aerosols and fluids into the vacuum system.

Features and Benefits

Liners

- 100% leak tested seamless liners minimize risk of contamination
- Permanently attached lids and a multifunctional lid tube help decrease the risk of fluid exposure and improve ease-of-use
- Integrated rotating lid tube promotes simple set-up
- Multiple liner options to meet clinical and safety needs
- Anti-foaming emulsion prevents froth formation and aerosolization, reducing the risk of contaminant exposure
- Non-Sterile, Disposable
- Latex Free and BPA Free



Reusable Suction Canister			Corresponding Disposable Liner		
MFR #	DESCRIPTION	SIZE	MFR #	DESCRIPTION	SIZE
16-43423-01	Reusable Suction Canister	1500 mL	16-43044-01	VAC-GARD® Pour Spout Suction Liner	1500 mL
16-43449-05	Reusable Suction Canister	1000 mL	16-43010-05	VAC-GARD® Pour Spout Suction Liner	1000 mL
16-43445-01	Reusable Suction Canister	2000 mL	16-43044-05	VAC-GARD® Pour Spout Suction Liner	2000 mL

ORDERING INFORMATION

MFR #	DESCRIPTION	SIZE	QTY
16-43203-01	Suction Canister, Disposable, Non-Sterile	1200 mL	1 ea, 8 ea/bx, 4 bx/cs
16-43208-01	Suction Canister, Disposable, Non-Sterile	800 mL	72 ea/cs
16-43208-05	Suction Canister, Disposable, Non-Sterile	800 mL	1 ea, 8 ea/bx, 4 bx/cs
16-43206-01	Suction Canister, Disposable, Non-Sterile	2000 mL	1 ea, 3 ea/bx, 4 bx/cs
16-43423-01	Suction Canister, Reusable, Non-Sterile	1500 mL	1 ea/bag, 10 bags/cs
16-43449-05	Suction Canister, Reusable, Non-Sterile	1000 mL	1 ea/bag, 10 bags/cs
16-43445-01	Suction Canister, Reusable, Non-Sterile	2000 mL	1 ea/bag, 10 bags/cs
16-43044-01	Pour Spout Suction Liner, Non-Sterile	1500 mL	10 ea/bag, 5 bags/cs
16-43010-05	Pour Spout Suction Liner, Non-Sterile	1000 mL	10 ea/bag, 5 bags/cs
16-43044-05	Pour Spout Suction Liner, Non-Sterile	2000 mL	10 ea/bag, 5 bags/cs

McKesson Medical-Surgical
 9954 Mayland Drive
 Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
 800.347.2456

Home Medical Equipment
 888.822.8111

Primary Care
 866.625.2679

McKesson Brands Support
 800.777.4908