Intermittent Catheters

McKesson intermittent catheters minimize urethral trauma and reduce the chance of bringing bacteria in to the bladder. The ready-to-use catheters are designed to help quickly and easily self-catheterize. They're offered in straight and Coudé tips and female or male/universal lengths. Available in the following French sizes: 10, 12, 14, 16, 18 Fr/FR

McKesson hydrophilic catheters are the preferred cathing system for those who experience the sensation of pain, pressure, or discomfort during the cathing process. Hydrophilic catheters come pre-lubricated and always optimally hydrated to promote excellent gliding properties that ensure easier, more comfortable insertion. Hydrophilic catheters have been shown to lead to better patient satisfaction with this type of technique by reducing the risk of urethral damage due to decreased friction. Available in the following French sizes: 10, 12, 14, 16, 18 Fr/FR

Intermittent Catheters

Intermittent Catheters

HYDROPHILIC | STERILE WATER PACKET

Features and Benefits

Uncoated Catheters

- Offered in straight or Coudé tips
- Funnel end for ease of handling (except MFR# 16-F614L)
- · Color-coded to indicate French size
- Available in 10, 12, 14, 16, 18 Fr/FR
- Non-coated
- Not made with natural rubber latex

Hydrophilic Catheters

- Offered in straight or Coudé tips
- Packaged with sterile water
- Coated with a hydrophilic layer consisting of PVP/ poly(vinyl pyrrolidone) and sodium chloride/NaCl
- Color-coded to indicate French size
- Available in 10, 12, 14, 16, 18 Fr/FR
- Hydrophilic
- Not made with natural rubber latex



14 Fr

FEMALE



Coudé and Straight tips



12 Fr

MER # 16-FH612

Funnel and Luer ends

MSKESSON

ORDERING INFORMATION

UNCOATED				
MFR #	DESCRIPTION	COLOR	QTY	
16-M1610	Intermittent Catheter, Uncoated, Male/Universal, Straight Tip, PVC, Sterile, 16", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F610	Intermittent Catheter, Uncoated, Female, Straight Tip, PVC, Sterile, 6", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1610C	Intermittent Catheter, Uncoated, Male/Universal, Coudé Tip, PVC, Sterile, 16", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1612	Intermittent Catheter, Uncoated, Male/Universal, Straight Tip, PVC, Sterile, 16", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F612	Intermittent Catheter, Uncoated, Female, Straight Tip, PVC, Sterile, 6", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1612C	Intermittent Catheter, Uncoated, Male/Universal, Coudé Tip, PVC, Sterile, 16", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1614	Intermittent Catheter, Uncoated, Male/Universal, Straight Tip, PVC, Sterile, 16", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F614	Intermittent Catheter, Uncoated, Female, Straight Tip, PVC, Sterile, 6", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1614C	Intermittent Catheter, Uncoated, Male/Universal, Coudé Tip, PVC, Sterile, 16", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F614L	Intermittent Catheter, Uncoated, Female, Luer End, Straight Tip, PVC, Sterile, 6", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-U1014	Intermittent Catheter, Uncoated, Universal, Straight Tip, PVC, Sterile, 10", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1616	Intermittent Catheter, Uncoated, Male/Universal, Straight Tip, PVC, Sterile, 16", 16 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F616	Intermittent Catheter, Uncoated, Female, Straight Tip, PVC, Sterile, 6", 16 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1616C	Intermittent Catheter, Uncoated, Male/Universal, Coudé Tip, PVC, Sterile, 16", 16 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-F618	Intermittent Catheter, Uncoated, Female, Straight Tip, PVC, Sterile, 6", 18 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1618C	Intermittent Catheter, Uncoated, Male/Universal, Coudé Tip, PVC, Sterile, 16", 18 Fr		1 ea, 30 ea/bx, 10 bx/cs	
16-M1618	Intermittent Catheter, Uncoated, Male/Universal, Straight Tip, PVC, Sterile, 16", 18 Fr		1 ea, 30 ea/bx, 10 bx/cs	

MSKESSON

ORDERING INFORMATION

	HYDROPHILIC		
MFR #	DESCRIPTION	COLOR	QTY
16-FH810	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Female, Straight Tip, PVC, Sterile, 8", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1610C	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Coudé Tip, PVC, Sterile, 16", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1610	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Straight Tip, PVC, Sterile, 16", 10 Fr		1 ea, 30 ea/bx, 10 bx/cs
16-FH612	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Female, Straight Tip, PVC, Sterile, 6", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1612C	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Coudé Tip, PVC, Sterile, 16", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1612	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Straight Tip, PVC, Sterile, 16", 12 Fr		1 ea, 30 ea/bx, 10 bx/cs
16-FH612	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Female, Straight Tip, PVC, Sterile, 6", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs
16-FH814	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Female, Straight Tip, PVC, Sterile, 8", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1614C	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Coudé Tip, PVC, Sterile, 16", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs
16-MH1614	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Straight Tip, PVC, Sterile, 16", 14 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1616C	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Coudé Tip, PVC, Sterile, 16", 16 Fr		1 ea, 30 ea/bx, 10 bx/cs
16-MH1616	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Straight Tip, PVC, Sterile, 16", 16 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1618C	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Coudé Tip, PVC, Sterile, 16", 18 Fr		1 ea, 30 ea/bx, 10 bx/cs
16- MH1618	Intermittent Catheter, Hydrophilic, Sterile Water Packet, Male/Universal, Straight Tip, PVC, Sterile, 16", 18 Fr		1 ea, 30 ea/bx, 10 bx/cs

MSKESSON

McKesson Medical-Surgical 9954 Mayland Drive Richmond, VA 23233 mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health 800.347.2456

Home Medical Equipment 888.822.8111

Primary Care 866.625.2679

McKesson Brands Support 800.777.4908

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2018 McKesson Medical-Surgical Inc.