

Nasopharyngeal Airways

The McKesson nasopharyngeal airways are made from a soft, flexible neoprene material that helps prevent damage to the nasal passages. Our nasopharyngeal airway helps maximize patient comfort while offering the best solution for latex sensitivity.

Features and Benefits

Robertazzi Nasopharyngeal Airways

- For quick and safe airway patency without intubation
- Flexible to aid in patient comfort
- Blunt tip
- Available in nine sizes to accommodate a variety of airway types
- Individually wrapped, sterile packaged
- Robertazzi nasopharyngeal airway kits are ideal for crash carts
- HCPCS A7044
- Not made with natural rubber latex



PVC Nasopharyngeal Airways

- For quick and safe airway patency without intubation
- Flexible to aid in patient comfort
- Beveled tip helps with gentle insertion
- Smooth PVC surface helps prevent damage to nasal passage
- Standard 15 mm flange fitting
- Individually wrapped - sterile packaged
- HCPCS A7044
- Not made with natural rubber latex



Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

ROBERTAZZI NASOPHARYNGEAL AIRWAY						
MFR #	DESCRIPTION	ATTRIBUTES	INNER DIAMETER	OUTER DIAMETER	LENGTH	QTY
159-75-20	Robertazzi Nasopharyngeal Airway, 20Fr	Blunt Tip, Sterile	5.0 mm	6.7 mm	115 mm	10 ea/bx
159-75-22	Robertazzi Nasopharyngeal Airway, 22Fr	Blunt Tip, Sterile	5.5 mm	7.3 mm	120 mm	10 ea/bx
159-75-24	Robertazzi Nasopharyngeal Airway, 24Fr	Blunt Tip, Sterile	6.0 mm	7.6 mm	120 mm	10 ea/bx
159-75-26	Robertazzi Nasopharyngeal Airway, 26Fr	Blunt Tip, Sterile	6.5 mm	8.6 mm	130 mm	10 ea/bx
159-75-28	Robertazzi Nasopharyngeal Airway, 28Fr	Blunt Tip, Sterile	7.0 mm	9.3 mm	135 mm	10 ea/bx
159-75-30	Robertazzi Nasopharyngeal Airway, 30Fr	Blunt Tip, Sterile	7.5 mm	10.0 mm	135 mm	10 ea/bx
159-75-32	Robertazzi Nasopharyngeal Airway, 32Fr	Blunt Tip, Sterile	8.0 mm	10.7 mm	145 mm	10 ea/bx
159-75-34	Robertazzi Nasopharyngeal Airway, 34Fr	Blunt Tip, Sterile	8.5 mm	11.2 mm	160 mm	10 ea/bx
159-75-36	Robertazzi Nasopharyngeal Airway, 36Fr	Blunt Tip, Sterile	9.0 mm	12.0 mm	175 mm	10 ea/bx

PVC NASOPHARYNGEAL AIRWAY					
MFR #	DESCRIPTION	ATTRIBUTES	INNER DIAMETER	LENGTH	QTY
159-PVC-14	PVC Nasopharyngeal Airway, 14Fr	Beveled Tip, Clear, Sterile	3.0 mm	90 mm	10 ea/bx
159-PVC-16	PVC Nasopharyngeal Airway, 16Fr	Beveled Tip, Clear, Sterile	3.5 mm	105 mm	10 ea/bx
159-PVC-18	PVC Nasopharyngeal Airway, 18Fr	Beveled Tip, Clear, Sterile	4.0 mm	115 mm	10 ea/bx
159-PVC-20	PVC Nasopharyngeal Airway, 20Fr	Beveled Tip, Clear, Sterile	5.0 mm	120 mm	10 ea/bx
159-PVC-22	PVC Nasopharyngeal Airway, 22Fr	Beveled Tip, Clear, Sterile	5.5 mm	125 mm	10 ea/bx
159-PVC-24	PVC Nasopharyngeal Airway, 24Fr	Beveled Tip, Clear, Sterile	6.0 mm	140 mm	10 ea/bx
159-PVC-26	PVC Nasopharyngeal Airway, 26Fr	Beveled Tip, Clear, Sterile	6.5 mm	145 mm	10 ea/bx
159-PVC-28	PVC Nasopharyngeal Airway, 28Fr	Beveled Tip, Clear, Sterile	7.0 mm	155 mm	10 ea/bx
159-PVC-30	PVC Nasopharyngeal Airway, 30Fr	Beveled Tip, Clear, Sterile	7.5 mm	160 mm	10 ea/bx
159-PVC-32	PVC Nasopharyngeal Airway, 32Fr	Beveled Tip, Clear, Sterile	8.0 mm	180 mm	10 ea/bx
159-PVC-34	PVC Nasopharyngeal Airway, 34Fr	Beveled Tip, Clear, Sterile	8.5 mm	185 mm	10 ea/bx
159-PVC-36	PVC Nasopharyngeal Airway, 36Fr	Beveled Tip, Clear, Sterile	9.0 mm	190 mm	10 ea/bx



WARNING: This product can expose you to chemicals including Chloroprene which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

ORDERING INFORMATION

ROBERTAZZI NASOPHARYNGEAL AIRWAY KITS				
MFR #	DESCRIPTION	ATTRIBUTES	CONTENTS	QTY
159-75-04	Robertazzi Nasopharyngeal Airway Kit	Blunt Tip, Sterile	1 Nasopharyngeal Airway, 20 Fr 1 Nasopharyngeal Airway, 24 Fr 1 Nasopharyngeal Airway, 28 Fr 1 Nasopharyngeal Airway, 32 Fr 4 Packs of Lubricating Jelly, Sterile	1 kit/bx
159-75-05	Robertazzi Nasopharyngeal Airway Kit	Blunt Tip, Sterile	1 Nasopharyngeal Airway, 28 Fr 1 Nasopharyngeal Airway, 30 Fr 1 Nasopharyngeal Airway, 32 Fr 1 Nasopharyngeal Airway, 34 Fr 1 Nasopharyngeal Airway, 36 Fr 5 Packs of Lubricating Jelly, Sterile	1 kit/bx
159-75-06	Robertazzi Nasopharyngeal Airway Kit	Blunt Tip, Sterile	1 Nasopharyngeal Airway, 22 Fr 1 Nasopharyngeal Airway, 24 Fr 1 Nasopharyngeal Airway, 26 Fr 1 Nasopharyngeal Airway, 28 Fr 1 Nasopharyngeal Airway, 30 Fr 1 Nasopharyngeal Airway, 32 Fr 6 Packs of Lubricating Jelly, Sterile	1 kit/bx
159-75-09	Robertazzi Nasopharyngeal Airway Kit	Blunt Tip, Sterile	1 Nasopharyngeal Airway, 20 Fr 1 Nasopharyngeal Airway, 22 Fr 1 Nasopharyngeal Airway, 24 Fr 1 Nasopharyngeal Airway, 26 Fr 1 Nasopharyngeal Airway, 28 Fr 1 Nasopharyngeal Airway, 30 Fr 1 Nasopharyngeal Airway, 32 Fr 1 Nasopharyngeal Airway, 34 Fr 1 Nasopharyngeal Airway, 36 Fr 9 Packs of Lubricating Jelly, Sterile	1 kit/bx



WARNING: This product can expose you to chemicals including Chloroprene which is known to the State of California to cause cancer. For more information go to www.P65Warnings.ca.gov.

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908