

McKESSON



Blood Collection Tube Holders

McKesson Prevent® HT blood collection tube holders, available with and without safety hinges, are compatible with common needle sizes to be used for routine blood collection procedures. The tube holders complete the blood collection assembly when combined with appropriate needles and evacuated blood collection tubes to reduce risk of back-end needle stick injuries and improve overall compliance with OSHA single-use holder guidelines.

FEATURES AND BENEFITS

Universal Features

- Non-stackable holders in bulk packaging
- Universal design with threaded female adapter compatible with common needle sizes and fits both 13 mm and 16 mm blood collection tubing
- Clear polypropylene construction allows unobstructed view of tube and blood as it is drawn
- Raised edge tube walls and tabbed end for secure grip
- Not made with natural rubber latex

Additional Features for

Tube Holders with Safety Hinge

- Holder-based safety hinge provides protection on both ends and prevents re-use to promote OSHA compliance
- Immediate one-handed safety feature activation, confirmed with an audible and tangible "click," reduces exposure time and helps protect from accidental needle stick injuries
- Used to reduce risk of back-end needle stick injuries and improve overall compliance with OSHA single-use holder guidelines



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

McKESSON



Blood Collection Tube Holders

MFR#	DESCRIPTION	For Use With	QTY
16-BCTHLDR	Tube Holders	Safety Needles	25/bags, 40 bags/case
16-BCTHLDRSF	Tube Holders with Safety Hinge	Non-Safety Needles	250/bags, 4 bags/case

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908