

Urine Specimen Container

McKesson urine specimen container provides convenient urine collection and sample transfer in a streamlined closed system. It allows urine to be transferred into one or more vacuum tubes without exposure to the sample. The closed system prevents healthcare professionals from having to open the collection cup and pour the specimen, reducing leakage, spills, and contamination. The cup features graduation marks in both ounces and milliliters: up to 4 ounces in increments of 0.5 ounces and up to 110 milliliters in increments of 10 milliliters. It includes a patient safety label to cover the fixed needle.

Features and Benefits

- Features a leak-resistant screw cap to minimize risk of contamination and exposure to specimen
- Closed system transfer of urine eliminates the need for pouring and thus avoids spills
- Works with most vacuum urine tubes
- Includes temperature strip to help prevent adulteration
- Needle is fixed to transfer device for healthcare staff and patient safety
- Frosted panel for patient information containing the following fields: Name, Room Number, Doctor, Date, Time
- Sterile, individually bagged
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	CAPACITY	QTY
16-UCC4	Urine Specimen Container with Integrated Transfer Device and Temperature Strip	4 fl oz	1 ea/bg, 200 bgs/cs

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.