


McKesson Brands Catheters

Easy to Use. Designed for Comfort. Wide Selection. Affordable Prices.



Why is the urology category important?

Why McKesson?

Every day, you face multiple challenges when it comes to running a business and serving your patients. McKesson Medical-Surgical is here to help. At McKesson, our focus is always on your success and meeting your patients' needs. That's why our expert team can help to make sure you have the right products, right when you need them.



Between
15%
and
25%

of hospitalized patients receive short-term indwelling catheters



14%
growth
(over previous year)



\$600
million

across all segments



An estimated
5%
of long-term care facility residents use catheters

Source: HIDA Home Care reports, CDC, 2017 and 2018; GHX reporting 2019



McKesson supports the health of your business by offering:



More products, more choice

Quality and selection with more national brands from a single source



Distribution you can count on

Reliable and flexible delivery options including patient home delivery



Solutions to enhance your business

Automated ordering solutions and resources to achieve seamless billing and reimbursements



Reduce readmissions

Reduce hospital readmissions with our comprehensive clinical programs and product offerings



99% fill rate

We take pride in ensuring the right products get to you and your patients



Educate clinical staff

Help your staff achieve patient and resident care goals with specialized clinical training and education

An introduction to urinary catheters

A urinary catheter is a hollow, flexible tube that collects urine from the bladder, leading to a drainage bag.



Common reasons for catheter use

- Diagnosis of urinary incontinence or urinary retention
- For a brief time following major surgery
- After diabetes complications
- Spinal injury
- After prostate surgeries
- Conditions like a stroke or multiple sclerosis

Catheter types


- Foley Balloon
- Intermittent
- External

| Foley HCPCS Codes | Intermittent HCPCS Codes | External HCPCS Codes |
|--|---|---|
| A4311-A4316 – are kits that include an indwelling catheter | A4351 – Intermittent straight tip with or without coating | A4326 – Male external catheter with collection chamber |
| A4338 – Indwelling catheter with two-way latex with coating (silicone, silicone elastomer, hydrophilic, etc.) | A4352 – Intermittent Coudé tip with or without coating | A4349 – Male external catheter |
| A4340 – Indwelling catheter, specialty type (Coudé, mushroom, wing, etc.) | A4353 – Intermittent with insertion supplies (this could be a closed system, introducer tip, etc.) | |
| A4344 – Indwelling catheter, two-way all silicone | | |
| A4346 – Indwelling three-way (continuous irrigation) | | |

An introduction to urinary catheters (continued)


You can be confident in the products you offer to your patients when you purchase medical supplies from McKesson. We combine the quality of leading manufacturers, with a fast, nationwide distribution network to help you run a smoother, more efficient business – all backed by a 100% satisfaction guarantee.

French Catheter Sizing



Actual Size

Catheter Sizing Continuum



Catheter sizes measured in French scale

All McKesson Brand catheters feature:



Polished eyes to help minimize urethral trauma



Color coded packaging to indicate French size



Variety of sizes and tips



Not made with DEHP, DINP, or BPA



Individually packaged for **single sterile use**

PVC Uncoated Catheters

McKesson intermittent catheters help minimize urethral trauma and reduce the chance of bringing bacteria into the bladder. These PVC uncoated catheters are ready-to-use and designed to help your patients self-catheterize quickly and easily.

Features

- Staggered eyes for maximum urine flow
- Polished eyes to help minimize urethral trauma
- Flexible funnel end with ridged design for grip and ease of handling
- Color-coded funnel and packaging label for easy identification of French size
- Packaging design includes loop cutout for hanging to maintain sterility
- Individually packaged for single, sterile use
- Not made with DEHP, DINP, BPA or natural rubber latex



| DESCRIPTION | SIZE | MFR # | ITEM # | ITEM # (BX) |
|---|-------|-----------|---------|-------------|
| McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 6" | 10FR | 16-F610 | 1020810 | 31661900 |
| | 12FR | 16-F612 | 1020812 | 31621900 |
| | 14FR | 16-F614 | 1020814 | 36641900 |
| | 16FR | 16-F616 | 1020818 | 56611900 |
| | 18 FR | 16-F618 | 1135715 | 36981900 |
| McKesson Intermittent Catheter, Uncoated PVC, Luer Tip, 6" | 14FR | 16-F614L | 1020816 | 35641900 |
| McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 10" | 14FR | 16-U1014 | 1135717 | 37331900 |
| McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 16" | 10FR | 16-M1610 | 1020809 | 69611900 |
| | 12FR | 16-M1612 | 1020811 | 62121900 |
| | 14FR | 16-M1614 | 1020813 | 66641900 |
| | 16FR | 16-M1616 | 1020817 | 69661900 |
| | 18 FR | 16-M1618 | 1135716 | 37071900 |
| McKesson Intermittent Catheter, Uncoated PVC, Coude Tip, 16" | 10FR | 16-M1610C | 1135718 | 75701900 |
| | 12FR | 16-M1612C | 1135719 | 75711900 |
| | 14FR | 16-M1614C | 1020815 | 32641900 |
| | 16FR | 16-M1616C | 1020819 | 23661900 |
| | 18 FR | 16-M1618C | 1135720 | 37081900 |

PVC Hydrophilic Catheters

McKesson hydrophilic catheters are pre-lubricated and hydrated to help promote excellent gliding properties. With easier, more comfortable insertion, hydrophilic catheters offer decreased friction. This leads to better patient satisfaction by helping reduce the risk of urethral damage.

Features

- Staggered eyes for maximum urine flow
- Polished eyes help minimize urethral trauma
- Flexible funnel end with ridged design for grip and ease of handling
- Color-coded funnel and packaging label for easy identification of French size
- Packaging design includes loop cutout for hanging to maintain sterility
- Individually packaged for single, sterile use
- Not made with DEHP, DINP, BPA or natural rubber latex



| DESCRIPTION | SIZE | MFR # | ITEM # | ITEM # (BX) |
|---|-------|-----------|---------|-------------|
| McKesson Intermittent Catheter, Hydrophilic, Straight Tip, 6" | 12FR | 16-FH612 | 1135721 | 36991900 |
| | 14FR | 16-FH614 | 1135722 | 75691900 |
| McKesson Intermittent Catheter, Uncoated PVC, Luer Tip, 6" | 10FR | 16-FH810 | 1135723 | 37021900 |
| | 14FR | 16-FH814 | 1135724 | 56391900 |
| McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 10" | 10FR | 16-MH1610 | 1135725 | 37091900 |
| | 12FR | 16-MH1612 | 1135726 | 37131900 |
| | 14FR | 16-MH1614 | 1135727 | 37171900 |
| | 16FR | 16-MH1616 | 1135728 | 37251900 |
| | 18 FR | 16-MH1618 | 1135729 | 37291900 |
| | 10FR | 16-M1610 | 1020809 | 69611900 |
| McKesson Intermittent Catheter, Hydrophilic, Coude Tip, 16" | 12FR | 16-M1612 | 1020811 | 62121900 |
| | 14FR | 16-M1614 | 1020813 | 66641900 |
| | 16FR | 16-M1616 | 1020817 | 69661900 |
| | 18 FR | 16-M1618 | 1135716 | 37071900 |
| | | | | |

Silicone Foley Catheters

A Foley is a catheter with two channels and one tube. The first channel drains urine into a collection bag, while the other channel has a valve on the outside end and connects to a balloon at the inside tip. The balloon is inflated with sterile water when it lies inside the bladder to prevent the catheter from slipping out.

A Foley catheter is meant to be in place for no longer than three months. Traditional straight catheters are inserted once and discarded when the bladder has been emptied.

McKesson 2-Way 100% Silicone Foley Catheters are ideal for patients with latex allergies and longer term catheterization. Silicone provides chemical and thermal stability, low surface tension and hydrophobicity.

Features

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra



| DESCRIPTION | SIZE | MFR # | ITEM # | ITEM # (BX) |
|---|------|-------------|---------|-------------|
| McKesson Foley Catheter, Silicone, Straight Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume | 12FR | 16-SFOL512 | 1152043 | 65121910 |
| | 14FR | 16-SFOL514 | 1152044 | 65141910 |
| | 16FR | 16-SFOL516 | 1152045 | 65161910 |
| | 18FR | 16-SFOL518 | 1152046 | 65181910 |
| | 20FR | 16-SFOL520 | 1152047 | 65201910 |
| | 22FR | 16-SFOL522 | 1152048 | 65221910 |
| | 24FR | 16-SFOL524 | 1152049 | 65241920 |
| McKesson Foley Catheter, Silicone, Straight Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume | 16FR | 16-SFOL3016 | 1152038 | 63161910 |
| | 18FR | 16-SFOL3018 | 1152039 | 63181910 |
| | 20FR | 16-SFOL3020 | 1152040 | 63201910 |
| | 22FR | 16-SFOL3022 | 1152041 | 63221910 |
| | 24FR | 16-SFOL3024 | 1152042 | 63241910 |
| McKesson Foley Catheter, Silicone, Coude Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume | 12FR | 16-SFOL512C | 1152031 | 15121910 |
| | 14FR | 16-SFOL514C | 1152032 | 15141910 |
| | 16FR | 16-SFOL516C | 1152033 | 15161910 |
| | 18FR | 16-SFOL518C | 1152034 | 15181910 |
| | 20FR | 16-SFOL520C | 1152035 | 15201910 |
| | 22FR | 16-SFOL522C | 1152036 | 15221910 |

Hydrophilic-Coated Latex Foley Catheters

McKesson hydrophilic-coated latex Foley catheters are retention type urological catheters made of silicone-coated natural latex material with the addition of a lubricious/hydrophilic coating to facilitate the insertion of the catheter.

Features

- Hydrophilic coating helps urethral insertion and helps reduce irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Contains natural rubber latex



| DESCRIPTION | SIZE | MFR # | ITEM # | ITEM # (BX) |
|---|------|--------------|---------|-------------|
| McKesson Foley Catheter, Hydrophilic Coated Latex, Straight Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume | 12FR | 16-LFOL512 | 1152068 | 51211910 |
| | 14FR | 16-LFOL514 | 1152069 | 51411910 |
| | 16FR | 16-LFOL516 | 1152070 | 16151910 |
| | 18FR | 16-LFOL518 | 1152071 | 18151910 |
| | 20FR | 16-LFOL520 | 1152072 | 52011910 |
| | 22FR | 16-LFOL522 | 1152073 | 52211910 |
| | 24FR | 16-LFOL524 | 1152074 | 52411910 |
| McKesson Foley Catheter, Hydrophilic Coated Latex, Straight Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume | 14FR | 16-LFOL3014 | 1152062 | 30141910 |
| | 16FR | 16-LFOL3016 | 1152063 | 31661910 |
| | 18FR | 16-LFOL3018 | 1152064 | 31861910 |
| | 20FR | 16-LFOL3020 | 1152065 | 32061910 |
| | 22FR | 16-LFOL3022 | 1152066 | 32261910 |
| | 24FR | 16-LFOL3024 | 1152067 | 32461910 |
| McKesson Foley Catheter, Hydrophilic Coated Latex, Coude Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume | 12FR | 16-LFOL512C | 1152055 | 15621910 |
| | 14FR | 16-LFOL514C | 1152056 | 51461910 |
| | 16FR | 16-LFOL516C | 1152057 | 51661910 |
| | 18FR | 16-LFOL518C | 1152058 | 51861920 |
| | 20FR | 16-LFOL520C | 1152059 | 52061910 |
| | 22FR | 16-LFOL522C | 1152060 | 52261920 |
| | 24FR | 16-LFOL524C | 1152061 | 52461910 |
| McKesson Foley Catheter, Hydrophilic Coated Latex, Coude Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume | 16FR | 16-LFOL3016C | 1152050 | 30161920 |
| | 18FR | 16-LFOL3018C | 1152051 | 30181910 |
| | 20FR | 16-LFOL3020C | 1152052 | 30201920 |
| | 22FR | 16-LFOL3022C | 1152053 | 22031920 |
| | 24FR | 16-LFOL3024C | 1152054 | 30241920 |

Silicone Coated Latex Foley Catheters

McKesson two-way siliconized latex Foley catheters are durable and flexible to help improve patient comfort. Silicone is typically resistant to chemical attack and insensitive to temperature changes.

Features

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Contains natural rubber latex



| DESCRIPTION | SIZE | MFR # | ITEM # | ITEM # (BX) |
|---|------|--------------|---------|-------------|
| McKesson Foley Catheter, Silicone Coated Latex, Straight Tip, 5cc Balloon Size, 5-15cc Inflation Volume | 12FR | 16-SLFOL512 | 1152082 | 16521910 |
| | 14FR | 16-SLFOL514 | 1152083 | 16541910 |
| | 16FR | 16-SLFOL516 | 1152084 | 65611920 |
| | 18FR | 16-SLFOL518 | 1152085 | 85611920 |
| | 20FR | 16-SLFOL520 | 1152086 | 25611910 |
| | 22FR | 16-SLFOL522 | 1152087 | 22561910 |
| | 24FR | 16-SLFOL524 | 1152088 | 42561910 |
| McKesson Foley Catheter, Silicone Coated Latex, Straight Tip, 30cc Balloon Size, 30-40cc Inflation Volume | 12FR | 16-SLFOL3012 | 1152075 | 16121910 |
| | 14FR | 16-SLFOL3014 | 1152076 | 16141910 |
| | 16FR | 16-SLFOL3016 | 1152077 | 16611910 |
| | 18FR | 16-SLFOL3018 | 1152078 | 16181910 |
| | 20FR | 16-SLFOL3020 | 1152079 | 16201910 |
| | 22FR | 16-SLFOL3022 | 1152080 | 16221910 |
| | 24FR | 16-SLFOL3024 | 1152081 | 16241910 |

Catheter Accessories

Your patients who use catheters may also need accessories such as surgical lubricants, holders and bags. McKesson Brands has a portfolio of accessories to help you meet your patient needs.



| DESCRIPTION | MFR # | E1 # | ITEM # (BX) | ITEM # (CS) | HCPCS | UOM |
|--|-----------------|---------|-------------|-------------|-------|------------------------|
| McKesson Foley Catheter Holders | MCKFOL | 1107740 | 77401910 | 77401900 | A4334 | 1/EA, 10EA/BX, 10BX/CS |
| McKesson Catheter Securement Device | 172-5445-9 | 1063193 | 44591908 | n/a | A4334 | 1/EA, 50EA/BX |
| McKesson Premium Catheter Securement Device | 18-3400LFC | 1146461 | 18341910 | 18341920 | A4334 | 100/BX |
| McKesson Premium Catheter Securement Device | 18-3400LFC-10PK | 1146463 | 18431910 | 18431920 | A4334 | 10/PK |
| McKesson Urinary Leg Bag, 500 mL | 4601 | 854582 | N/A | 15481908 | A4358 | 1/EA, 48EA/CS |
| McKesson Urinary Leg Bag, 750 mL | 4602 | 854583 | N/A | 25481908 | A4358 | 1/EA, 48EA/CS |
| McKesson Urinary Leg Bag, 1000 mL | 4605 | 911680 | N/A | 47221908 | A4358 | 1/EA, 48EA/CS |
| McKesson Urinary Drain Bag with Cover | 16-2880C | 1112820 | N/A | 28801920 | A4358 | 1/EA, 20EA/CS |
| McKesson Brand Lubricating Jelly, 4 oz | 16-8919 | 1066700 | 16191400 | 16191460 | A4332 | 1/EA, 12EA/BX, 6BX/CS |
| McKesson Brand Lubricating Jelly, 3g Packets | 16-8942 | 1066698 | 16421400 | 16421460 | A4332 | 1/EA, 144EA/BX, 6BX/CS |
| McKesson Brand Lubricating Jelly, 5g Packets | 16-8946 | 1066699 | 16461400 | 16461460 | A4332 | 1/EA, 144EA/BX, 6BX/CS |

NEVER COMPROMISE ON CARESM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you needed, when you need them – making it even easier for you to provide the best care for your patients.

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

McKesson Medical-Surgical

9954 Mayland Dr.
Richmond, VA 23233
mms.mckesson.com/mckesson-brand

Primary Care

866.625.2679

Post-Acute and Home Health

800.347.2456

Home Medical Equipment

888.822.8111

McKesson Brands Support

800.777.4908

McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

All trademarks and registered trademarks are the property of their respective owners.

© 2020 McKesson Medical-Surgical Inc.

2019-1098