


McKesson Brands Catheters

Easy to Use. Designed for Comfort. Wide Selection. Affordable Prices.



Why is the urology category important?

Why McKesson?

Every day, you face multiple challenges when it comes to running a business and serving your patients. McKesson Medical-Surgical is here to help. At McKesson, our focus is always on your success and meeting your patients' needs. That's why our expert team can help to make sure you have the right products, right when you need them.



Between
15%
and
25%

of hospitalized patients receive short-term indwelling catheters



14%
growth
(over previous year)



\$600
million

across all segments



An estimated
5%
of long-term care facility residents use catheters

Source: HIDA Home Care reports, CDC, 2017 and 2018; GHX reporting 2019



McKesson supports the health of your business by offering:



More products, more choice

Quality and selection with more national brands from a single source



Distribution you can count on

Reliable and flexible delivery options including patient home delivery



Solutions to enhance your business

Automated ordering solutions and resources to achieve seamless billing and reimbursements



Reduce readmissions

Reduce hospital readmissions with our comprehensive clinical programs and product offerings



99% fill rate

We take pride in ensuring the right products get to you and your patients



Educate clinical staff

Help your staff achieve patient and resident care goals with specialized clinical training and education

An introduction to urinary catheters

A urinary catheter is a hollow, flexible tube that collects urine from the bladder, leading to a drainage bag.



Common reasons for catheter use

- Diagnosis of urinary incontinence or urinary retention
- For a brief time following major surgery
- After diabetes complications
- Spinal injury
- After prostate surgeries
- Conditions like a stroke or multiple sclerosis

Catheter types


- Foley Balloon
- Intermittent
- External

Foley HCPCS Codes	Intermittent HCPCS Codes	External HCPCS Codes
A4311-A4316 – are kits that include an indwelling catheter	A4351 – Intermittent straight tip with or without coating	A4326 – Male external catheter with collection chamber
A4338 – Indwelling catheter with two-way latex with coating (silicone, silicone elastomer, hydrophilic, etc.)	A4352 – Intermittent Coudé tip with or without coating	A4349 – Male external catheter
A4340 – Indwelling catheter, specialty type (Coudé, mushroom, wing, etc.)	A4353 – Intermittent with insertion supplies (this could be a closed system, introducer tip, etc.)	
A4344 – Indwelling catheter, two-way all silicone		
A4346 – Indwelling three-way (continuous irrigation)		

An introduction to urinary catheters (continued)


You can be confident in the products you offer to your patients when you purchase medical supplies from McKesson. We combine the quality of leading manufacturers, with a fast, nationwide distribution network to help you run a smoother, more efficient business – all backed by a 100% satisfaction guarantee.

French Catheter Sizing



Actual Size

Catheter Sizing Continuum



14 French is the most common size.

Catheter sizes measured in French scale

All McKesson Brand catheters feature:



Polished eyes to help minimize urethral trauma



Color coded packaging to indicate French size



Variety of sizes and tips



Not made with DEHP, DINP, or BPA



Individually packaged for **single sterile use**

PVC Uncoated Catheters

McKesson intermittent catheters help minimize urethral trauma and reduce the chance of bringing bacteria into the bladder. These PVC uncoated catheters are ready-to-use and designed to help your patients self-catheterize quickly and easily.

Features

- Staggered eyes for maximum urine flow
- Polished eyes to help minimize urethral trauma
- Flexible funnel end with ridged design for grip and ease of handling
- Color-coded funnel and packaging label for easy identification of French size
- Packaging design includes loop cutout for hanging to maintain sterility
- Individually packaged for single, sterile use
- Not made with DEHP, DINP, BPA or natural rubber latex



DESCRIPTION	SIZE	MFR #	ITEM #	ITEM # (BX)
McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 6"	10FR	16-F610	1020810	31661900
	12FR	16-F612	1020812	31621900
	14FR	16-F614	1020814	36641900
	16FR	16-F616	1020818	56611900
	18 FR	16-F618	1135715	36981900
McKesson Intermittent Catheter, Uncoated PVC, Luer Tip, 6"	14FR	16-F614L	1020816	35641900
McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 10"	14FR	16-U1014	1135717	37331900
McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 16"	10FR	16-M1610	1020809	69611900
	12FR	16-M1612	1020811	62121900
	14FR	16-M1614	1020813	66641900
	16FR	16-M1616	1020817	69661900
	18 FR	16-M1618	1135716	37071900
McKesson Intermittent Catheter, Uncoated PVC, Coude Tip, 16"	10FR	16-M1610C	1135718	75701900
	12FR	16-M1612C	1135719	75711900
	14FR	16-M1614C	1020815	32641900
	16FR	16-M1616C	1020819	23661900
	18 FR	16-M1618C	1135720	37081900

PVC Hydrophilic Catheters

McKesson hydrophilic catheters are pre-lubricated and hydrated to help promote excellent gliding properties. With easier, more comfortable insertion, hydrophilic catheters offer decreased friction. This leads to better patient satisfaction by helping reduce the risk of urethral damage.

Features

- Staggered eyes for maximum urine flow
- Polished eyes help minimize urethral trauma
- Flexible funnel end with ridged design for grip and ease of handling
- Color-coded funnel and packaging label for easy identification of French size
- Packaging design includes loop cutout for hanging to maintain sterility
- Individually packaged for single, sterile use
- Not made with DEHP, DINP, BPA or natural rubber latex



DESCRIPTION	SIZE	MFR #	ITEM #	ITEM # (BX)
McKesson Intermittent Catheter, Hydrophilic, Straight Tip, 6"	12FR	16-FH612	1135721	36991900
	14FR	16-FH614	1135722	75691900
McKesson Intermittent Catheter, Uncoated PVC, Luer Tip, 6"	10FR	16-FH810	1135723	37021900
	14FR	16-FH814	1135724	56391900
McKesson Intermittent Catheter, Uncoated PVC, Straight Tip, 10"	10FR	16-MH1610	1135725	37091900
	12FR	16-MH1612	1135726	37131900
	14FR	16-MH1614	1135727	37171900
	16FR	16-MH1616	1135728	37251900
	18 FR	16-MH1618	1135729	37291900
	10FR	16-M1610	1020809	69611900
McKesson Intermittent Catheter, Hydrophilic, Coude Tip, 16"	12FR	16-M1612	1020811	62121900
	14FR	16-M1614	1020813	66641900
	16FR	16-M1616	1020817	69661900
	18 FR	16-M1618	1135716	37071900

Silicone Foley Catheters

A Foley is a catheter with two channels and one tube. The first channel drains urine into a collection bag, while the other channel has a valve on the outside end and connects to a balloon at the inside tip. The balloon is inflated with sterile water when it lies inside the bladder to prevent the catheter from slipping out.

A Foley catheter is meant to be in place for no longer than three months. Traditional straight catheters are inserted once and discarded when the bladder has been emptied.

McKesson 2-Way 100% Silicone Foley Catheters are ideal for patients with latex allergies and longer term catheterization. Silicone provides chemical and thermal stability, low surface tension and hydrophobicity.

Features

- Eliminates latex allergy concerns
- Smooth eyes improve patient comfort
- Reinforced tip provides a stylet for optimized insertion efficacy
- 100% silicone material has non-adherent characteristics and reduces encrustation and calcifications
- Symmetrical balloon provides a tight seal at the neck of the urethra



DESCRIPTION	SIZE	MFR #	ITEM #	ITEM # (BX)
McKesson Foley Catheter, Silicone, Straight Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume	12FR	16-SFOL512	1152043	65121910
	14FR	16-SFOL514	1152044	65141910
	16FR	16-SFOL516	1152045	65161910
	18FR	16-SFOL518	1152046	65181910
	20FR	16-SFOL520	1152047	65201910
	22FR	16-SFOL522	1152048	65221910
	24FR	16-SFOL524	1152049	65241920
McKesson Foley Catheter, Silicone, Straight Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume	16FR	16-SFOL3016	1152038	63161910
	18FR	16-SFOL3018	1152039	63181910
	20FR	16-SFOL3020	1152040	63201910
	22FR	16-SFOL3022	1152041	63221910
	24FR	16-SFOL3024	1152042	63241910
McKesson Foley Catheter, Silicone, Coude Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume	12FR	16-SFOL512C	1152031	15121910
	14FR	16-SFOL514C	1152032	15141910
	16FR	16-SFOL516C	1152033	15161910
	18FR	16-SFOL518C	1152034	15181910
	20FR	16-SFOL520C	1152035	15201910
	22FR	16-SFOL522C	1152036	15221910

Hydrophilic-Coated Latex Foley Catheters

McKesson hydrophilic-coated latex Foley catheters are retention type urological catheters made of silicone-coated natural latex material with the addition of a lubricious/hydrophilic coating to facilitate the insertion of the catheter.

Features

- Hydrophilic coating helps urethral insertion and helps reduce irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Contains natural rubber latex



DESCRIPTION	SIZE	MFR #	ITEM #	ITEM # (BX)
McKesson Foley Catheter, Hydrophilic Coated Latex, Straight Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume	12FR	16-LFOL512	1152068	51211910
	14FR	16-LFOL514	1152069	51411910
	16FR	16-LFOL516	1152070	16151910
	18FR	16-LFOL518	1152071	18151910
	20FR	16-LFOL520	1152072	52011910
	22FR	16-LFOL522	1152073	52211910
	24FR	16-LFOL524	1152074	52411910
McKesson Foley Catheter, Hydrophilic Coated Latex, Straight Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume	14FR	16-LFOL3014	1152062	30141910
	16FR	16-LFOL3016	1152063	31661910
	18FR	16-LFOL3018	1152064	31861910
	20FR	16-LFOL3020	1152065	32061910
	22FR	16-LFOL3022	1152066	32261910
	24FR	16-LFOL3024	1152067	32461910
McKesson Foley Catheter, Hydrophilic Coated Latex, Coude Tip, 16", 5cc Balloon Size, 5-15cc Inflation Volume	12FR	16-LFOL512C	1152055	15621910
	14FR	16-LFOL514C	1152056	51461910
	16FR	16-LFOL516C	1152057	51661910
	18FR	16-LFOL518C	1152058	51861920
	20FR	16-LFOL520C	1152059	52061910
	22FR	16-LFOL522C	1152060	52261920
	24FR	16-LFOL524C	1152061	52461910
McKesson Foley Catheter, Hydrophilic Coated Latex, Coude Tip, 16", 30cc Balloon Size, 30-40cc Inflation Volume	16FR	16-LFOL3016C	1152050	30161920
	18FR	16-LFOL3018C	1152051	30181910
	20FR	16-LFOL3020C	1152052	30201920
	22FR	16-LFOL3022C	1152053	22031920
	24FR	16-LFOL3024C	1152054	30241920

Silicone Coated Latex Foley Catheters

McKesson two-way siliconized latex Foley catheters are durable and flexible to help improve patient comfort. Silicone is typically resistant to chemical attack and insensitive to temperature changes.

Features

- Silicone coating helps reduce urethral irritation
- Smooth beveled symmetrical eyes minimize tissue damage
- Large oval eyes aid in maximum drainage
- Firm catheter tip designed to improve ease of insertion
- Symmetrical balloon provides a tight seal at the neck of the urethra
- Contains natural rubber latex



DESCRIPTION	SIZE	MFR #	ITEM #	ITEM # (BX)
McKesson Foley Catheter, Silicone Coated Latex, Straight Tip, 5cc Balloon Size, 5-15cc Inflation Volume	12FR	16-SLFOL512	1152082	16521910
	14FR	16-SLFOL514	1152083	16541910
	16FR	16-SLFOL516	1152084	65611920
	18FR	16-SLFOL518	1152085	85611920
	20FR	16-SLFOL520	1152086	25611910
	22FR	16-SLFOL522	1152087	22561910
	24FR	16-SLFOL524	1152088	42561910
McKesson Foley Catheter, Silicone Coated Latex, Straight Tip, 30cc Balloon Size, 30-40cc Inflation Volume	12FR	16-SLFOL3012	1152075	16121910
	14FR	16-SLFOL3014	1152076	16141910
	16FR	16-SLFOL3016	1152077	16611910
	18FR	16-SLFOL3018	1152078	16181910
	20FR	16-SLFOL3020	1152079	16201910
	22FR	16-SLFOL3022	1152080	16221910
	24FR	16-SLFOL3024	1152081	16241910

Catheter Accessories

Your patients who use catheters may also need accessories such as surgical lubricants, holders and bags. McKesson Brands has a portfolio of accessories to help you meet your patient needs.



DESCRIPTION	MFR #	E1 #	ITEM # (BX)	ITEM # (CS)	HCPCS	UOM
McKesson Foley Catheter Holders	MCKFOL	1107740	77401910	77401900	A4334	1/EA, 10EA/BX, 10BX/CS
McKesson Catheter Securement Device	172-5445-9	1063193	44591908	n/a	A4334	1/EA, 50EA/BX
McKesson Premium Catheter Securement Device	18-3400LFC	1146461	18341910	18341920	A4334	100/BX
McKesson Premium Catheter Securement Device	18-3400LFC-10PK	1146463	18431910	18431920	A4334	10/PK
McKesson Urinary Leg Bag, 500 mL	4601	854582	N/A	15481908	A4358	1/EA, 48EA/CS
McKesson Urinary Leg Bag, 750 mL	4602	854583	N/A	25481908	A4358	1/EA, 48EA/CS
McKesson Urinary Leg Bag, 1000 mL	4605	911680	N/A	47221908	A4358	1/EA, 48EA/CS
McKesson Urinary Drain Bag with Cover	16-2880C	1112820	N/A	28801920	A4358	1/EA, 20EA/CS
McKesson Brand Lubricating Jelly, 4 oz	16-8919	1066700	16191400	16191460	A4332	1/EA, 12EA/BX, 6BX/CS
McKesson Brand Lubricating Jelly, 3g Packets	16-8942	1066698	16421400	16421460	A4332	1/EA, 144EA/BX, 6BX/CS
McKesson Brand Lubricating Jelly, 5g Packets	16-8946	1066699	16461400	16461460	A4332	1/EA, 144EA/BX, 6BX/CS

NEVER COMPROMISE ON CARESM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you needed, when you need them – making it even easier for you to provide the best care for your patients.

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

McKesson Medical-Surgical

9954 Mayland Dr.
Richmond, VA 23233
mms.mckesson.com/mckesson-brand

Primary Care

866.625.2679

Post-Acute and Home Health

800.347.2456

Home Medical Equipment

888.822.8111

McKesson Brands Support

800.777.4908

McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

All trademarks and registered trademarks are the property of their respective owners.

© 2020 McKesson Medical-Surgical Inc.

2019-1098