

Safety Hypodermic Needles and Syringes with Safety Hypodermic Needles

Prevent® safety hypodermic needles are designed with a needle locking mechanism that helps to significantly reduce accidental needle sticks. Once used, the needle can be clipped into a sheath and locked into place. An audible click indicates device activation. In addition to the broad range of needle sizes, McKesson also offers a selection of Prevent syringes with safety hypodermic needles that include 1cc and 3cc syringes. Choose Prevent needles and syringes for safety, ease of use and standardization.

FEATURES AND BENEFITS

- Dual, reinforced needle locking mechanism helps reduce risk of needlestick injury
- Needle point style: beveled tip
- Simple, intuitive, no complicated mechanisms
- ISO color-coded packaging and needle hubs enable easy identification of gauge size
- Selection includes thin wall, ultra thin wall and regular wall needles
- Hub and needle both lock with audible and tactile feedback
- Needles are compatible with conventional luer slip and luer lock syringes
- Syringe material: plastic
- Needle material: stainless steel
- Sterile
- Non-pyrogenic
- Single use only
- Not made with natural rubber latex



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.



Safety Hypodermic Needles and Syringes with Safety Hypodermic Needles

MFR#	DESCRIPTION	QTY
102-SN3C22105S3	Syringes with Safety Hypodermic Needles, 3cc, 22G x 1 1/2 in	100 ea/bx, 4 bx/cs
102-SN3C221S3	Syringes with Safety Hypodermic Needles, 3cc, 22G x 1 in	100 ea/bx, 4 bx/cs
102-SN3C231S3	Syringes with Safety Hypodermic Needles, 3cc, 23G x 1 in	100 ea/bx, 4 bx/cs
102-SN3C251S3	Syringes with Safety Hypodermic Needles, 3cc, 25G x 1 in	100 ea/bx, 4 bx/cs
102-SN3C2558S3	Syringes with Safety Hypodermic Needles, 3cc, 25G x 5/8 in	100 ea/bx, 4 bx/cs
102-SNT1C2558S3	Syringes with Safety Hypodermic Needles, 1cc, 25G x 5/8 in	100 ea/bx, 4 bx/cs
102-SNT1C2705S3	Syringes with Safety Hypodermic Needles, 1cc, 27G x 1/2 in	100 ea/bx, 4 bx/cs
102-N18105S3	Safety Hypodermic Needles, 18G x 1 1/2 in	100 ea/bx, 8 bx/cs
102-N181S3	Safety Hypodermic Needles, 18G x 1 in	100 ea/bx, 8 bx/cs
102-N201S3	Safety Hypodermic Needles, 20G x 1 in	100 ea/bx, 8 bx/cs
102-N21105S3	Safety Hypodermic Needles, 21G x 1 1/2 in	100 ea/bx, 8 bx/cs
102-N211S3	Safety Hypodermic Needles, 21G x 1 in	100 ea/bx, 8 bx/cs
102-N22105S3	Safety Hypodermic Needles, 22G x 1 1/2 in	100 ea/bx, 8 bx/cs
102-N221S3	Safety Hypodermic Needles, 22G x 1 in	100 ea/bx, 8 bx/cs
102-N23105S3	Safety Hypodermic Needles, 23G x 1 1/2 in	100 ea/bx, 8 bx/cs
102-N231S3	Safety Hypodermic Needles, 23G x 1 in	100 ea/bx, 8 bx/cs
102-N25105S3	Safety Hypodermic Needles, 25G x 1 1/2 in	100 ea/bx, 8 bx/cs
102-N251S3	Safety Hypodermic Needles, 25G x 1 in	100 ea/bx, 8 bx/cs
102-N2558S3	Safety Hypodermic Needles, 25G x 5/8 in	100 ea/bx, 8 bx/cs
102-N2705S3	Safety Hypodermic Needles, 27G x 1/2 in	100 ea/bx, 8 bx/cs
102-N3005S3	Safety Hypodermic Needles, 30G x 1/2 in	100 ea/bx, 8 bx/cs