Amsino is proud to offer Infrared Forehead Thermometer



FAST, ACCURATE & COMPACT

Measures temperature in 1 second. Use for forehead reading. Large-screen LCD display, backlit, can be clearly displayed under any light

NON-CONTACT

The infrared technology reads temperature without any physical touch

FEVER FUNCTION

Alerts with beeps, vibration, and yellow color display in case of high body temperature; used in hospitals, home, and health clinics

USEFUL FOR ALL

Safe for use among newborns, infants, kids, and adults alike. It is a perfect choice for Moms and Dads who want to keep a close eye on the body temperature of babies, children and loved ones without disturbing their sleep

EASY & CONVENIENT

User-friendly design with LCD display; can switch between Fahrenheit and Celsius temperature units

Product Details

- · Measurement place: Forehead
- Display: LCD screen with 4-digits and special icons
- Display resolution: 0.1 °F (0.1 °C)
- Human body temperature mode: 89.6 109.2 °F (32.0 42.9 °C)
- Accuracy: ±0.4 °F (±0.2 °C) for the range 93.2 to 109.4 °F (34 to 43 °C)
- Small size: 5.9" x 1.5" x 1.7"
- Battery requirements: 2 x 1.5V AAA batteries
- Case packaging: 22.6"x 13.8" x 6.1"
- Case Net Weight: 5.9 lbs Case Gross Weight: 13.7 lbs
- Quantities/Unit: 40 pcs/case
- Warranty: One year manufacturer warranty

Item Code	Description	Quantity/Unit
YT-1	Infrared Forehead Thermometer	40/Case

This device is manufactured consistent with FDA 21 CFR Part 820, ISO 13485:2016 Quality Management System and is in compliance with ASTM E1965-98, ASTM E1112-00, ASTM E825-98, ANSI/AAMI ES60601-1:2005(R2012), ANSI/AAMI/IEC 60601-1-2:2014, ANSI/AAMI/IEC 60601-1-11:2015, ANSI/AAMI/IEC -62304:2006/A1:2016, ANSI/AAMI/ISO 10993-1, ISO 80601-2-56 standards. The product is CE marked. This device is not 510(k) cleared. This device can be legally distributed and sold in the US during COVID 19 in accordance with FDA's "Enforcement Policy for Clinical Electronic Thermometers During the Coronavirus Disease 2019 (COVID 19) Public Health Emergency" April 2020.

Call Amsino sales reps for pricing, lead time and any other order information. To view our complete product line, visit us at **www.amsino.com**



Amsino International, Inc.

708 Corporate Center Drive, Pomona, CA 91768 800.632.6746 • fax: 909.626.3888 • info@amsino.com Amsino management cooperates with Yuwell to bring new products to market to improve patient care and safety. Amsino is Yuwell's exclusive distributor in USA.

Product(s) described may not be licensed or available for sale in all countries. Please see the Instructions for Use/Operator's Manual for a complete listing of the indications, contraindications, warnings and precautions. The symbol @ indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. © Copyright 2020, Amsino International, Inc. All rights reserved.