

Centrifuge Tubes

McKesson centrifuge tubes are high-quality, translucent tubes molded from polypropylene (PP) material featuring a conical bottom and flat-top, leak-resistant blue caps. The tubes are printed with black graduation marks and have a large white writing area.

Features and Benefits

- Black printed graduations in 0.5 mL increments for easy measurements
- Large white printed writing space for convenient labeling
- Autoclavable tubes for continued use (without cap)
- Included screw caps are made of high density polyethylene (HDPE), flat-top and leak-resistant
- Not made with natural rubber latex



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

CENTRIFUGE TUBES

MFR #	DESCRIPTION	ATTRIBUTES	SIZE	QTY
177-6284	Polypropylene (PP), Non-Sterile	Blue Screw Cap, Graduations, White Writing Area	15 mL	500/cs
177-6285	Polypropylene (PP), Sterile	Blue Screw Cap, Graduations, White Writing Area	15 mL	25/bg, 20 bg/cs
177-6286	Polypropylene (PP), Racked, Sterile	Blue Screw Cap, Graduations, White Writing Area	15 mL	50/rack, 10 racks/cs
177-6288	Polypropylene (PP), Sterile	Blue Screw Cap, Graduations, White Writing Area	50 mL	25/bg, 20 bg/cs
177-6289	Polypropylene (PP), Racked, Sterile	Blue Screw Cap, Graduations, White Writing Area	50 mL	25/rack, 20 racks/cs
177-6295PL	Polypropylene (PP), Sterile	Red Screw Cap, Graduations, White Writing Area	15 mL	25/rack, 25 racks/cs
177-6296PL	Polypropylene (PP), Sterile	Red Screw Cap, Graduations, White Writing Area	15 mL	25/rack, 25 racks/cs
177-6297PL	Polypropylene (PP), Sterile	Red Screw Cap, Graduations, White Writing Area	50 mL	25/rack, 25 racks/cs
177-6298PL	Polypropylene (PP), Sterile	Red Screw Cap, Graduations, White Writing Area	50 mL	25/rack, 25 racks/cs

McKesson Medical-Surgical
 9954 Mayland Drive
 Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
 800.347.2456

Home Medical Equipment
 888.822.8111

Primary Care
 866.625.2679

McKesson Brands Support
 800.777.4908