

McKESSON



Shielded IV Catheter

McKesson Prevent® R shielded IV catheters are designed to be inserted into the patient's vascular system for short term use to sample blood, monitor blood pressure, or administer fluids intravenously. It has a push-button safety design to allow one-handed needle retraction into the safety shield for reduced risk of accidental needle-stick injuries.

FEATURES AND BENEFITS

- Push-button safety design allows one-handed needle retraction into the safety shield for reduced risk of accidental needle-stick injuries
- ISO-compliant, color-coded catheter hubs and packaging for easy identification of gauge size
- Flash chamber provides visual confirmation that catheter has entered the vein
- Soft, flexible, kink-resistant material further softens inside to conform to the shape of the vessel
- Radiopaque materials allow use during x-rays or other forms of radiation
- Non-pyrogenic and not made with natural rubber latex or DEHP



Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2018 McKesson Medical-Surgical Inc.

McKESSON



MFR#	McKB MFR#	DESCRIPTION	HUB TYPE	NOTCHED NEEDLE Y/N	SIZE	QTY
380267	192-PB14175	Shielded IV Catheter	Straight Hub	No	14G x 1.75 in	50 ea/bx, 4 bx/cs
380254	192-PB16116	Shielded IV Catheter	Straight Hub	No	16G x 1.16 in	50 ea/bx, 4 bx/cs
380244	192-PB18116	Shielded IV Catheter	Straight Hub	No	18G x 1.16 in	50 ea/bx, 4 bx/cs
380233	192-PB201	Shielded IV Catheter	Straight Hub	Yes	20G x 1 in	50 ea/bx, 4 bx/cs
380234	192-PB20116	Shielded IV Catheter	Straight Hub	Yes	20G x 1.16 in	50 ea/bx, 4 bx/cs
380223	192-PB221	Shielded IV Catheter	Straight Hub	Yes	22G x 1 in	50 ea/bx, 4 bx/cs
380211	192-PB245	Shielded IV Catheter	Straight Hub	Yes	24G x 0.5 in	50 ea/bx, 4 bx/cs
381412	192-PB2475	Shielded IV Catheter	Straight Hub	Yes	24G x 0.75 in	50 ea/bx, 4 bx/cs
380333	192-PB201W	Shielded IV Catheter	Winged Hub	Yes	20G x 1 in	50 ea/bx, 4 bx/cs
380334	192-PB20116W	Shielded IV Catheter	Winged Hub	Yes	20G x 1.16 in	50 ea/bx, 4 bx/cs
380323	192-PB221W	Shielded IV Catheter	Winged Hub	Yes	22G x 1 in	50 ea/bx, 4 bx/cs
380312	192-PB2475W	Shielded IV Catheter	Winged Hub	Yes	24G x 0.75 in	50 ea/bx, 4 bx/cs

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908