

McKESSON

PREVENT® G

Safety Winged Blood Collection Sets **GLIDE**

McKesson's Prevent G safety winged blood collection sets feature a single-use, sterile, winged needle bonded to flexible tubing with a Luer adapter. The sets are sterile and non-pyrogenic. Each blood collection set includes an easy-to-use integrated safety device that reduces accidental needle sticks. Since you don't have to change your needle insertion technique, it is easy to make the switch to Prevent.

FEATURES AND BENEFITS

- One-hand technique activates safety mechanism
- Compact design enhances control and flexibility
- Wing design permits improved positioning during utilization
- Flexible tubing assists in handling
- Built-in safety device
- Audible click confirms that needle is securely locked
- MFR #s 16-7BC2134, 16-7BC2334, and 16-BC2334H contain a 7-inch tube; MFR #s 16-BC2134, 16-BC2334, 16-BC2534, and 16-BC2134H have a 12-inch tube
- MFR #s 16-BC2134H and 16-BC2334H include a pre-attached tube holder
- MFR #s 16-BC2334NL comes without luer adapter
- Not made with DEHP
- Not made with natural rubber latex



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2021 McKesson Medical-Surgical Inc.

McKESSON



Safety Winged Blood Collection Sets **GLIDE**

MFR#	DESCRIPTION	QTY
16-BC2134	Safety Winged Blood Collection Set Glide, 21G x 3/4" x 12"	50 ea/bx, 4 bx/cs
16-BC2334	Safety Winged Blood Collection Set Glide, 23G x 3/4" x 12"	50 ea/bx, 4 bx/cs
16-BC2534	Safety Winged Blood Collection Set Glide, 25G x 3/4" x 12"	50 ea/bx, 4 bx/cs
16-7BC2134	Safety Winged Blood Collection Set Glide, 21G x 3/4" x 7"	50 ea/bx, 4 bx/cs
16-7BC2334	Safety Winged Blood Collection Set Glide, 23G x 3/4" x 7"	50 ea/bx, 4 bx/cs
16-BC2134H	Safety Winged Blood Collection Set Glide, with Tube Holder, 21G x 3/4" x 12"	25 ea/bx, 8 bx/cs
16-BC2334H	Safety Winged Blood Collection Set Glide, with Tube Holder, 23G x 3/4" x 7"	25 ea/bx, 8 bx/cs
16-BC2334NL	Safety Winged Blood Collection Set Glide, without Luer Adapter, 23G x 3/4" x 12"	50 ea/bx, 4 bx/cs



MFR# 16-BC2134



MFR# 16-BC2334



MFR# 16-BC2534