

IV Administration Sets

McKesson IV administration sets connect a IV solution bag to a patient that contain ports (needle and needle-free) that allow the the infusion of medication into the solution. The sets are designed with slide and roller clamps to control the flow rate.

Features and benefits

- Gravity sets are used to connect an IV solution bag to a patient
- The needle-free connectors (Y ports) allow the infusion of medication into the solution (MFR MS9212NF has 2 needle-free connectors)
- Male luer lock connection
- Easy-to-use slide and roller clamps help control the flow rate (slide clamps with MFR MS9212NF and MS921NF)
- Provides a drip rate of 20 drops/mL
- Air-vented spike with cap
- Includes standard bore tubing, see chart for lengths
- Storage conditions: do not store above 104°F (40°C)
- Not made with natural rubber latex or DEHP



ORDERING INFORMATION

MFR #	DESCRIPTION	ATTRIBUTES	TUBING LENGTH	QTY
MS9212NF	IV Administration Set, Standard Bore	Two Needle-Free Connectors, Air Vented, with Slide Clamp, Sterile	80 inch	1/pk, 50 pks/bx, 6 bx/cs
MS921NF	IV Administration Set, Standard Bore	Needle-Free Connector, Air Vented, with Slide Clamp, Sterile	75 inch	1/pk, 50 pks/bx, 6 bx/cs
MS919PNF	IV Administration Set, Standard Bore	Needle-Free Connector, Air Vented, Sterile	97 inch	1/pk, 50 pks/bx, 6 bx/cs

WARNING: This product can expose you to Bisphenol A (BPA), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

Satisfaction Guaranteed

For complete details, please visit mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.