


Dermal Curettes

McKesson dermal curettes are designed for scraping skin lesions and growths such as warts, basal cell carcinomas, and melanomas. The loop is sharp on the inside and blunt on the outside so the lesion can be removed while the surrounding skin is not damaged. A variety of tip sizes are available to accommodate varying lesion sizes.

Features and Benefits

- Stainless steel blade
- Textured handle provides better comfort and control
- Curette sizes: 2mm - 7mm
- Length: 5.5 inches
- Individually packaged
- Sterile, disposable
- Not made with natural rubber latex



ORDERING INFORMATION

MFR #	DESCRIPTION	SIZE	QTY
16-DC2	Dermal Curette, Disposable, Stainless Steel, Sterile	2mm	1 ea, 50 ea/bx, 4 bx/cs
16-DC3	Dermal Curette, Disposable, Stainless Steel, Sterile	3mm	1 ea, 50 ea/bx, 4 bx/cs
16-DC4	Dermal Curette, Disposable, Stainless Steel, Sterile	4mm	1 ea, 50 ea/bx, 4 bx/cs
16-DC5	Dermal Curette, Disposable, Stainless Steel, Sterile	5mm	1 ea, 50 ea/bx, 4 bx/cs
16-DC7	Dermal Curette, Disposable, Stainless Steel, Sterile	7mm	1 ea, 50 ea/bx, 4 bx/cs

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.