

Rollators

McKesson rollators are ideal when arm strength is limited. Rollators, unlike standard walkers, do not require users to lift the unit with each step and they also include brakes to maintain control. With convenient storage and seating, our rollators allow users to transport personal items without requiring the use of hands and provide the ability for users to sit or rest. Choose McKesson's rollators as a mobility solution for the enhanced maneuverability across various terrain that will keep users moving smoothly and safely.

Features and benefits

Lightweight rollators

- Aluminum frame construction safely supports users up to 300 lbs
- 13-pound, lightweight design with folding capabilities for easy storage and portability
- Ergonomically-designed handles are easy-to-grip, reduce hand pressure and are adjustable to accommodate users' height and preferences
- Oversized padded seat and backrest are made of strong polyurethane material for durability, have a contoured design for ergonomic comfort, and provide cut and split resistance to prevent rubbing and irritation
- Hinged padded backrest can be folded up, down or removed
- Deluxe loop-lock system operates like a bicycle hand brake to slow down or temporarily stop and easily locks into stationary position by pushing handle into down-position
- Swiveling, 6-inch caster wheels with soft-grip tires are non-marking and add stability and control for indoor and outdoor use
- Serrated-edged brakes provide firm hold
- Includes convenient zippered storage pouch on the bottom of the seat and spacious vinyl carry pouch below for transporting personal items with added privacy and security
- 14-inch seat width; 20-inch seat-to-floor height
- Not made with natural rubber latex
- HCPCS: E0143, E0156



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Rollators

Features and benefits

Durable rollators

- Steel frame construction safely supports up to 300 lbs
- Designed with folding hinge for easy storage and portability
- Ergonomically-designed handles are easy-to-grip, reduce hand pressure and are adjustable to accommodate users' height and preferences
- Lever allows for independent handle height and angle adjustment
- Oversized padded seat and backrest are made of strong polyurethane material for durability, have a contoured design for ergonomic comfort, and provide cut and split resistance to prevent rubbing and irritation
- Hinged padded backrest can be folded up, down or removed
- Deluxe loop-lock system operates like a bicycle hand brake to slow down or temporarily stop and easily locks into stationary position by pushing handle into down-position
- Swiveling, 7 1/2-inch, non-marking casters are ideal for indoor and outdoor use
- Includes convenient cloth storage pouch for personal items
- 12-inch seat width; 21-inch seat-to-floor height
- Not made with natural rubber latex
- HCPCS: E0143, E0156



Rollators

Features and benefits

Heavy-Duty rollators

- Exceptionally strong steel frame construction safely supports bariatric users up to 400 lbs
- Designed with folding hinge for easy storage and transport
- Ergonomically-designed handles are easy-to-grip, reduce hand pressure and are adjustable to accommodate users' height and preferences
- Self-threading knob allows easy, tool-free handle height adjustment
- Oversized padded seat and backrest are made of strong polyurethane material for durability, have a contoured design for ergonomic comfort, and provide cut and split resistance to prevent rubbing and irritation
- Hinged padded backrest can be folded up, down or removed
- Deluxe loop-lock system, made of internal aluminum casting, operates like a bicycle hand brake to slow down or temporarily stop and easily locks into stationary position by pushing handle into down-position
- Swiveling, 8-inch caster wheels with soft-grip tires are non-marking and add stability and control for indoor and outdoor use
- Seat folds up to conveniently access the spacious vinyl carry pouch below the seat for transporting personal items with added privacy and security
- 19-inch seat width; 22-inch seat-to-floor height
- Not made with natural rubber latex
- HCPCS: E0149, E0156



ORDERING INFORMATION

ROLLATORS					
MFR #	DESCRIPTION	FRAME	COLOR	WEIGHT CAPACITY	QTY
146-R726RD	Lightweight	Aluminum	Red	300 lbs	1 ea
146-R726BL	Lightweight	Aluminum	Blue	300 lbs	1 ea
146-R726GR	Lightweight	Aluminum	Green	300 lbs	1 ea
146-R726BK	Lightweight	Aluminum	Black	300 lbs	1 ea
146-R726SL	Lightweight	Aluminum	Silver	300 lbs	1 ea
146-RTL10261RD	Lightweight	Aluminum	Red	300 lbs	1 ea
146-R728RD	Lightweight	Aluminum	Red	300 lbs	1 ea
146-R728BL	Lightweight	Aluminum	Blue	300 lbs	1 ea
146-10257RD-1	Durable	Steel	Red	300 lbs	1 ea
146-10257BL-1	Durable	Steel	Blue	300 lbs	1 ea
146-10216RD-1	Heavy-Duty	Steel	Red	400 lbs	1 ea
146-10216BL-1	Heavy-Duty	Steel	Blue	400 lbs	1 ea
146-10215RD-1	Bariatric	Steel	Red	500 lbs	1 ea
146-10215BL-1	Bariatric	Steel	Blue	500 lbs	1 ea