

Wheelchairs

McKesson wheelchairs are ideal for patients who have difficulty walking due to illness, injury or disability. McKesson offers a comprehensive selection of wheelchairs that promote patient independence, improve mobility and enhance quality of life. With durable construction, easy-to-adjust footrests and lightweight, easy-to-clean upholstery, McKesson wheelchairs provide long-lasting performance and comfort. Dual axle options provide easy transitions to hemi-level seat height to allow users to propel the wheelchair with their feet. Choose a McKesson wheelchair based on arm style, footrest preferences, reclining options, size and weight capacity.

Features and Benefits

Standard Wheelchairs

- Carbon steel frame with a textured, silver-vein finish for enhanced visual appeal and ease of maintenance
- Embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Urethane tires are mounted on composite wheels to provide durability and a smooth ride over most surfaces
- 8-inch front casters are adjustable in three positions
- Convenient carry pouch for transporting personal items
- Not made with natural rubber latex



Lightweight Wheelchairs

- Carbon steel frame with a textured, silver-vein finish for enhanced visual appeal and ease of maintenance
- Nylon upholstery is lightweight, and breathable for user's lasting comfort
- Urethane tires are mounted on mag-style, composite wheels to provide durability and a smooth ride over most surfaces
- Precision-sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- 8-inch front caster wheels provide a smooth ride over most surfaces
- Flipback attachable armrests for patient comfort and easy transfer
- Not made with natural rubber latex



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Wheelchairs

Features and Benefits

Heavy-Duty Wheelchairs

- Increased height, width and durable steel construction yields heavy-duty chair ideal for bariatric use
- 450 lb weight capacity
- Carbon steel frame with a triple-coated chrome finish provides for an attractive, chip-proof appearance
- Double-embossed vinyl upholstery is durable, lightweight, attractive and easy-to-clean
- Urethane tires are mounted on composite wheels to provide durability and a smooth ride over most surfaces
- Precision-sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Push-to-lock wheel locks
- Not made with natural rubber latex



Reclining Wheelchairs

- Hydraulic reclining mechanism allows for adjustments up to 180°
- Carbon steel frame with a textured, silver-vein finish for enhanced visual appeal and ease of maintenance
- Nylon upholstery is lightweight, and breathable for user's lasting comfort
- Urethane tires are mounted on composite wheels to provide durability and a smooth ride over most surfaces
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Single axle design
- Designed with wheels set back on frame and rear anti-tippers to prevent toppling
- Includes headrest extension with cushioned head immobilizer
- Not made with natural rubber latex



ORDERING INFORMATION

STANDARD, DUAL AXLE WHEELCHAIRS							
MFR #		ARM STYLE	FOOTREST OR LEGREST TYPE	SEAT WIDTH	WEIGHT CAPACITY	HCPCS	QTY
146-SSP216DDA-SF	⚠	Detachable Desk Arms	Swing-Away Footrests	16 in	250 lbs.	K0001 or K0002	1 ea
146-SSP216DDA-ELR	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests	16 in	250 lbs.	K0001 or K0002	1 ea
146-SSP218DDA-ELR	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests	18 in	300 lbs.	K0001 or K0002	1 ea
146-SSP218FA-SF	⚠	Fixed Arms	Swing-Away Footrests	18 in	300 lbs.	K0001 or K0002	1 ea
146-SSP218FA-ELR	⚠	Fixed Arms	Swing-Away, Elevating Footrests	18 in	300 lbs.	K0001 or K0002	1 ea
146-SSP218DDA-SF	⚠	Detachable Desk Arms	Swing-Away Footrests	18 in	300 lbs.	K0001 or K0002	1 ea
146-SSP220DDA-SF	⚠	Detachable Desk Arms	Swing-Away Footrests	20 in	350 lbs.	K0001 or K0002+E2201	1 ea
146-SSP220DDA-ELR	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests	20 in	350 lbs.	K0001 or K0002+E2201	1 ea
LIGHTWEIGHT, DUAL AXLE WHEELCHAIRS							
146-K316DDA-SF	⚠	Flip-Back Detachable Desk Arm	Swing-Away Footrests	16 in	300 lbs.	K0003	1 ea
146-K316DDA-ELR	⚠	Flip-Back Detachable Desk Arm	Swing-Away, Elevating Footrests	16 in	300 lbs.	K0003	1 ea
146-K318DDA-SF	⚠	Flip-Back Detachable Desk Arm	Swing-Away Footrests	18 in	300 lbs.	K0003	1 ea
146-K318DDA-ELR	⚠	Flip-Back Detachable Desk Arm	Swing-Away, Elevating Footrests	18 in	300 lbs.	K0003	1 ea
146-K320DDA-SF	⚠	Flip-Back Detachable Desk Arm	Swing-Away Footrests	20 in	300 lbs.	K0003+E2201	1ea
146-K320DDA-ELR	⚠	Flip-Back Detachable Desk Arm	Swing-Away, Elevating Footrests	20 in	300 lbs.	K0003+E2201	1 ea
HEAVY DUTY, DUAL AXLE WHEELCHAIRS							
146-STD22ECDDA-SF	⚠	Detachable Desk Arms	Swing Away Footrests	22 in	450 lbs.	K0007	1 ea
146-STD22ECDDA-ELR	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests	22 in	450 lbs.	K0007+K0195 or E0990	1 ea
146-STD22ECDFA-SF	⚠	Detachable Full Arms	Swing Away Footrests	22 in	450 lbs.	K0007	1ea
146-STD24ECDDA-SF	⚠	Detachable Desk Arms	Swing Away Footrests	24 in	450 lbs.	K0007	1ea
146-STD24ECDDA-ELR	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests	24 in	450 lbs.	K0007+K0195 or E0990	1 ea
RECLINING, SINGLE AXLE WHEELCHAIRS							
146-SSP18RBDDA	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests with Metal Foot Plates	18 in	300 lbs.	K0001+E1226 +K0195 or E0990	1 ea
146-SSP20RBDDA	⚠	Detachable Desk Arms	Swing-Away, Elevating Footrests with Metal Foot Plates	20 in	350 lbs.	K0001+E1226	1 ea
ACCESSORIES							
146-STDS1040S		Cup Holder	N/A	N/A	N/A	N/A	1 ea, 6 ea/cs



WARNING: This product can expose you to styrene, which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov.