

Universal securement device

The McKesson Grip-Lok universal securement device improves patient comfort and secures a variety of catheters, lines, and tubes. Its unique wrap-around design provides security around the entire tube, unlike tape, and delivers a physical barrier between the skin and the tube. Use the universal securement device for lines and tubes sized 9-24 FR (3-8 mm) in diameter.

Features and benefits

- Allows healthcare professionals and caregivers to open and close the flap to reposition and access the line without removing the universal securement device from the skin
- Attaches and removes quickly and easily with its hook-and-loop design; alcohol can be used to assist with removal
- Delivers superior securement and protects against removal from accidental pulls
- Meets INS guidelines for engineered securement
- Store at 41°-104° F (5°- 40° C)
- Helps reduce the risk of dislodgement of catheters, drains, and lines
- Reduces the risk of allergic reactions and skin irritation through its hypoallergenic, breathable, and latex-free materials
- Keeps device close to the skin with its low-profile design
- Made in USA
- Not made with natural latex



ORDERING INFORMATION

MFR #	DESCRIPTION	SIZE	QTY
18-3300M	Universal Securement Device, Hypoallergenic Acrylic Adhesive, Hook-And-Loop Closure, Sterile	Medium, 3.5 in (8.9 cm)	1 ea, 100 ea/bx, 10 bx/cs

Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.