

## Overbed tables

McKesson overbed tables are a convenient solution for users with reduced mobility. Commonly used in healthcare facilities, overbed tables allow patients to eat or work while in bed and enable healthcare professionals to conveniently access medical equipment from the tabletop. From bedsidings to living spaces, the walnut woodgrain top blends naturally in most décor settings, while the height-adjustable, wheeled design provides transitional use between beds, recliners, and chairs. Low-pressure laminate tabletops are easy-to-clean and feature a large surface area for versatile use, including dining, writing, reading, and using a laptop or tablet. The H-base design is optimal for wheelchair users, allowing the 20-inch clearance required for most manual wheelchairs to roll up to the overbed table. The U-base design minimizes interference with other equipment, accommodating side positioning without floor posts or items under the bed getting in the way.

### Features and benefits

- Conveniently rolls under beds and chairs
- Height-adjustable, wheeled design allows transitional use between beds, recliners, and chairs
- Spring-assisted lift enables smooth gliding height adjustment from seated position and requires minimal strength or mobility
- Lever easily adjusts height from 29.13 to 42.25 inches (74 to 107.3 centimeters)
- Versatile, large table surface measures 28 x 15 inches (71.1 x 38.1 centimeters)
- Four dual-wheel nylon casters swivel for convenient positioning and portability
- Rear-locking wheels keep table in place
- Supports up to 50 lbs
- Not made with natural latex



## ORDERING INFORMATION

MFR #	DESCRIPTION	ATTRIBUTES	QTY
16-11610	Non-Tilt Overbed Table, H-Base	Adjustable Height, Walnut, Chrome	1 ea/cs
16-11620	Non-Tilt Overbed Table, U-Base	Adjustable Height, Walnut, Chrome	1 ea/cs

### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.