

# Clear the way for healing + powerful absorption

## DURABLE PERFORMANCE

**DURAMAX<sup>◊</sup> S** non-occlusive **super-absorbent** dressing allows for a moist environment due to its **breathable properties**.

Protection of periwound skin.\*

The two layers sandwiching the **absorbent** central core strengthen the dressing and keep the polymers in place in the event of damage.

**DURAMAX<sup>◊</sup> S Dressing is comfortable for patients**

**MINIMIZES PAIN:** The cross-linked silicone wound contact layer does not adhere to the wound bed, minimizing pain upon removal or dressing changes.

**MINIMIZES LEAKAGE:** can help to protect periwound skin from wound exudate and excessive moisture.

## MAXIMUM ABSORPTION

**DURAMAX<sup>◊</sup> S Silicone Super Absorbent Dressing manages excess exudate**

With its **retention** capacity, **DURAMAX<sup>◊</sup> S Dressing** prevents maceration and irritation of the wound bed and periwound tissue due to exudate.<sup>1</sup>

**DURAMAX<sup>◊</sup> S Dressing** absorbs and locks in viscous exudate.<sup>1</sup>

**DURAMAX<sup>◊</sup> S Dressing** can absorb and retain fluid while under standard compression.<sup>1</sup>



\* Can remain in place for up to 7 days.

# + A new partner in wound healing

**DURAMAX<sup>®</sup> S Silicone Super Absorbent Dressing** is indicated for malodorous wounds with moderate to heavy exudate levels, including diabetic foot ulcers, venous leg and arterial ulcers, first- and second-degree burns, post-operative wounds and traumatic wounds.

## Unique 6-layer structure

### Polypropylene and polyethylene film

- Permeable to water vapor

### Non-woven polypropylene layers

### Microperforated silicone contact layer<sup>1</sup>

- Atraumatic removal
- Repositionable
- Comfortable for patients



### Super-absorbent polymer core

- Absorbs and locks in exudate, even under compression<sup>1</sup>

Order no.	Size	Each/unit	Unit/case	HCPCS
66023131	4" x 4"	10	10	A6196
66023132	6" x 8"	10	10	A6197
66023134	8" x 10"	10	10	A6198
66023135	8" x 16"	10	10	A6198

Please see Instructions for Use for indications, contraindications, warnings, precautions and other important information.

Reference: 1. Test report No. JCSQ202212511. Mandatory information. DURAMAX S, super-absorbent silicone dressing.

Advanced Wound Management  
Smith+Nephew, Inc.  
Fort Worth, TX 76109 USA  
[www.smith-nephew.com](http://www.smith-nephew.com)

Customer Care Center  
1 800 441-8227  
T 817 900-4000  
F 817 900-4100

©2022 Smith+Nephew, Inc.  
®Trademark of Smith+Nephew  
SCPE44-30185-1022