

Silicone adhesive foam dressings

Silicone adhesive foam dressings: Highly conformable dressing for wounds with light to moderate exudate. Perforated, soft silicone gel adhesive provides gentle, but firm, adhesion and minimizes pain and trauma to the wound and skin upon removal. The dressing may be lifted and repositioned without losing its adherent properties. Breathable, yet moisture- and bacteria-resistant outer film layer seals the wound margins. The absorbent inner hydrocellular foam pad maintains a moist wound environment to help promote healing. McKesson silicone foam dressings are an excellent choice when fragile skin is a concern.

Indications: Wide range of acute and chronic wounds like venous leg ulcers, pressure ulcers, diabetic ulcers and traumatic wounds.

Contraindications: Patients with a known hypersensitivity to the product itself or its components.

Change Frequency: Change dressing as directed by a healthcare professional.



Satisfaction Guaranteed

For complete details, please visit
mms.mckesson.com/mckesson-brands.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

ORDERING INFORMATION

MFR #	DESCRIPTION	QTY	HCPCS	EXUDATE LEVEL*
4841	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 9 in x 9 in	1/ea, 5/bx, 10 bx/cs		●●○
4842	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20 bx/cs	A6212	●●○
4843	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20 bx/cs	A6212	●●○
4844	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20 bx/cs	A6212	●●○
4845	Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Sacral, Gentle Adhesive, Sterile, 7 in x 7 in	1/ea, 10/bx, 20 bx/cs	A6213	●●○
4862	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20 bx/cs	A6209	●●○
4863	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20 bx/cs	A6209	●●○
4864	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20 bx/cs	A6210	●●○
4865	Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Sacral, Gentle Adhesive, Sterile, 7 in x 7 in	1/ea, 10/bx, 20 bx/cs	A6211	●●○
4872	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20 bx/cs	A6212	●○○
4873	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20 bx/cs	A6212	●○○
4874	Lite Hydrocellular Foam Dressings, Silicone Adhesive, With Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20 bx/cs	A6212	●○○
4892	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 3 in x 3 in	1/ea, 10/bx, 20 bx/cs	A6209	●○○
4893	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 4 in x 4 in	1/ea, 10/bx, 20 bx/cs	A6209	●○○
4894	Lite Hydrocellular Foam Dressings, Silicone Adhesive, Without Border, Gentle Adhesive, Sterile, 6 in x 6 in	1/ea, 10/bx, 20 bx/cs	A6210	●○○

EXUDATE LEVEL KEY: ●○○ = Light ●●○ = Moderate ●●● = Heavy

*Exudate level indicates highest level of exudate that the dressing can accommodate