

## Sterilization wraps dual layer

McKesson sterilization wraps provide effective material to wrap instruments and equipment for autoclave sterilization. Designed to deliver excellent barrier properties against both air and waterborne bacteria, the quality SMS construction also offers a high degree of wet or dry tensile strength. Available in various sizes, the sterilization wraps feature soft, cloth-like sheets with two colors bonded together making it easier and more efficient to use.

### Features and benefits

- For ethylene oxide (EO) gas sterilization processes
- 100% polypropylene (SMS) construction featuring spunbond polypropylene on the inner and outer layers and a meltblown polypropylene middle layer
- Protects against both aqueous and airborne contamination
- Consistently provides an excellent barrier maximizing shelf-life and performance
- Dual-layer, two-color design features opposing white and blue sides bonded together
- Improves ease and efficiency required in tray assessments
- Options available varying by grade and size
- Feature 400, 500, and 600-grade options to fit the strength and protection required for the intended task
- Disposable design for single use only
- Non-sterile
- Not made with natural rubber latex



### Satisfaction Guaranteed

For complete details, please visit [mms.mckesson.com/mckesson-brands](https://mms.mckesson.com/mckesson-brands).

### Never Compromise on Care™

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

## ORDERING INFORMATION

MFR #	DESCRIPTION	COLOR	GRADE	SIZE	QTY
16-BW3036	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	400 Grade	36 in x 36 in	36/bg, 2 bgs/cs
16-BW3024	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	500 Grade	24 in x 24 in	60/bg, 2 bgs/cs
16-BW5036	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	500 Grade	36 in x 36 in	36/bg, 2 bgs/cs
16-BW6048	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	600 Grade	48 in x 48 in	24/bg, 1 bg/cs
16-BW6045	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	600 Grade	45 in x 45 in	24/bg, 2 bgs/cs
16-BW3024	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	400 Grade	24 in x 24 in	60/bg, 2 bgs/cs
16-BW6036	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	600 Grade	36 in x 36 in	36/bg, 2 bgs/cs
16-BW5048	Sterilization Wraps Dual Layer, 100% Polypropylene (SMS)	Blue/White	500 Grade	48 in x 48 in	24/bg, 1 bg/cs