

**MATERIAL SAFETY DATA SHEET 0050901400.03.01**

Coloplast® Filtrodor  
Coloplast A/S

Version 1  
1/2003

ANSI Z400.1-1988

Page 1 of 3

**1. Identification of the product and of the company**

**PRODUCT NAME**  
Coloplast® Filtrodor

**PRODUCT CODE(s)**  
0509

**Manufacturer**  
Coloplast A/S  
Høttedam 1  
3050 Humlebæk  
Denmark  
+45 49 11 11 11

**Contact person**  
Product Manager  
770/281/8400

[telephone in Denmark during working hours 8.30 a.m. – 4.30 p.m.]

**2. Composition**

**2.1 OSHA Regulatory Status**

This product contains Copper, CAS 7440-50-8, which is classified as hazardous under OSHA regulations (TWA 1 mg/m3 dusts). Due to its containment in a polymer matrix, the substance poses no immediate hazard.

**2.2 Other components:**

**CAS no.:**

Filters:

Activated carbon	7440-44-0
Copper	7440-50-8
Polyurethane foam	68400-67-9
Siliconized film	—

Needle:

Stainless Steel —

Packaging

Corrugated Cardboard	—
Paper	—
Polyethylene (PE)	9002-88-4
Aluminium foil	—

Ingredients not identified are proprietary and/or non-hazardous. All chemical ingredients appear on the EPA TSCA inventory.

**3. Hazards identification**

The product consists primarily of polymer materials. Product poses no immediate hazard, but can give rise to toxic fumes in a fire.

**4. First aid measures**

<i>Skin contact</i>	This product is intended to be in contact with the skin as directed in the instruction for use
<i>Eye contact</i>	Not relevant
<i>Inhalation</i>	Not relevant
<i>Ingestion</i>	Not relevant
<i>Information</i>	Show this Safety Data Sheet to a physician or emergency ward

**5. Fire fighting measures**

<b>Flash point (°F)</b>	Not relevant
<b>Risks</b>	Smoke, soot, CO/CO <sub>2</sub> , and other toxic fumes. Avoid contact with heated or molten Product. Do not breathe vapors from heated product.
<b>Extinguishing media</b>	Use water spray, carbon dioxide, dry chemical powder or appropriate foam.
<b>Protective measures</b>	Keep personnel removed from and upwind of fire. Wear full firefighting turn-out gear (full Bunker gear), and respiratory protection (SCBA).

I-401322 C-0509  
E-244987 C-0509 D6m

**MATERIAL SAFETY DATA SHEET 0050901400.03.01**Coloplast® Filtrador  
Coloplast A/SVersion 1  
1/2003

ANSI Z400.1-1998

Page 2 of 3

**6. Accidental release measures**

Protection of persons	No specific measures.
Protection of environment	See section 13
Clean-up operation	No specific measures.

**7. Handling and storage***Handling*  
See section 8.*Storage*

Avoid direct sunlight. Store products dry and at room temperatures to maintain product integrity.

**8. Exposure controls/personal protection**

Eye/face protection	None required
Skin protection	None required
Respiratory protection	None required

**9. Physical and chemical properties**

Appearance	Round Filters and single packed Needle
Flash point (°F)	Not relevant
Boiling point (°F)	Not relevant
Solubility	Not relevant
Odor	None

**10. Stability and reactivity**

The product is stable under ambient conditions. High temperatures to be avoided.

**11. Toxicological information**

This product is a medical device and has been tested according to the requirements of medical devices. According to current knowledge this product is considered non-toxic. For further information please contact Coloplast A/S at the telephone number found in section 1.

**12. Ecological information**

Ecotoxicological data have not been determined specifically for this product. Based on toxicity data on the ingredients no ecotoxicological effects are expected. However the product is not biodegradable and discharge to the environment should be avoided.

For further information please contact Coloplast A/S at the telephone number found in section 1.

**13. Disposal considerations**

If this product as supplied becomes a waste, it does not meet the criteria of a hazardous waste as defined under the Resource Conservation and Recovery Act (RCRA) 40 CFR 261.

Under normal private use the product may be disposed of together with other household waste per RCRA 40 CFR 261.4.B1.

The disposal should always be in compliance with Federal, State and local regulations.

**14. Transport Information**

Not dangerous goods (ADR 2003)

#20 p.3

**MATERIAL SAFETY DATA SHEET 0050901400.03.01**

Coloplast® Filtrodor  
Coloplast A/S

Version 1  
1/2003

ANSI Z400.1-1998

Page 3 of 3

**15. Regulatory information**

This product is a medical device and is regulated under the Federal Food Drug and Cosmetic Act. This MSDS is supplied as an additional service. Medical devices do not require an MSDS for hazard communication as stated in 29 CFR 1910.1200.

California Prop65: None

Further information see section 2

**16. Other information**

THE ABOVE INFORMATION HAS BEEN COMPILED FROM SOURCES BELIEVED TO BE RELIABLE AND IS ACCURATE TO THE BEST OF OUR KNOWLEDGE. HOWEVER, COLOPLAST CORP. CANNOT GIVE ANY GUARANTEES REGARDING INFORMATION FROM OTHER SOURCES AND EXPRESSLY DOES NOT MAKE ANY WARRANTIES, NOR ASSUMES ANY LIABILITY, FOR ITS USE.