MSKESSON

Nasal Cannulas

McKesson nasal cannulas are used for delivery of lower flows of oxygen and are ideal for long-term oxygen therapy. They are available in a variety of tube lengths and configurations. All cannulas are designed for single patient use. The line also includes foam ear cushions, which can be added to the tubing for additional comfort. Ear cushions slide easily onto tubing and put less pressure directly on ear to help prevent chafing and tissue injury. O2/CO2 Sampling Cannulas offer accurate, easy to use cannulas for CO2 monitoring applications.

Features and Benefits

- · Made of injection-molded soft plastic for patient comfort
- Anatomical tips and over-the-ear placement provide an easy fit helping to reduce the risk of irritation and pressure sores
- "No Crush" lumen tubing helps maintain open flow of oxygen
- Universal oxygen connectors ensure compatibility with all types of oxygen concentrators
- Nasal prongs are available in straight, contoured, and extra comfort styles
- Straight prongs help ensure consistent delivery of oxygen
- Contour prongs have a soft feel to help prevent skin shearing
- Over-the-ear comfort design helps minimize friction and aid in the reduction of pressure ulceration
- Ear cushions are disposable
- O2/CO2 Sampling cannulas in various lengths, sizes and configurations for CO2 monitoring while administering conscious sedation or post-op recovery
- Not made with natural rubber latex

■ Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.



Never Compromise on Care SM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2020 McKesson Medical-Surgical Inc.

MSKESSON

ORDERING INFORMATION

MFR#	DESCRIPTION	QTY
	NASAL CANNULAS	
32636	Adult, Straight Tip w/Stabilizer Slide (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
16-3302E	Adult, Straight Tip, Universal Connector (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
16-331E	Pediatric, Extra Comfort Plus, Contour Tip (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
16-3318E	Adult, Contoured Tip, Universal Connector (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
32637	Adult, Contoured Tip, (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
32639	Adult, Contoured Tip, (14' Kink Resistant Tubing)	1 ea, 5 ea/cs
32640	Adult, Contoured Tip, (25' Kink Resistant Tubing)	1 ea, 5 ea/cs
32638	Adult, Over-the-Ear Design, Extra Comfort Plus (7' Kink Resistant Tubing)	1 ea, 5 ea/cs
16-0055-50	Adult, Contoured Tip, (50' Crush Resistant Tubing)	1 ea, 20 ea/cs
	NASAL CANNULA EAR CUSHIONS	
32648	Ear Cushions, Adult	1 pair/pk, 50 pk/cs
32649	Adult Cannula, Contour Tip w/Ear Cushions pre-attached (7' Kink Resistant Tubing)	1 ea, 25 ea/cs
	CO ₂ SAMPLING NASAL CANNULAS	
16-0503	Adult, 7' O ₂ / CO ₂ Sampling Line w/Male Luer	1 ea, 25 ea/cs
16-0504	Adult, 7' O ₂ / CO ₂ Sampling Line w/Female Luer	1 ea, 25 ea/cs
16-0581	Adult, 7' $\rm O_2$ Line, 2" $\rm CO_2$ Sampling Line w/Male Luer, w/2" Female x Female Luer Adaptor	1 ea, 25 ea/cs
16-0538	Adult, 10' O ₂ / CO ₂ Sampling Line w/Male Luer	1 ea, 25 ea/cs
16-0539	Adult, 10' O ₂ / CO ₂ Sampling Line w/Female Luer	1 ea, 25 ea/cs
16-0542	Adult, 14' O ₂ / CO ₂ Sampling Line w/Male Luer, w/2" Female x Female Luer Adaptor	1 ea, 25 ea/cs
16-0441	Pediatric, 7' O ₂ / CO ₂ Sampling Line w/Male Luer, w/2" Female x Female Luer Adaptor	1 ea, 25 ea/cs
	ORAL/NASAL CANNULAS	
16-0430	Adult, 10' O ₂ / CO ₂ Crush Resistant Tubing w/Male Luer, w/2" Pigtail w/Female Luers	1 ea, 25 ea/cs