

Loop Electrodes

McKesson Argent® loop electrodes measure 3/32-inches in diameter and feature a stainless steel shaft. They are shipped in a patented tray to maintain their shape and size and have safety grip insulators. Each electrode features tungsten wire for superior shape and wire integrity throughout the excision procedure. This durable construction minimizes thermal damage and yields excellent histopathology results. They are packaged in individual sterile packs.

Features and Benefits

- » Constructed of durable tungsten wire for superior shape and integrity throughout the procedure
- » Standard 3/32-inch stainless-steel shaft designed to fit most electrosurgical pencils
- » Individually packaged in protective tray to prevent damage during shipping
- » Not made with natural rubber latex



MFR #	Description	Specs	QTY
454	Loop Electrodes, LLETZ, Sterile, Single Patient Use	10 mm x 10 mm	1 ea/pouch, 5 pouches/cs
456	Loop Electrodes, LLETZ, Sterile, Single Patient Use	15 mm x 10 mm	1 ea/pouch, 5 pouches/cs
466	Loop Electrodes, LLETZ, Sterile, Single Patient Use	20 mm x 10 mm	1 ea/pouch, 5 pouches/cs
458	Loop Electrodes, LLETZ, Sterile, Single Patient Use	20 mm x 15 mm	1 ea/pouch, 5 pouches/cs
457	Loop Electrodes, LLETZ, Sterile, Single Patient Use	20 mm x 8 mm	1 ea/pouch, 5 pouches/cs

Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit within 12 months from the purchase date. McKesson Brands' satisfaction guarantee is limited to products purchased directly from McKesson. If a McKesson Brands product is not purchased directly from McKesson, any return request or request for refund must be directed to the selling party.

Never Compromise on CareSM

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Electrosurgical Loop Electrodes

McKesson Medical-Surgical
9954 Mayland Drive
Richmond, VA 23233
mms.mckesson.com/mckesson-brands

Post-Acute Care and Home Health
800.347.2456

Home Medical Equipment
888.822.8111

Primary Care
866.625.2679

McKesson Brands Support
800.777.4908