

# Electrosurgical Electrodes

McKesson Argent® coated electrodes are coated with PTFE (polytetrafluoroethylene). Unlike standard electrodes that accumulate discarded blood and tissue, these coated electrodes resist build up and can be easily wiped clean. McKesson Argent coated electrodes are available in needle or blade styles, are for single patient use and are not made of natural rubber latex. McKesson Argent needle electrodes are available standard, extended and straight styles with options for coated tips. The tips are coated with PTFE and move more effectively through the tissue. Any debris can easily be removed with a damp sponge or cloth.

## Features and Benefits

- » Coated with PTFE (polytetrafluoroethylene) to resist eschar build up
- » Available in needle or blade styles
- » Available in standard, extended and straight styles
- » Single patient use
- » Not made with natural rubber latex



**Satisfaction Guaranteed**  
 If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

Description	MFR #	QTY
Blade Electrodes, Extended, Coated, 6"	22-ES04T	1ea/pouch, 12 pouches/bx
Blade Electrodes, Coated, 2.5"	22-ES58T	1ea/pouch, 12 pouches/bx
Blade Electrodes, Modified, Coated, 2.5"	22-ES59T	1ea/pouch, 12 pouches/bx
Blade Electrodes, Extended, Modified, Coated, 6"	22-ES39T	1ea/pouch, 12 pouches/bx
Needle Electrodes, Straight, Coated, 2.75"	22-ES02T	1ea/pouch, 12 pouches/bx
Needle Electrodes, Extended, Coated, 6"	22-ES03T	1ea/pouch, 12 pouches/bx
Needle Electrodes, Modified, Coated, 2.75"	22-ES38T	1ea/pouch, 12 pouches/bx
Needle Electrodes, Extended, Modified, Coated, 6"	22-ES40T	1ea/pouch, 12 pouches/bx
Needle Electrodes, Standard, 2.84"	468	1ea, 50 ea/bx