

## Spinal Needles

McKesson spinal needles feature two options for conducting spinal procedures that maximize patient comfort. The regular spinal needles help to reduce trauma by providing a precise fit between stylet and cannula bevel. Their smooth insertion resistance helps the caregiver to identify tissue layers. The epidural option comes with detachable wings for added flexibility. Their wide, metric marks and metal echogenic stylets offer maximum visibility, accuracy and ease of placement. McKesson Spinal Needles are color coded for easy identification and available in a wide range of colors, walls, points and sizes.

### Features and Benefits

- Color-coded stylets for easy gauge identification and large, clear hubs for enhanced tactile feel
- Computerized grinding for spinal needles to help ensure sharpness, consistency and durability
- Precision point grinding for epidural needles to minimize coring and trauma
- Not made with natural rubber latex



### Satisfaction Guaranteed

If you are not completely satisfied with any McKesson Brands product, you may return it for a full refund or credit.

### Never Compromise on Care<sup>SM</sup>

With McKesson Brands, you'll never have to compromise. We deliver quality products at competitive prices, backed by our 100% satisfaction guarantee. We have decades of experience in the market, and we offer you support and work to help ensure you have the products you need, when you need them—making it even easier for you to provide the best care for your patients.

The product information contained in this document, including the product images and additional product materials, was collected from various supplier sources. All product claims and specifications are those of the product suppliers and have not been independently verified by McKesson Medical-Surgical or its affiliates ("McKesson"). McKesson is not responsible for errors or omissions in the product information.

The properties of a product may change or be inaccurate following the posting or printing of the product information in the document, either in the print or online version. Caution should be exercised when using or purchasing any products from McKesson's online or print documents by closely examining the product packaging and the labeling prior to use.

Due to product changes, information listed in this document is subject to change without notice. This information is placed solely for your convenience in ordering and McKesson disclaims all responsibility for its completeness and accuracy, whether or not the inaccuracy or incompleteness is due to fault or error by McKesson.

© 2018 McKesson Medical-Surgical Inc.

## ORDERING INFORMATION

MFR #	DESCRIPTION	QTY
4628V2	Spinal Needles, Thin Wall, Pink, Quincke Point, 18G x 3½"	1 ea, 25/bx, 4 bx/cs
4627V2	Spinal Needles, Thin Wall, Yellow, Quincke Point, 20G x 6"	1 ea, 25/bx, 4 bx/cs
4629V2	Spinal Needles, Regular Wall, Yellow, Quincke Point, 20G x 3½"	1 ea, 25/bx, 4 bx/cs
4630V2	Spinal Needles, Regular Wall, Black, Quincke Point, 22G x 3"	1 ea, 25/bx, 4 bx/cs
4631V2	Spinal Needles, Regular Wall, Black, Quincke Point, 22G x 3½"	1 ea, 25/bx, 4 bx/cs
4633V2	Spinal Needles, Regular Wall, Black, Quincke Point, 22G x 5"	1 ea, 25/bx, 4 bx/cs
4634V2	Spinal Needles, Regular Wall, Black, Quincke Point, 22G x 7"	1 ea, 25/bx, 4 bx/cs
4632V2	Spinal Needles, Regular Wall, Orange, Quincke Point, 25G x 3½"	1 ea, 25/bx, 4 bx/cs
3076V2	Epidural Needles, Extra Thin Wall, Pink, Tuohy Point, Detachable Wing, 18G x 3½"	1 ea, 25/bx, 4 bx/cs
3077V2	Epidural Needles, Regular Wall, Yellow, Tuohy Point, Detachable Wing, 20G x 3½"	1 ea, 25/bx, 4 bx/cs



**\* WARNING:** Cancer and Reproductive Harm - For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

McKesson Medical-Surgical  
9954 Mayland Drive  
Richmond, VA 23233  
[mms.mckesson.com/mckesson-brands](http://mms.mckesson.com/mckesson-brands)

Post-Acute Care and Home Health  
800.347.2456

Home Medical Equipment  
888.822.8111

Primary Care  
866.625.2679

McKesson Brands Support  
800.777.4908