

Protekt® ProTemp Infrared Non-Contact Thermometer



Temperature taking has never been easier with the Non-Contact Thermometer. The Protekt ProTemp provides convenience, accuracy and precision when taking a temperature. With the push of a button, instantly and accurately take a person's temperature without the need of any physical contact. Simply point the thermometer at the patient's forehead, click "Scan" and their temperature instantly will be displayed.

- Non-touch eliminates cross contamination
- Hygienic and easy to use
- High temperature alert
- Fast and accurate, 1 second response time
- FDA & CE approved

FEATURES AND BENEFITS:

- Designed to take the body temperature of a person without the interference from the heat of the surrounding environment.
- Never wake the patient to take their temperature.
- Audible alarm when the temperature is too high. [$\leq 100.4^{\circ}\text{F}$ (38°C)]
- Memory for 32 previous measurements.
- Back-lit LCD screen.
- Works with 2 AA batteries (not included)

- Data can be displayed in either Celsius or Fahrenheit.
- Body-Mode Temperature Range: $89.6^{\circ}\text{F} \sim 109.2^{\circ}\text{F}$ ($32^{\circ}\text{C} \sim 42.9^{\circ}\text{C}$)
- Measuring Distance: $1.2\text{in} \sim 2\text{in}$ ($3\text{cm} \sim 5\text{cm}$)
- Accuracy: $\pm 0.6^{\circ}\text{F}$ (0.3°C)
- Automatic shut off.
- Practical and easy to use.
- 6 Month limited warranty.

40011	Non-Contact Infrared Thermometer
-------	----------------------------------



270 Washington Street, Mount Vernon, NY 10553
Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005
www.proactivemedical.com • sales@proactivemedical.com