

# SEATING AND POSITIONING

 **Proactive**  
medical products



270 Washington Street Mount Vernon, NY 10553  
Tel: 855-237-7622 • Fax: 845-205-6005  
[www.proactivemedical.com](http://www.proactivemedical.com)

## PROTEKT® FOAM CUSHION (E2601)

The Protekt® Foam Cushion contours to the patient's body, which reduces pressure on vulnerable areas, and assists with the prevention and management of pressure ulcers. The high-performance, breathable stretch nylon cover, designed for enhanced patient comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 1 YEAR NON-PRORATED WARRANTY



70001	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x2"
70002	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x2"
70003	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x2"
71001	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"
71002	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"
71003	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x18"x3"
71004	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"
71008	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x3"
72001	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"
72002	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"
72003	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"
72060	Foam Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"

## PROTEKT® BARIATRIC FOAM CUSHION (E2602)

71005	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x3"
71020	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x3"
71006	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x3"
71012	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x3"
71007	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x3"
71017	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x20"x3"
71019	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x18"x3"
71010	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x20"x3"
71026	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x18"x3"
71016	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x20"x3"
72004	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"
72062	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x4"
72005	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"
72106	Foam Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x4"
72006	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x18"x4"
72107	Foam Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x20"x4"
72063	Foam Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	28"x18"x4"
72007	Foam Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	28"x20"x4"
72064	Foam Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	30"x18"x4"
72046	Foam Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	30"x20"x4"

## PROTEKT® FOAM WEDGE CUSHION

The Protekt® Foam Wedge Cushion offers all the features of the Protekt® Foam Cushion, with the added benefit of a wedge shape for posture support. This design enhances safety for patients who are prone to sliding down and potentially falling while seated. It also provides support for patients who require frequent hip repositioning. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 1 YEAR NON-PRORATED WARRANTY



76149	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"-1.5"
76150	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"-1.5"
76152	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"-1.5"
76056	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"-2"
76161	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"-2"
76162	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"-2"
76053	Foam Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"-2"

## PROTEKT® BARIATRIC FOAM WEDGE CUSHION

76168	Foam Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x16"x4"-2"
76054	Foam Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"-2"
76055	Foam Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"-2"
76062	Foam Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x4"-2"

## PROTEKT® GEL CUSHION (E2603)

The Protekt® Gel-Cushion is designed to redistribute pressure away from vulnerable contact points, improve patient comfort, and protect against friction and shear. Ideal for use in standard or bariatric wheelchairs and geri chairs. All gel-bladders are dual chamber, which ensure even pressure distribution to support and protect areas prone to pressure ulcers. This design ideally maintains a cool surface temperature, which helps prevent skin breakdown. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 18 MONTH NON-PRORATED WARRANTY

73001	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x2"
73002	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x2"
73003	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x2"
74001	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"
74002	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"
74003	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x18"x3"
74004	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"
74014	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x3"
75001	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"
75002	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"
75003	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"
75007	Gel Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"

**Dual Chamber Gel-Bladder** is designed to eliminate gel migration and regulate body temperature.

**Non-contoured Seating Surface** enables easy patient transfers.



**Easy-Snap Buckle** conveniently secures the cushion to any wheelchair surface for added safety.

**Breathable Stretch Cover**

- Encourages healing and comfort by removing moisture and reducing shear.
- Anti-microbial. Anti-bacterial, water resistant & flame retardant (CA117).

## PROTEKT® BARIATRIC GEL CUSHION (E2603)

74005	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x3"
74019	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x3"
74006	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x3"
74008	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x3"
74007	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x3"
74020	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x20"x3"
74015	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x18"x3"
74021	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x20"x3"
74016	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x18"x3"
74017	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x20"x3"
75004	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"
75008	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x4"
75005	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"
75009	Gel Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x4"
75006	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x18"x4"
75010	Gel Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x20"x4"
75035	Gel Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	28"x18"x4"
75016	Gel Cushion w/Stretch Nylon Cover	750 lb wt. cap.	28"x20"x4"
75031	Gel Cushion w/Stretch Nylon Cover	750 lb wt. cap.	30"x18"x4"
75017	Gel Cushion w/Stretch Nylon Cover	750 lb wt. cap.	30"x20"x4"

## PROTEKT® GEL WEDGE CUSHION



\* 18 MONTH NON-PRORATED WARRANTY

The Protekt® Gel Wedge Cushion offers all the features of the Protekt® Gel Cushion, with the added benefit of a wedge shape for posture support. This design enhances safety for patients who are prone to sliding down and potentially falling while seated. It also provides support for patients who require frequent hip repositioning. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

76149G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"-1.5"
76150G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"-1.5"
76152G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"-1.5"
76056G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"-2"
76161G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"-2"
76162G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"-2"
76053G	Gel Wedge Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"-2"

76168G	Gel Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x16"x4"-2"
76054G	Gel Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"-2"
76055G	Gel Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"-2"
76062G	Gel Wedge Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x4"-2"

## PROTEKT® BARIATRIC GEL WEDGE CUSHION

## PROTEKT® FOAM COCCYX CUSHION

The Protekt® Foam Coccyx Cushion offers all the features of the Protekt® Foam Cushion, with the added benefit of a coccyx cut-out that is specifically designed for patients at risk for pressure ulcers in the coccyx area. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 1 YEAR NON-PRORATED WARRANTY



70001CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x2"
70002CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x2"
70003CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x2"
71001CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"
71002CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"
71003CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x18"x3"
71004CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"
71008CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x3"
72001CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"
72002CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"
72003CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"
72060CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"

## PROTEKT® BARIATRIC FOAM COCCYX CUSHION

71005CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x3"
71020CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x3"
71006CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x3"
71012CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x3"
71007CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x3"
71017CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	26"x20"x3"
71019CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x18"x3"
71010CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	28"x20"x3"
71026CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x18"x3"
71016CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	30"x20"x3"
72004CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"
72062CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x4"
72005CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"
72106CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x4"
72006CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x18"x4"
72107CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	600 lb wt. cap.	26"x20"x4"
72063CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	28"x18"x4"
72007CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	28"x20"x4"
72064CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	30"x18"x4"
72046CX	Foam Coccyx Cushion w/ Stretch Nylon Cover	750 lb wt. cap.	30"x20"x4"

## PROTEKT® FOAM WEDGE CUSHION WITH POMMEL

The Protekt® Foam Wedge Cushion with Pommel offers all the features of the Protekt® Foam Wedge Cushion, with the added benefit of a pommel to promote posture support and proper hip positioning. The pommel also assists with hip abduction to prevent hip rotation. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 1 YEAR NON-PRORATED WARRANTY



76149P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"-1.5"
76150P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"-1.5"
76152P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"-1.5"
76056P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"-2"
76161P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"-2"
76162P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"-2"
76053P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"-2"

76168P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	500 lb wt. cap.	22"x16"x4"-2"
76054P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"-2"
76055P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"-2"
76062P	Foam Wedge Cushion w/ Pommel w/Stretch Nylon Cover	500 lb wt. cap.	26"x18"x4"-2"

## PROTEKT® BARIATRIC FOAM WEDGE CUSHION W/ POMMEL

## PROTEKT® GEL COCCYX CUSHION

The Protekt® Gel Coccyx Cushion offers all the features of the Protekt® Gel Cushion, with the added benefit of a coccyx cut-out that is specifically designed for patients at risk for pressure ulcers in the coccyx area. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 18 MONTH NON-PRORATED WARRANTY



73001CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x2"
73002CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x2"
73003CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x2"
74001CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"
74002CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"
74003CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	18"x18"x3"
74004CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"
74014CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x3"
75001CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"
75002CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"
75003CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"
75007CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"

## PROTEKT® BARIATRIC GEL COCCYX CUSHION

74005CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x3"
74019CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x3"
74006CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x3"
74008CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x3"
74007CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x3"
74020CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	26"x20"x3"
74015CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	28"x18"x3"
74021CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	28"x20"x3"
74016CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	30"x18"x3"
74017CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	30"x20"x3"
75004CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"
75008CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x4"
75005CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"
75009CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x4"
75006CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	26"x18"x4"
75010CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	600 lb wt. cap.	26"x20"x4"
75035CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	750 lb wt. cap.	28"x18"x4"
75016CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	750 lb wt. cap.	28"x20"x4"
75031CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	750 lb wt. cap.	30"x18"x4"
75017CX	Gel Cushion w/ Coccyx w/ Stretch Nylon Cover	750 lb wt. cap.	30"x20"x4"

## PROTEKT® GEL WEDGE CUSHION WITH POMMEL



\* 18 MONTH NON-PRORATED WARRANTY

The Protekt® Gel Wedge Cushion with Pommel offers all the features of the Protekt® Gel Wedge Cushion, with the added benefit of a pommel to promote posture support and proper hip positioning. The pommel assists with hip abduction to prevent hip rotation. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

76149GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"-1.5"
76150GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"-1.5"
76152GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"-1.5"
76056GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x4"-2"
76161GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x4"-2"
76162GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x4"-2"
76053GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x4"-2"

76168GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	22"x16"x4"-2"
76054GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x4"-2"
76055GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x4"-2"
76062GP	Gel Wedge Cushion w/ Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x4"-2"

## PROTEKT® BARIATRIC GEL WEDGE CUSHION W/ POMMEL

## PROTEKT® SUPREME MOLDED FOAM CUSHION (E2605)



\* 18 MONTH NON-PRORATED WARRANTY

The Protekt® Supreme Molded Foam Cushion is a high-quality molded foam, specially contoured to maximize stabilization, promote proper positioning, prevent skin breakdown, and enhance patient comfort. Lateral and medial leg contours promote hip alignment and postural symmetry, while a raised front contour controls forward sliding. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

76001SP	Contour Molded Foam Cushion	275 lb wt. cap.	16"x16"x3"
76002SP	Contour Molded Foam Cushion	275 lb wt. cap.	18"x16"x3"
76003SP	Contour Molded Foam Cushion	275 lb wt. cap.	18"x18"x3"
76004SP	Contour Molded Foam Cushion	275 lb wt. cap.	20"x16"x3"
76008SP	Contour Molded Foam Cushion	275 lb wt. cap.	20"x18"x3"

## PROTEKT® BARIATRIC SUPREME MOLDED FOAM CUSHION (E2606)

76005SP1	Contour Molded Foam Cushion	500 lb wt. cap.	22"x16"x3"
76005SP	Contour Molded Foam Cushion	500 lb wt. cap.	22"x18"x3"
76006SP	Contour Molded Foam Cushion	500 lb wt. cap.	24"x18"x3"
76007SP	Contour Molded Foam Cushion	500 lb wt. cap.	26"x18"x3"

## PROTEKT® ANTI-THRUST FOAM CUSHION



\* 18 MONTH NON-PRORATED WARRANTY

The Protekt® Anti-Thrust Foam Cushion is a high-quality positioning cushion that utilizes an anti-thrust barrier to eliminate forward and sacral sitting problems. High-density foam increases comfort and relieves pressure. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

76206	Protekt Anti-Thrust Foam Cushion	16"x16"x5"x3"
76201	Protekt Anti-Thrust Foam Cushion	18"x16"x5"x3"
76202	Protekt Anti-Thrust Foam Cushion	20"x16"x5"x3"
76205	Protekt Anti-Thrust Foam Cushion	20"x18"x5"x3"
76208	Protekt Anti-Thrust Foam Cushion	22"x18"x5"x3"
76209	Protekt Anti-Thrust Foam Cushion	24"x18"x5"x3"
76210	Protekt Anti-Thrust Foam Cushion	26"x18"x5"x3"

## PROTEKT® ANTI-THRUST CUSHION WITH GEL



\* 18 MONTH NON-PRORATED WARRANTY

The Protekt® Anti-Thrust Cushion with Gel is a high-quality positioning cushion that utilizes an anti-thrust barrier to eliminate forward and sacral sitting problems. All gel-bladders are dual chamber, which ensure even pressure distribution to support and protect areas prone to high pressure. This design ideally maintains a cool surface temperature, which helps prevent skin breakdown. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

76206G	Protekt Anti-Thrust Cushion with Gel	16"x16"x5"x3"
76201G	Protekt Anti-Thrust Cushion with Gel	18"x16"x5"x3"
76202G	Protekt Anti-Thrust Cushion with Gel	20"x16"x5"x3"
76205G	Protekt Anti-Thrust Cushion with Gel	20"x18"x5"x3"
76208G	Protekt Anti-Thrust Cushion with Gel	22"x18"x5"x3"
76209G	Protekt Anti-Thrust Cushion with Gel	24"x18"x5"x3"
76210G	Protekt Anti-Thrust Cushion with Gel	26"x18"x5"x3"


## THE PROTEKT® ULTRA GEL-CUSHION WITH VGPT™ TECHNOLOGY (E2607)

The Propekt® Ultra Gel-Cushion with VGPT™ offers the ultimate level in patient skin protection and posture support. Clinically designed to evenly distribute a patient's weight, the Visco Gel Polymer Technology (VGPT™) is a revolutionary support surface that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam. It is constructed with a dual-chamber gel bladder and two layers of foam -- a contoured, molded foam base and a layer of gel-infused memory foam for added skin protection and comfort. The dual-chamber gel bladder ensures even gel distribution to areas vulnerable to pressure ulcers and maintains a cool surface temperature to prevent skin breakdown. A contoured surface with tapered leg abductors supports proper hip positioning and provides greater surface area to optimize skin protection. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection.

\* 18 MONTH NON-PRORATED WARRANTY



77001	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	16"x16"x3"
77002	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	18"x16"x3"
77008	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	18"x18"x3"
77003	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	20"x16"x3"
77009	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	275 lb wt. cap.	20"x18"x3"



## PROTEKT® ULTRA GEL-CUSHION WITH VGPT™ TECHNOLOGY (E2608)

77004	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	22"x18"x3"
77010	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	22"x20"x3"
77005	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	24"x18"x3"
77011	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	24"x20"x3"
77006	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	26"x18"x3"
77012	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	26"x20"x3"
77014	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	28"x20"x3"
77013	Molded Foam Base w/ Gel Bladder, Visco Gel-Foam and Pommel w/ Stretch Nylon Cover	500 lb wt. cap.	30"x20"x3"

## PROTEKT® O2 DRY FLOTATION AIR CUSHION (E2624/E2625)

The Protekt® O2 cushion is an excellent choice for the prevention and treatment of pressure ulcers. The soft, interconnected air cells allow for deep immersion, providing constant low pressure across the seating surface. Air-adjustable 4" or 2" cells provide optimal pressure redistribution for skin tissue protection and pressure ulcer management, while dual-valve compartments address seating challenges, such as pelvic obliquities and sacral seating, and promote proper positioning. It has a nylon breathable cover with non-skid bottom and safety straps.

\* 5 YEAR NON-PRORATED WARRANTY



78001	Dry Flotation Air Cushion	500 lb wt. cap.	16"x16"x2"
78002	Dry Flotation Air Cushion	500 lb wt. cap.	18"x16"x2"
78007	Dry Flotation Air Cushion	500 lb wt. cap.	18"x18"x2"
78003	Dry Flotation Air Cushion	500 lb wt. cap.	20"x16"x2"
78011	Dry Flotation Air Cushion	500 lb wt. cap.	20"x18"x2"
78012	Dry Flotation Air Cushion	500 lb wt. cap.	22"x18"x2"
78013	Dry Flotation Air Cushion	500 lb wt. cap.	22"x20"x2"
78004	Dry Flotation Air Cushion	500 lb wt. cap.	16"x16"x4"
78005	Dry Flotation Air Cushion	500 lb wt. cap.	18"x16"x4"
78008	Dry Flotation Air Cushion	500 lb wt. cap.	18"x18"x4"
78006	Dry Flotation Air Cushion	500 lb wt. cap.	20"x16"x4"
78009	Dry Flotation Air Cushion	500 lb wt. cap.	20"x18"x4"
78014	Dry Flotation Air Cushion	500 lb wt. cap.	22"x18"x4"
78015	Dry Flotation Air Cushion	500 lb wt. cap.	22"x20"x4"

## PROTEKT® GERI-CHAIR FOAM OVERLAY WITH VGPT™ TECHNOLOGY


The Protekt® Geri Chair Foam Overlay contours to the patient's body, which reduces pressure on vulnerable areas and enhances patient comfort. The Visco Gel Polymer Technology (VGPT™) offers a gel-infused memory foam top layer that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial, and soft to the touch. The multi-directional stretch design alleviates pressure on bony prominences and provides optimal skin protection. The optional dual-chamber gel bladder ensures even gel distribution to areas vulnerable to pressure ulcers and maintains a cool surface temperature to prevent skin breakdown.

\* 1 YEAR NON-PRORATED WARRANTY

**Dual Chamber Gel-Bladder** is designed to eliminate gel migration and regulate body temperature (optional)



VGPT™ reduces occurrence of pressure sores by channeling heat away from the body significantly faster than traditional foam, virtually eliminating perspiration.



(Pictured with gel option)

90011	Geri-Chair Overlay w/ Gel Infused Visco Memory Foam and Gel Bladder	275 lb wt. cap.	19"x72"x3"
90012	Geri-Chair Overlay w/ Gel Infused Visco Memory Foam and Gel Bladder	275 lb wt. cap.	24"x72"x3"
90011F	Geri-Chair Overlay w/ Gel Infused Visco Memory Foam	500 lb wt. cap.	19"x72"x3"
90012F	Geri-Chair Overlay w/ Gel Infused Visco Memory Foam	500 lb wt. cap.	24"x72"x3"

## PROTEKT® SEAT RELIEF ALTERNATING PRESSURE CUSHION WITH VGPT™ TECHNOLOGY



The Protekt™ Seat Relief Alternating Pressure Seat Cushion System has four alternating bladders that provide the patient a superior pressure redistribution seating surface. The patented Visco Gel Polymer Technology (VGPT™) foam topper cools the seating surface and maximizes the pressure redistribution of the air bladders. The Protekt® Seat Relief redistributes pressure over the entire contact interface and stimulates capillary blood flow for the prevention and management of pressure ulcers. The water-resistant stretch cover is anti-microbial and fire-retardant. The non-skid bottom keeps the cushion in place ensuring proper patient positioning.

\* 1 YEAR NON-PRORATED WARRANTY

80120	Alternating Pressure Cushion w/Pump	275 lb wt. cap.	16"x16"x3"
80121	Alternating Pressure Cushion w/Pump	275 lb wt. cap.	18"x16"x3"
80124	Alternating Pressure Cushion w/Pump	275 lb wt. cap.	18"x18"x3"
80122	Alternating Pressure Cushion w/Pump	275 lb wt. cap.	20"x16"x3"

## PROTEKT® AIRE GERI-CHAIR ALTERNATING PRESSURE OVERLAY



Protekt® Aire Geri Chair Alternating Pressure Overlay is an alternating pressure recliner geri chair air mattress system. This adaptable system allows patients, who are at high risk for pressure ulcers, to be seated in a variety of types of chairs. Patients are no longer confined to bed due to risks of sitting without appropriate pressure redistribution. The Protekt Aire Geri Chair Overlay System can be integrated with a geriatric chair, recliner, or standard high-back wheelchair. Thirty-one bladders gently alternate, providing low interface pressures, in addition to comfortably supporting the patient. The static function provides safe egress or ingress from the system.

\* 1 YEAR NON-PRORATED WARRANTY

80018	Geri-Chair Alternating Pressure Geri-Chair Overlay with Pump - 20"	275 lb wt. cap.	69"x20"x3"
80019	Geri-Chair Alternating Pressure Geri-Chair Overlay with Pump - 24"	500 lb wt. cap.	69"x24"x3"

# PROTEKT® GAP FILLERS ENTRAPMENT PROTECTION

## ENSURE THAT ALL GAPS AROUND THE MATTRESS ARE FILLED

Protekt® Gap Fillers are custom-made to fit any gap and protect the resident from bed entrapment. They are ideal for mitigating the danger of entrapment in Zones 3 and 7. By closing the gaps around the perimeter of the mattress, the risk of entrapment of the head, neck, chest or limbs that may occur at the headboard, footboard and rails is eliminated.

### General features of the Protekt® Gap Fillers:

- Fill the gaps between the mattress and headboard / footboard / side rails
- Help prevent head, neck, chest and limb entrapment in Zones 3 and 7
- Easy-to-clean vinyl cover
- Custom-made to any size ensuring a perfect fit every time
- Meets CA #117
- 1 Year Non-Prorated Warranty



### Protekt® Rail Guard

- Prevents entrapment along the side of the mattress
- Fills the gap between the mattress and the rail
- Attaches securely to the rail with buckles
- Works with most rails when raised or lowered

58300	Protekt® Rail Guard	Half – 36"x6"x3"
58301	Protekt® Rail Guard	Full – 70"x6"x3"

### Protekt® Mattress Extender

- Safely accommodates taller patients
- Prevents entrapment along the head/foot board
- Fills the gap between the mattress and the foot board
- Optional buckles secure the extender to the head/foot board

65010	Protekt® Mattress Extender	36"x6"x3"
65031	Protekt® Mattress Extender	36"x6"x4"
65032	Protekt® Mattress Extender	36"x6"x5"



### Protekt® Wedge Rail Guards (pair)

- Prevent entrapment along the side of the mattress
- Fill the gap between the mattress and the rail
- Double as a side rail gap filler and a patient positioning device
- Attaches securely to the rail with buckles

58400	Protekt® Wedge Rail Guards (pair)	Half – 36" (pr)
58401	Protekt® Wedge Rail Guards (pair)	Full – 70" (pr)

### Protekt® Roll Control Bolster

- Alternative to side rails.
- Prevents the patient from falling out of bed.
- Easily attaches to bed with quick-release buckles and can be positioned to accommodate individual patient needs
- Provides secure positioning for side-lying patients.
- Not considered a restraint.

58500	Protekt® Roll Control Bolster	Double 34"Lx7"Hx8"D
-------	-------------------------------	---------------------



## PROTEKT® HEEL LOFT PILLOW WITH VGPT™



60150	Heel Loft Pillow with VGPT™	24"x16"x7"x5.5"
60150G	Heel Loft Pillow with VGPT™ with Gel Bladder	24"x16"x7"x5.5"

The Protekt® Heel Loft Pillow completely off-loads pressure on the heels by comfortably elevating the feet and suspending the heels. In addition to inhibiting the development of pressure ulcers, it also alleviates lower back pain by relaxing the spine. The Visco Gel Polymer Technology (VGPT™) offers a gel-infused memory foam top layer that provides unprecedented pressure redistribution, rapid heat dissipation, and vast increases in breathability compared to traditional foam redistribution. The high-performance, breathable stretch nylon cover, designed for pressure relief and comfort, is water-proof, anti-microbial and soft to the touch.

- Completely removes pressure from the heels
- Prevents the formation of pressure ulcers
- Helps heal existing pressure ulcers
- Breathable stretch nylon cover reduces shearing
- Non-skid bottom prevents sliding
- Optional gel bladder adds superior comfort and cooling
- 1 Year Non-Prorated Warranty

## PROTEKT® POSITION WEDGE



50030	30° Positioning Wedge 7"x12"x18"
50031	30° Positioning Wedge 7"x12"x24"
50032	30° Positioning Wedge 7"x12"x34"

The Protekt® Positioning Wedge is used to comfortably reposition immobile patients. By reducing pressure on vulnerable areas, it assists with the prevention and management of pressure ulcers. The wedge is designed with a 30-degree lateral incline to help with proper positioning and to relieve pressure on the hip bone and coccyx. The high-performance breathable stretch nylon cover, designed for pressure relief and comfort, is fluid-proof, anti-microbial and soft to the touch.

- Recommended incline angle of 30° helps relieve pressure on hip bone and coccyx
- High-density polyurethane foam construction
- Breathable stretch nylon cover reduces shearing
- Non-skid bottom prevents sliding
- 1 Year Non-Prorated Warranty

## PROTEKT FOOTREST EXTENDER/LEGREST PAD

The Protekt® Footrest Extender provides a base of support and elevates the legs one or two inches, improving seat alignment, posture, and comfort while alleviating foot drop symptoms.

- Protects lower leg from contact with wheelchair
- Prevents feet from slipping behind footrests
- Mitigates foot drop with removable extender that raises the footrest by 1, 2 or 3 inches
- Secures easily with quick-release straps
- 1 Year Non-Prorated Warranty



47001	Footrest Extender/Leg Rest Pad w/1" Pad, 16"-20"
47002	Footrest Extender/Leg Rest Pad w/2" Pad, 16"-20"
47003	Footrest Extender/Leg Rest Pad w/3" Pad, 16"-20"
47004	Footrest Extender/Leg Rest Pad w/1" Pad, 22"-24"
47005	Footrest Extender/Leg Rest Pad w/2" Pad, 22"-24"
47006	Footrest Extender/Leg Rest Pad w/3" Pad, 22"-24"

## PROTEKT LEGREST PAD

The Protekt® Legrest Pad provides a cushion support behind the calves and prevents the feet from getting caught under the wheelchair by blocking them from slipping behind the footrest.

- Protects lower leg from contact with wheelchair
- Prevents feet from slipping behind footrests
- Secures easily with quick-release straps
- 1 Year Non-Prorated Warranty



47020	Leg Rest Pad Only 16"-20"
47024	Leg Rest Pad Only 22"-24"



---

270 Washington Street Mount Vernon, NY 10553  
Tel: 855-237-7622 • Fax: 845-205-6005  
sales@proactivemedical.com • www.proactivemedical.com