

Be Independent.™



# CURVE GRAB BAR™

BECAUSE YOUR GRAB BAR CAN PIVOT TO YOU



PRODUCT #9000



WATCH A DEMO





Mounts to Wall for Strong Support



Ergonomic Ladder Grip



Pivots and Locks Against Wall



Handrail for Additional Support

## FEATURES & BENEFITS

- **Promotes Safety & Mobility** - Protects against falls and makes standing or sitting easy.
- **Adaptable** - Can be pivoted to the position most convenient for the user.
- **Horizontal Bar** - Provides extra stability.
- **Space-Saving Design** - Pivots and locks against the wall when not in use.

## SPECIFICATIONS

Materials: Rust resistant zinc plated steel

5 Locking Positions

Product Length: 40"

Length of Curve Handle: 18.75"

Length of Horizontal Base Bar: 14"

Length of Extension From the Wall at 90° Position: 12.5"

Diameter: 1.25"

Weight Capacity: 300 lbs

Weight of Product: 12 lbs

Package Dimensions: 29" L x 14.25" W x 3.25" H; 12 lbs

Case of 4 Dimensions: 29" L x 14.25" W x 13" H; 50 lbs

## EASY ASSEMBLY (No Additional Tools Required)



Insert Top and Bottom Mounts



Attach Horizontal Support Bar



Secure to the Wall