

home assistance products



GUARDIAN[®]

Walking Aids

Bath Safety

Home Comfort





GUARDIAN *Select* Your customers have the choice of enhanced features on products with the Guardian Select designation.

GUARDIAN *Plus+* Easily find products that have a weight capacity of 500 lbs. and greater.

Medline Industries is one of the leading manufacturers of home medical equipment and is dedicated to providing innovative, high quality products including the Guardian® brand.

When it comes to walking aids and personal care products, no company is more committed to the needs of health care providers and their patients, than Medline Industries. The Guardian® brand has set the standard for uncompromising quality and product selection for over sixty years.

UNCOMPROMISING QUALITY

Guardian® is committed to delivering the highest level of quality in every item that bears our name. The Guardian® mark is your assurance that each product has been built to meet our stringent standards because your patients' independence and safety are dependent on our assistive devices.

UNMATCHED PRODUCT SELECTION

From pediatrics to bariatrics, the Guardian® product offering is unmatched. We provide a wide range of products to meet your customers' needs.

WARRANTED TO LAST

All Guardian® DME products offer the peace-of-mind of a generous Lifetime Limited Warranty against product defects.

Ambulatory Aids

- 01 Guardian Pediatric Walking Aids
- 02 Guardian Crutches
- 03 Guardian Single Point & Quad Canes
- 04 Guardian Rolling Walkers
- 06 Guardian Folding Walkers
- 08 Guardian Walker Accessories
- 09 Guardian Rubber Accessories
- 10 Specifications

Bath Safety and Home Assistance Products

- 12 Guardian Bath Safety
- 14 Guardian Bathing Products
- 15 Guardian Grab Bars
- 16 Guardian Home Assistance Products
- 17 Guardian Toilet Assist Aids
- 18 Guardian Commodes
- 20 Specifications

Home Comfort Products

- 21 Medline Home Comfort Products



Model G07782

STRIDER® PRESCRIPTION WALKER

- Accommodates most physical requirements with the optional platform attachments.
- Handles adjust without tools to an infinite number of positions to provide stability and posture. 90° bend standard (Tyke models are 45°).
- Extra wide base provides stability, safety and balance.
- Outfitted with 5" front wheels to easily roll over irregular surfaces.

Model:

G07782 (Tyke-Red) G07783 (Tweeners-Green)



Model G07749

CHILDREN'S FOLDING WALKER

- Lightweight anodized aluminum "A" frame design provides greater balance.
- Height-adjustable "grow with you" design.
- Comfortable vinyl hand grips designed for a child's hand.
- Noise-free silencer on footpieces.
- Folds easily for travel.

Model:

G07749 (Child-Brite) G07750 (Toddler-Red)



Model G07751

CHILDREN'S NON-FOLDING WALKER

- 3" wheels maneuver over normal floor obstacles: cracks, carpet, or outside play areas.
- Large rear rubber tips provide long braking action.
- Lightweight, sturdy aluminum frame with proportionally large width to length for wide stance and greater stability
- Height-adjustable.

Model:

G07751 (Tyke)



Model G05163

CHILDREN'S FOREARM CRUTCHES

- Telescoping parts are manufactured with internal bushings and external lock nuts to provide totally silent usage.
- Vinyl-coated, tapered, contoured arm cuffs mold to user for extra comfort
- Heavy-wall, high strength aluminum tubing.
- Black components hide signs of wear and scratches
- Available to fit 'lil tykes from 2 feet tall.

Model:

G05162 (Youth-Brite), G05163 (Child-Brite)



Model G53314-8

CHILDREN'S AXILLARY CRUTCHES

- Push-button adjusted footpiece with fitting guide on crutch footpiece provides accurate fitting.
- 1" adjustment increments on both handle and height provide a precise fit, reducing fatigue.
- Latex free rubber components

Model:

G53314-8

Guardian® Pediatric Walking Aids

Guardian® has a walking aid to meet everyone's needs, even the smallest child. Our Strider Prescription Walkers have been a favorite of therapists. The infinite configurations accommodate the special needs of children with disabilities. All pediatric walking aids feature "grow-with-you" height adjustability and rugged durability.



Guardian® Crutches

Crutches are designed for users that cannot bear weight on one lower limb. Axillary crutches are generally limited to short-term use due to a lower extremity injury. Forearm crutches provide stability and independence for those requiring assistance in ambulation on a long term basis. Guardian crutches are available in a variety of sizes and are easily adjustable to provide a precise fit.

GUARDIAN^{Select} RED•DOT® CRUTCHES

- Clip-lock adjusted handle and push-button adjusted height allows for 10 second adjustments.
- 1" adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue.
- High density latex free rubber cushions and hand grips provide long wear and greater comfort.
- How-to-Fit and Use guides included as part of risk management program.
- 300 lb. weight capacity.

Model:

G90-214-8 (Tall) G91-214-8 (Adult) G92-214-8 (Youth)

QUIK-FIT® CRUTCHES

- Bolt adjusted handle and push-button adjusted footpieces with fitting scale allow for accurate fitting.
- 1" adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue.
- Latex-free rubber components.
- How-to-Fit and Use guides included.
- 300 lb. weight capacity.

Model:

G50314-8 (Tall) G51314-8 (Adult)
G52314-8 (Youth) G53314-8 (Child)

GUARDIAN^{Plus+} BARIATRIC CRUTCHES

- Heavy-duty steel construction provides durability and supports up to 550 lbs.
- Push-button adjusted footpieces with fitting scale allow for accurate, fast fitting.
- 1" adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue.

Model:

G60314B (Tall Adult) G61314B (Adult)

FOREARM CRUTCHES

- Telescoping parts are manufactured with internal bushings and external lock nuts to provide totally silent usage.
- Vinyl-coated, tapered, contoured arm cuffs mold to user for extra comfort
- Heavy-wall, high strength aluminum tubing.
- Black components hide signs of wear and scratches

Model:

G05160 (Tall Adult) G05161 (Adult) G05162 (Youth)



Model G91-214-8



Model G51314-8



Model G61314B



Model G05161



G30341

G30372

G05360

STANDARD CANES

- Broad range of styles and colors.
- Engineered and tested to ensure performance and quality.
- Soft foam handgrip enhances comfort.
- Locking silencer provides rattle free use.
- Solid brass detent button resists wear and breakage.
- Height easily adjusts.
- 250 lb. weight capacity

Quad Canes

Large Base Quad:
Large chrome-plated, 4 point steel base provides maximum stability.
G30340 (Brite)
G30341 (Black)

Small Base Quad:
4 point support for stability and small base for easy clearance.
G30345 (Brite)
G30346 (Black)

Single Point Canes

Standard Handle:
Allows user to retain cane over an arm while opening doors or performing other tasks.
G30370 (Brite)
G30371 (Black)
G30372 (Bronze)

Offset Handle with Wrist Strap:
Offers great stability and a convenient wrist strap.
G05360 (Black)
G05362 (Brite)
G05364 (Bronze)
G05365 (Bronze Slim)

GUARDIAN Select CANES

All of the features of the standard Guardian cane plus:

- 300 lb. weight capacity.
- Handle inserts and reinforced metal inserts at stress points provide great strength and durability.
- Contoured vinyl handgrip is easier to grasp and offers enhanced durability for long-term use.

GUARDIAN Plus+ BARIATRIC QUAD

• Small steel base supports weights up to 500 lbs.
G30343B (Black)

Quad Canes

Select Quadlite:
Molded resin base is strong yet light-weight, offering a slight flex for self-leveling and stability.
G05845S (Brite Adult)
G05846S (Brite Youth)

Select Large Base Quad:
For maximum stability.
G05340S (Brite)

Select Small Base Quad:
Convenient size for use on stairs.
G05345S (Brite)

Single Point Canes

Select Offset handle with Sure Grip:
Offers an anatomically correct design that relieves palm fatigue and provides extra comfort for the long term user.
G05555 (Brite)

Select Standard Handle:
G05370 (Brite)

Select Offset Handle:
G05355 (Brite)
G05356 (Black)



G30343B

G05845S

G05356

G05555

WOODEN CANES

- Attractive, classic style.
- Available in a variety of wood color finishes.

Model:

G05390B (Black) G05390M (Mahogany) G05390W (Walnut)

FOLDING CANE

- Folds for convenient storage and travel.
- Comes in attractive black finish.
- Adjustable height.

Model:

G30375



G05390B

G05390M

G05390W

G30375

Guardian® Canes

For people needing balance and support, Guardian canes are available in a selection of designs and hand grip styles.

WALKING ASSISTANT®

7 FEATURES IN 1 CANE

The Walking Assistant is the newest evolution in walking canes. Never before has a cane offered this much assistance to make every day activities and mobility easier. A grabber picks up items from the floor, a bright LED lights your path, a strong magnet picks up metal items, a non-slip grip allows cane to rest on surfaces, a wrist strap with handy clip for keys make this cane *Intelligent Support.*

Model:

G05940 (Tall 5'10"+)
G05950 (Med 5'5" - 5'10")
G05960 (Short 5'5" under)



Guardian® Folding Walkers

Guardian Easy Care Walkers are the most imitated design, but never duplicated! We build quality features into our walkers for enhanced stability and durability because someone's health is dependent on this product. Compare these quality enhancements:

- Reinforced tube connectors at high-stress areas.
- Solid brass detent buttons that are stronger to support the full weight of the user.
- Double steel cross frame for extra rigidity
- Double side cross braces for extra rigidity
- Screws, that can be retightened, connect side cross braces
- Deburred tubing to prevent cutting through rubber tips
- Tight tolerance between frame and footpieces for rattle-free use
- Durable folding mechanism is designed to last.



Model G30755P

GUARDIAN EASY CARE® WALKERS

Engineered with quality differences that offer maximum stability and longer life.

- Dual folding mechanism is easy to operate with palms or fingers. Ideal for those with limited dexterity.
- Lightweight aluminum frame with steel cross brace provides durability.
- Folds to only 4" for easy transport and storage.
- 300 lb. weight capacity.

Model:

G30755P (Adult) G30756P (Youth)

Models with 5" Wheels:

G30757W (Adult) G30758W (Youth)



Model G07755

GUARDIAN *Select* RED•DOT® WALKERS

- The familiar Red•Dot dual folding mechanism releases simply by depressing the tab with finger, arm or palm.
- Lightweight aluminum frame, with the exclusive Guardian steel "K" cross brace, is chrome plated for added durability.
- Folds to only 4 1/2" for easy transport and storage.
- Preferred by institutions across the nation.
- 350 lb. weight capacity.

Model:

G07755 (Adult) G07756 (Youth)



Model G30760P

ONE BUTTON FOLDING WALKER

- One-button folding offers quick, convenient operation and is ideal for those who are weaker on one side.
- Lightweight aluminum frame supports up to 300 lbs. yet weighs only 5 lbs.
- Folds to 4 1/2" for ease of transport and storage.

Model:

G30760P (Adult) G30761P (Youth)



Model G07767

GUARDIAN *Select* EXTRA-DUTY WALKERS

- Same quality features as Red•Dot with extra width.
- Lightweight – weighs only 7 1/2 pounds.
- 400 lb. weight capacity.
- 19” between hand grips.

Model:

G07767 (Adult) G07768 (Youth)



Model G30754B

GUARDIAN *Plus+* BARIATRIC WALKER

- Reinforced aluminum frame with dual steel cross brace provides durability.
- Extra wide frames with up to 22” between hand grips accommodate larger users.

Model:

G30754B (650 lb. Weight Capacity)



Model G07780

STRIDER® PRESCRIPTION WALKER

- Accommodates most physical requirements with the optional platform attachments.
- Handle adjusts without tools to an infinite number of positions to provide stability and posture. 90° bend is standard.
- Extra wide base provides stability, safety and balance.
- Outfitted with 5" wheels to easily roll over irregular surfaces.

Model:

G07780 (Adult-black) G07781 (Middler-blue)



Model G07770

SIDESTEPPER HEMI WALKER

- Dual grip cane/hemi-walker offers extra assistance when getting in and out of chairs.
- Provides support and stability to those with limited or no dexterity in one arm or hand.
- Lightweight frame is easy to lift and maneuver.
- Wide base provides more stability than a quad cane.

Model:

G07770 (Adult)
G07771 (Youth)





Envoy[®] Rolling Walker Series

For people interested in mobility and support with style, the Guardian[®] Envoy Rolling Walkers are a fashionable alternative. Easy-to-use hand brakes, large padded seats for resting, removable padded backrest for support and comfort. Easy-to-fold mechanisms make transport or storage convenient.



Model G07981T

ENVOY 380 - 3 WHEEL WALKER

- Powder-coated steel frame for long life.
- Tri-wheel design with 8" wheels is easy to maneuver.
- Loop style brakes control speed and easily lock wheels.
- Folds to only 9" for easy transport and storage.
- Tote bag, basket and tray accessories are available.
- Weight capacity 250 lbs.

Model:

- G07981T (titanium color)
- G07913 (tray and basket)
- G07912 (tote bag)



Model G07886G

ENVOY 460 - ECONOMY ROLLING WALKER

- Weighing 15.5 lbs., the lightweight, anodized aluminum frame with 6" wheels loads into vehicles with ease.
- Loop style brakes control speed and easily lock wheels while seated.
- Padded seat and backrest provide comfort.
- Store personal items securely in basket underneath seat.
- Weight capacity 275 lbs.

Model:

- G07886G (Green) G 07886B (Blue) G07886RD (Red)



Model G07887R

ENVOY 480 - ROLLING WALKER

- Lightweight, anodized aluminum frame has height adjustable handgrips.
- Loop style brakes control speed and easily lock wheels while seated.
- Deluxe padded seat and curved padded backrest provide enhanced support and comfort.
- Large 8" wheels are great for outdoor terrain.
- Large basket and seat are standard.
- Weight capacity 300 lbs.

Models:

- G07887B (Blue) G07887R (Rose)



Model G07889T



Model G07890TB

**ENVOY 480HD
HEAVY DUTY ROLLING WALKER**

- Extra wide design, with 21.5" between handles, accommodates larger users.
- Anodized aluminum frame has height adjustable handgrips.
- Loop style brakes control speed and easily lock wheels while seated.
- Deluxe padded seat and curved padded backrest provide enhanced support and comfort.
- Large 8" wheels are great for outdoor terrain.
- Basket and seat are standard.
- Weight capacity 375 lbs.

Model:

G07889T (titanium color)

**GUARDIAN Plus+ ENVOY 485B
BARIATRIC ROLLING WALKER**

- Reinforced steel frame offers greater durability and support.
- Extra width between handles (22.25") accommodates larger users.
- Handles point forward to promote better posture and walking gait.
- Loop style brakes control speed and easily lock wheels while seated.
- Curved padded backrest provides comfort and support.
- Large 8" wheels are great for outdoor terrain.
- Basket and seat are standard.
- Weight capacity 500 lbs.

Model:

G07890TB (titanium color)





FRONT WHEEL ATTACHMENTS

- Front wheel footpieces not only allow you to easily maneuver over surfaces which provide difficulty for walker users, but they also provide a more natural walking stride.
- Available in a choice of 3" and 5" inline wheels or 360° casters which allow you to change directions with ease.

Note: Wheel attachments are to be used on the front legs of the walker only.



5" Wheels

Model:
G07722-8 (1" dia.)
G07722G (1" retail)
G07723 (7/8" dia.)
G07722-8B (bariatric)

5" Swivel Caster

Model: G07811
(1" dia. retail)

3" Swivel Caster

Model: G07825-8
(1" dia. retail)

3" Wheels

Model:
G07731-8 (1" dia.)
G07731G (1" retail)
G07729 (7/8" dia.)

Guardian® Walker Accessories

Guardian offers a full line of walker accessories, replacement parts and attachments to enhance the life, flexibility and enjoyment of your walker. Refer to the accessory selection guide for accessories available to fit your walker.

REAR BRAKE ATTACHMENTS

- Provides safe, gradual rear leg braking action and smooth forward movement when used with front wheel attachments.
- Brake firmly engages the rubber tip on the floor when weight is placed on the walker.
- Choice of easy rolling wheel brakes or Glidebrakes.



Glide Caps

Models:
G07904 (retail)
G07900CS
(bulk 24 pr)

Glidebrakes

Model: G07732G
(1" dia. retail)

Wheel Brakes

Model: G07832-8
(1" dia. retail)

GLIDE CAPS

- Included with all wheel footpieces.
- Allows walkers to easily and quietly glide across most non-abrasive surfaces.
- Fits tips with 1 1/2" O.D.

FOOTPIECE EXTENSIONS

- Extend the height of the walker by up to 7".

Models:

G07708 with 3" wheels, 1" dia.
G07724 with 5" wheels, 7/8" dia.
G07725 with 5" wheels, 1" dia.
G07712 with rubber tips, 1" dia.

WALKER ENHANCEMENT ACCESSORIES



Folding Flip Tray
Model G07850

Walker Seat
Model G07872

Basket Organizer
Model G07878

Walker Basket
Model G07715
G07713 (for one-button)

Front Pouch
Model G07741

Platform Attachment
Model G07702 (Adult)
Model G07703 (Youth)



Model G00018

CRUTCH CUSHIONS

- Designed with thick cushioning for extra comfort.
- Durable Latex-free rubber.
- Fits most aluminum and wooden crutches.
- Attractive retail packaging.

Model:

G00013 (Tan) G00018 (Gray)



Model G00017

CRUTCH HANDGRIPS

- Durable Latex-free rubber provides a comfortable sturdy grip.
- Fits most aluminum and wood crutches with bolt mountings.
- Attractive retail packaging.

Model:

G00016 (Tan) G00017 (Gray) G00015 (Red Dot)

CRUTCH TIPS

- Large surface area provides superior tracking.
- Heavy-duty latex-free rubber tip flexes for cushioning and shock absorption.
- Reinforced with metal insert for extended wear.
- Fits aluminum and wood crutches.

Model:

G00500 (Tan 7/8") G00502 (Gray 7/8")
 G00840 (Tan 7/8") G00841 (Black 7/8") G00842 (Gray 7/8")



Model G00502

QUAD CANE TIPS

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.
- Special tread for skid resistance and secure traction.
- Four tips per package.

Model:

G02161 (Black 1/2") G02162 (Gray 1/2")
 G02171 (Black 5/8") G02172 (Gray 5/8")



Model G02162

CANE TIPS

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.
- Two tips per package.
- Attractive retail packaging.

Model:

G30561-PC (Black 5/8") G30571-PC (Black 3/4")
 G0572-PC (Gray 3/4") G30573-PC (Brown 3/4")
 G30581-PC (Black 7/8")



Model G30571-PC

UTILITY TIPS

- Heavy-duty rubber tip is reinforced with metal insert for longer wear.
- Fit walkers, commodes and transfer benches.
- Four tips per package.

Model:

G01182 (Gray 3/4") G01202 (Gray 1") G01212 (Gray 1 1/8")



Model G01212

Guardian[®] Rubber Accessories

You should replace the tips on your walker, cane or crutches whenever they appear to show wear to maintain a safe gripping surface. Remove the existing tip to use the sizing guide on the replacement tip cartons.

General Guidelines for Replacement Guardian Tips

Model	Tip Size
Aluminum Canes with 7/8" Handle	3/4"
Aluminum Canes with 3/4" Handle	5/8"
Wood Canes with 7/8" Handle	7/8"
Small Base Quad Canes	1/2"
Large Base Quad Canes	5/8"
Axillary Crutches	7/8"
Forearm Crutches	7/8"
Walkers	1 1/8"

ROLLING WALKERS

Model	Height Range of Handles	Seat Height	Width Between Handles	Overall Width	Overall Length	Weight	Max. Weight Capacity
Envoy 380	31 $\frac{1}{2}$ "-36 $\frac{3}{4}$ " (81-93 cm)	N/A	14 $\frac{1}{2}$ " (37 cm)	25 $\frac{1}{2}$ " (65 cm)	24 $\frac{1}{2}$ " (62 cm)	14 $\frac{1}{2}$ lbs. (6.6 kg)	250 lbs. (113 kg)
Envoy 460	34"-38" (86-96 cm)	21 $\frac{3}{4}$ " (55 cm)	18 $\frac{1}{2}$ " (47 cm)	24" (61 cm)	23" (58 cm)	15 $\frac{1}{2}$ lbs. (7 kg)	275 lbs. (125 kg)
Envoy 480	30 $\frac{3}{4}$ "-36" (78-91 cm)	20" (51 cm)	18" (46 cm)	24 $\frac{1}{2}$ " (62 cm)	30 $\frac{1}{2}$ " (77 cm)	20 $\frac{1}{2}$ lbs. (9.3 kg)	300 lbs. (136 kg)
Envoy 480HD	31"-36 $\frac{1}{2}$ " (79-93 cm)	20" (51 cm)	21 $\frac{1}{2}$ " (55 cm)	28 $\frac{1}{2}$ " (72 cm)	30 $\frac{1}{2}$ " (77 cm)	22 $\frac{1}{2}$ lbs. (10 kg)	375 lbs. (170 kg)
Envoy 485B	32 $\frac{1}{4}$ " - 34 $\frac{1}{4}$ " (82-87 cm)	22" (56 cm)	22 $\frac{1}{4}$ " (56 cm)	26 $\frac{1}{4}$ " (67 cm)	30 $\frac{1}{2}$ " (77 cm)	25 lbs. (11.5 kg)	500 lbs. (227 kg)

FOLDING WALKERS

MODEL	SIZE	HEIGHT ADJUSTS	UNIT WEIGHT	BASE WIDTH	INSIDE GRIP WIDTH	BASE DEPTH, OPEN	BASE DEPTH, FOLDED	MAX. WEIGHT CAPACITY
Easy Care	Adult	32 $\frac{1}{2}$ "-39 $\frac{1}{2}$ " (82-100 cm)	6 $\frac{1}{2}$ lbs. (3 kg)	24" (61 cm)	16 $\frac{1}{2}$ " (42 cm)	17" (43 cm)	4" (10 cm)	300 lbs. (136 kg)
	Youth	27 $\frac{1}{4}$ "-34 $\frac{1}{4}$ " (69-87 cm)	6 lbs. (2.7 kg)	24" (61 cm)	16 $\frac{1}{2}$ " (42 cm)	16 $\frac{1}{4}$ " (41 cm)	4" (10 cm)	300 lbs. (136 kg)
Red Dot	Adult	32 $\frac{1}{2}$ "- 39 $\frac{1}{2}$ " (82-100 cm)	6 $\frac{1}{2}$ lbs. (3 kg)	24" (61 cm)	16 $\frac{1}{2}$ " (42 cm)	17" (43 cm)	4 $\frac{1}{2}$ " (11 cm)	350 lbs. (159 kg)
	Youth	27 $\frac{1}{4}$ "- 34 $\frac{1}{4}$ " (69-87 cm)	6 lbs. (2.7 kg)	24" (61 cm)	16 $\frac{1}{2}$ " (42 cm)	16 $\frac{1}{4}$ " (41 cm)	4 $\frac{1}{2}$ " (11 cm)	350 lbs. (159 kg)
One Button	Adult	31 $\frac{1}{2}$ "- 38 $\frac{1}{2}$ " (80-98 cm)	5 $\frac{3}{4}$ lbs. (2.6 kg)	24" (61 cm)	17 $\frac{1}{4}$ " (44 cm)	18 $\frac{1}{2}$ " (47cm)	4 $\frac{1}{2}$ " (11 cm)	300 lbs. (136 kg)
	Youth	27 $\frac{1}{4}$ "- 34 $\frac{1}{4}$ " (69-87 cm)	5 lbs. (2.3 kg)	23" (58 cm)	17 $\frac{1}{4}$ " (44 cm)	16 $\frac{1}{2}$ " (42 cm)	4 $\frac{1}{2}$ " (11 cm)	300 lbs. (136 kg)
Strider	Adult	28"-39" (71-99 cm)	9 lbs. (4.1 kg)	26 $\frac{1}{2}$ " (67 cm)	adjustable	27 $\frac{1}{2}$ " (70 cm)	8" (20 cm)	250 lbs. (113 kg)
	Middler	24"- 36" (61-91 cm)	8 lbs. (3.6 kg)	28" (71 cm)	adjustable	25 $\frac{1}{2}$ " (65 cm)	8" (20 cm)	200 lbs. (91 kg)
Extra-Duty	Adult	32 $\frac{1}{4}$ "- 39 $\frac{1}{4}$ " (82-99 cm)	7 $\frac{1}{2}$ lbs. (3.4 kg)	26 $\frac{1}{4}$ " (67 cm)	19" (48 cm)	17 $\frac{1}{2}$ " (44 cm)	4 $\frac{1}{2}$ " (11 cm)	400 lbs. (181 kg)
	Youth	27 $\frac{1}{4}$ "- 34 $\frac{1}{4}$ " (69-87 cm)	7 $\frac{1}{2}$ lbs. (3.4 kg)	26 $\frac{1}{4}$ " (67 cm)	19" (48 cm)	17" (43 cm)	4 $\frac{1}{2}$ " (11 cm)	400 lbs. (181 kg)
Bariatric	Adult	32"- 39" (81-99 cm)	8 $\frac{1}{2}$ lbs. (3.8 kg)	28" (71 cm)	22" (56 cm)	18" (46 cm)	4 $\frac{1}{2}$ " (11 cm)	650 lbs. (295 kg)
Sidestepper	Adult	32"- 36" (81-91 cm)	3 $\frac{3}{4}$ lbs. (1.7 kg)	17 $\frac{1}{4}$ " (44 cm)	N/A	15 $\frac{3}{4}$ " (40 cm)	3" (7.6 cm)	250 lbs. (113 kg)
	Youth	28 $\frac{1}{2}$ "- 32" (72-81 cm)	3 $\frac{3}{4}$ lbs. (1.7 kg)	17 $\frac{1}{4}$ " (44 cm)	N/A	15 $\frac{3}{4}$ " (40 cm)	3" (7.6 cm)	250 lbs. (113 kg)

SINGLE POINT & QUAD CANES

MODEL	SIZE	APPROXIMATE USER HEIGHT	HEIGHT ADJUSTS	UNIT WEIGHT	BASE/SHAFT WIDE	BASE/SHAFT NARROW	REPLACEMENT TIP	MAXIMUM WEIGHT CAPACITY (STD/SELECT)
Quad Canes	Adult	4'10" - 6'4"	29" - 38"	2.5 lbs.	8" x 12"	6" x 8"	G02162 (small base)	250 lbs. / 300 lbs.
		147 - 193 cm	74 - 97 cm	1.1 kg	20 - 30 cm	15 - 20 cm	G02172 (large base)	113 kg / 136 kg
	Youth	3'10" - 5'4"	25" - 34"	2 lbs.	8" x 12"	6" x 8"	G02162 (small base)	250 lbs. / 300 lbs.
		117 - 163 cm	63 - 86 cm	.9 kg	20 - 30 cm	15 - 20 cm	G02172 (large base)	113 kg / 136 kg
Single Point Canes	Adult	4'10" - 6'4"	31" - 40"	.75 lb.	7/8" dia.	3/4" dia.	G30571-PC/72 (7/8")	250 lbs. / 300 lbs.
		147 - 193 cm	79 - 100 cm	.3 kg	2.2 cm dia.	1.9 cm dia.	G30561-PC (3/4")	113 kg / 136 kg
Wooden Canes	Adult	Up to 6'2"	Up to 36"	.75 lbs.	7/8" dia.	N/A	G30581-PC	300 lbs.
		188 cm	91 cm	.3 kg	2.2 cm dia.			136 kg
Folding Cane	Adult	5' - 6'	33" - 37"	.75 lbs.	3/4" dia.	N/A	30561	200 lbs.
		152 - 183 cm	84 - 94 cm	.3 kg	1.9 cm dia.			91 kg

AXILLARY CRUTCHES

MODEL	SIZE	APPROXIMATE USER HEIGHT	HEIGHT ADJUSTS	MAX. WEIGHT CAPACITY
Red•Dot/ Quik-fit	Tall Adult	5'10" - 6'6"	52" - 60"	300 lbs.
		178 - 198 cm	132 - 152 cm	136 kg
Red•Dot/ Quik-fit	Adult	5'1" - 5'9"	44" - 52"	300 lbs.
		155 - 175 cm	112 - 132 cm	136 kg
Red•Dot/ Quik-fit	Youth	4'6" - 5'2"	37" - 45"	300 lbs.
		132 - 157 cm	94 - 114 cm	136 kg
Quik-fit	Child	4'0" - 4'6"	31 $\frac{1}{2}$ " - 37 $\frac{1}{2}$ "	200 lbs.
		122 - 137 cm	80 - 95 cm	91 kg
Bariatric	Tall Adult	5'10" - 6'6"	52" - 60"	550 lbs.
		178 - 198 cm	132 - 152 cm	249 kg
Bariatric	Adult	5'1" - 5'9"	44" - 52"	550 lbs.
		155 - 175 cm	112 - 132 cm	249 kg

FOREARM CRUTCHES

SIZE	APPROXIMATE USER HEIGHT	HEIGHT ADJUSTS	UNIT WEIGHT	ADJUSTMENT CUFF TO HANDLE	MAX. WEIGHT CAPACITY
Tall Adult	5'10" - 6'6"	35" - 41"	4 $\frac{3}{4}$ lbs.	9 $\frac{1}{2}$ " - 11"	300 lbs.
	178 - 198 cm	89 - 104 cm	2 kg	24 - 28 cm	136 kg
Adult	5' - 6'2"	30" - 36"	4 $\frac{3}{4}$ lbs.	9 $\frac{1}{2}$ " - 11"	300 lbs.
	152 - 183 cm	76 - 91cm	2 kg	24 - 28 cm	136 kg
Youth	4'2" - 5'2"	25" - 31"	4 lbs.	8 $\frac{1}{2}$ " - 10"	300 lbs.
	127 - 158 cm	64 - 79 cm	1.8 kg	22 - 25 cm	136 kg
Child	3'2" - 4'6"	19" - 27"	2 $\frac{1}{2}$ lbs.	7 $\frac{1}{4}$ " - 8 $\frac{1}{2}$ "	200 lbs.
	97 - 137 cm	48 - 69 cm	1.1 kg	24 - 28 cm	91 kg

Specifications

PEDIATRIC WALKERS

MODEL	SIZE	APPROXIMATE USER HEIGHT	HEIGHT ADJUSTS	UNIT WEIGHT	BASE WIDTH	INSIDE GRIP WIDTH	BASE DEPTH OPEN	BASE DEPTH FOLDED	MAX. WEIGHT CAPACITY
Strider	Tweeners	3' - 5'2"	24"-36"	7 1/4 lbs.	24 1/4"	adjustable	25 1/2"	9 1/2"	175 lbs.
		91 - 157 cm	61-91 cm	3.3 kg	62 cm		65 cm	24 cm	80 kg
	Tyke	2'2" - 3'6"	14"-19.5"	5 1/2 lbs.	22"	adjustable	20"	9 1/2"	90 lbs.
		66 - 107 cm	36-50 cm	2.5 kg	56 cm		51 cm	24 cm	41 kg
Children's Folding	Child	4' - 4'8"	24 1/2" - 28 1/2"	5 lbs.	22"	12 3/4"	16 3/4"	8"	200 lbs.
		122 - 142 cm	62 - 72cm	2.3 kg	56 cm	32 cm	41 cm	20 cm	91 kg
	Toddler	3'4" - 4'	20 1/2" - 24 1/2"	5 lbs.	18 1/2"	10 3/4"	13 1/4"	8"	175 lbs.
		102 - 122 cm	52 - 62cm	2.3 kg	47 cm	27 cm	34 cm	20 cm	80 kg
Tyke Non-Folding	Tyke	2'8" - 3'6"	15 1/2" - 20 1/2"	5 lbs.	17"	13"	15 1/2"	N/A	80 lbs.
		81 - 107 cm	39 - 52cm	2.3 kg	43cm	33cm	39cm		36 kg

WALKER ACCESSORY SELECTION GUIDE

Model	Description	Size	Red•Dot				Easy Care				Single	Strider				Pediatric		Bari
			G07755	G07756	G07767	G07768	G30755P	G30757W	G30756P	G30758W	G30760P/G30761P	G07780	G07781	G07783	G07782	G07749/G07750	G07751	G30754B
G07731-8	3" (7.62 cm) Fixed Wheels - pair	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07731G	3" (7.62 cm) Fixed Wheels - pair Retail Pkg	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07729	3" (7.62 cm) Fixed Wheels - pair	7/8" (2.22 cm)														•		
G07708	3" (7.62 cm) Fixed Wheels Extended - set	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07722-8	5" (12.7 cm) Fixed Wheels - pair	1" (2.54 cm)	•	•	•	•	•	Inc	Inc	•	Inc							
G07722G	5" (12.7 cm) Fixed Wheels - pair Retail Pkg	1" (2.54 cm)	•	•	•	•	•	Inc	Inc	•	Inc							
G07722-8B	5" (12.7 cm) Fixed Wheels bariatrics - pair	1" (2.54 cm)															•	
G07723	5" (12.7 cm) Fixed Wheels - pair	7/8" (2.22 cm)										Inc	Inc					
G07724	5" (12.7 cm) Fixed Wheels Extended - set	7/8" (2.22 cm)														•		
G07725	5" (12.7 cm) Fixed Wheels Extended - set	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07825-8	3" (7.62 cm) Swivel Wheels - pair	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07820-8	3" (7.62 cm) Swivel Wheels w/Glidebrakes	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07811	5" (12.7 cm) Swivel Wheels - pair	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07732G	GlideBrakes - pair	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07832-8	Wheel Brakes - pair	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07900CS	Glide Caps - Bulk 24 prs/case	1 1/2" (3.81 cm)	•	•	•	•	•	•	•	•	•	•	•	•	•		•	
G07904	Glide Caps - 2 pairs/retail pkg	1 1/2" (3.81 cm)	•	•	•	•	•	•	•	•	•	•	•	•	•		•	
G07712	Extension Footpieces - set of 4	1" (2.54 cm)	•	•	•	•	•	•	•	•	•							
G07796	Platform Attachment	1" (2.54 cm)									•							
G07797	Platform Attachment	7/8" (2.22 cm)										•	•					
G07798	Platform Attachment	3/4" (1.9 cm)														•		
G07702	Platform Attachment, Adult	1" (2.54 cm)	•	•			•	•										
G07703	Platform Attachment, Youth	1" (2.54 cm)	•	•					•	•								
G07741	Carrying Pouch - Front		•	•	•	•	•	•	•	•	•	•	•					
G07850	Fliptray		•	•			•	•	•	•								
G07715	Carrying Basket		•	•			•	•	•	•		•	•					
G07713	Carrying Basket - One Button										•							
G07872	Walker Seat		•	•			•	•	•	•								

Guardian® Bath Safety

Guardian® offers a complete line of bath safety products to help meet all of your customers' needs. With a commitment to delivering the highest level of quality, each Guardian product is engineered for safety and durability. These bath safety products offer the ideal choice for comfort and convenience.

GUARDIAN EASY CARE® SHOWER CHAIR/STOOL

A lightweight and economical choice for those who require a basic chair with height adjustability and an extra wide seating surface.

- Sturdy, lightweight, anodized aluminum frame resists rusting.
- Legs adjust in 1" increments with just a push of a button to ideal height between 13" and 20".
- Available with a back for extra support and leverage, or backless for easier access in bathing.
- Textured, blow-molded polyethylene wide track seat and backrest with built-in hand-holds provide extra safety.
- Rubber feet grip even wet surfaces.
- Available in attractive retail packaging.
- Weight capacity 250 lbs.

Model:

G30400-4 (Chair assembled-Gray)

G30401-4 (Stool assembled-Gray)

G30402 (Chair unassembled-Retail-White)

G30403 (Stool unassembled-Retail-White)



Model G30402



Model G0403



Model G30406B

GUARDIAN *Plus+* BARIATRIC SHOWER CHAIR

- Frame is reinforced to support weights up to 500 lbs.
- Extra wide seating surface accommodates larger users.
- Seat height adjusts from 15"-19" for a comfortable fit.
- Rubber tips grip bathtub surface.

Model:

G30406B



Model G98027

GUARDIAN *Select* PADDED SHOWER CHAIR

- Extra thick padding on seat and backrest provides additional comfort for long-term use.
- Extra deep 16 1/2" seat for greater stability and comfort.
- Legs adjust in 1" increments to ideal height between 18" and 22".
- Rubber tips grip bathtub surface.
- 300 lb. weight capacity.

Model:

G98027



Model G98308A

NON-PADDED TRANSFER BENCH

- Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage.
- Lightweight, anodized aluminum frame is height adjustable in 1/2" increments.
- Backrest attaches for right or left facing tubs easily without tools.
- Side arm and seat back handle provide extra stability and leverage when transferring.
- Weight capacity 300 lbs.

Model:
G98308A (Assembled)



Model G98338F

PADDED TRANSFER BENCH

- Vinyl covered foam backrest and seating surface offer a highly resilient cushioned surface for comfort.
- Lightweight, folding aluminum frame is height adjustable in 1/2" increments.
- Suction cup footpieces secure the bench to the tub floor.
- Backrest adjusts for right or left facing tubs.
- Side arm and seat handle provide extra stability and leverage when transferring.
- A cost effective solution with the added feature of padding.
- Weight capacity 300 lbs.

Model:
G98338F (Assembled)



Model G98309

GUARDIAN^{Select} NON-PADDED TRANSFER BENCH

- Suction cup footpieces secure the bench to tub floor
- Lightweight, anodized aluminum frame is height adjustable in 1/2" increments to accommodate most tubs.
- Backrest adjusts without tools for right or left facing tubs.
- Built-in soap dish and hand-held shower slot for user convenience.
- Weight capacity 300 lbs.

Model:
G98309 (Assembled)
G98309G (Unassembled-Retail)



Model G98013A

PADDED TRANSFER BENCH WITH COMMODE OPENING

- A commode opening affords the user greater comfort and better access for cleaning.
- Thick, high density foam cushions and backrest for comfortable seating that won't bottom out.
- Suction cup footpieces secure the bench to tub floor.
- Lightweight, folding aluminum frame is height adjustable in 1/2" increments to accommodate most tubs.
- Weight capacity 300 lbs.

Model:
G98013A (Assembled)



Model G98310B

GUARDIAN^{Plus+} BARIATRIC TRANSFER BENCH

- 500 lb. weight capacity accommodates larger users.
- Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage.
- Aluminum frame is height adjustable in 1/2" increments.
- Backrest adjusts without tools for right or left facing tubs.
- Side arm and seat handle provide extra stability and leverage when transferring.
- Safely provides assistance to those who have difficulty stepping over bathtub walls.
- Backrest and armrest installed with detent buttons for tool free assembly.

Model:
G98310B (Assembled)

Guardian® Bathing Products

ERGONOMIC HAND HELD SHOWER

- Easily adapts to existing shower arm, delivering standard or pulsating spray.
- Extra long 86" (218 cm) hose for use while seated.
- Built-in water saver. Uses less than 2.5 gallons (9.5 L) per minute.
- On/Off control and massage are conveniently positioned at your fingertips.
- Diverter valve allows you to also use existing shower head.

Model:

G98035 (Hand held shower with diverter valve.)



Model G98035

DELUXE TUB GRAB BAR

- Clamp design, including stainless steel components, stays secure after repeated use.
- Patented step through design provides an unobstructed entry/exit path allowing better hand positioning.
- Bi-level handgrip for maximum flexibility and safety.
- Tool free installation and large rubber-lined pads won't mar tub surface.
- Adjusts to fit tubs 2.75"-6.25" (7-16 cm) thick.

Model:

G98006



Model G98006

BATHTUB SAFETY RAILS

- An economical choice for a safe grip around the tub.
- White powder-coated steel with foam-lined clamps that won't mar tub surface.
- Mounting bracket adjusts 3.5" to 6" (9-16 cm) to fit most standard tubs.

Model:

G98080 (8"H x 11"W - 20 x 27 cm)



Model G98080

Guardian® Grab Bars



Model G81060KE

G81000 SERIES – CHROME PLATED WALL GRAB BARS

- Adjustable oval flange rotates to best position for secure mounting.
- 1" (2.5 cm) diameter; 18 gauge chrome plated steel.
- Knurled finish provides a sure grasp even when wet.
- 12"-36" comes in attractive retail packaging.

Model	Length
G81020KE	12" (30 cm)
G81030KE	16" (41 cm)
G81040KE	18" (46 cm)
G81050KE	24" (61 cm)
G81060KE	32" (81 cm)
G81070KE	36" (91 cm)




Model G87560

G87500 SERIES – WHITE FINISH WALL GRAB BAR

- White painted steel matches most decors.
- 1" (2.5 cm) diameter; 18 gauge steel with 2.5" (6 cm) flange.
- White screw covers included.
- 12"-36" comes in attractive retail packaging.

Model	Length
G87520	12" (30 cm)
G87530	16" (41 cm)
G87540	18" (46 cm)
G87550	24" (61 cm)
G87560	32" (81 cm)
G87570	36" (91 cm)



Model G98120

SAFE-T-POLE

- Floor-to-ceiling pole that safely assists standing, sitting, climbing stairs, or transferring.
- No tools required for installation. Screw mechanism with padded bracket fits securely without marring surface.
- White powder-coated, heavy gauge steel is attractive and durable.

Model: G98120

Diameter:	1.5" (4 cm)
Height Adjusts:	7' 8.25" - 8' 2.5" (234-260 cm)
Weight:	14 lbs (6.4 kg)
Color:	White



Guardian® Home Assistance Products

SITZ BATH

- To ease discomfort from hemorrhoids and other disorders.
- Contour design for comfort during use.
- Complete with vinyl bag and tubing for convenient home use.
- Easy to follow instructions included.
- Plastic bowl conveniently fits over most toilets.
- Easy to clean with household cleaners.

Model: G903-PC



Model G903-PC

MALE URINAL

- Ideal for use when it's difficult to get out of bed.
- Snap-on cap prevents spills and odor.
- Easy to measure — imprinted in ounces and cc's.
- Easy to clean with household cleaners.

Model: G900



Model G900

FEMALE URINAL

- Contoured design for a complete and secure seal during use.
- Designed to prevent accidental spills.
- Easy to use in several positions.
- Easy to measure — imprinted in ounces and cc's.

Model: G901-PC



Model G901-PC

HOSPITAL STYLE BEDPAN

- Ideal for use when it's difficult to get out of bed.
- Contoured "hospital" style for comfortable use.
- Molded from polyethylene plastic for light weight and odor resistance.
- Easy to clean with household cleaners.

Model: G904-PC



Model G904-PC

FRACTURE STYLE BEDPAN

- Fracture bedpans are for use with immobile or fracture patients.
- Convenient for use with the assistance of a caregiver.
- Contour shape for extra patient comfort
- May be used as a female urinal.
- Molded from polyethylene plastic for light weight and odor-resistance.
- Easy to clean with household cleaners.

Model: G905



Model G905



Model G30260



Model G30270A



Model G30250



Model G30300

GUARDIAN Select LOCKING RAISED TOILET SEATS

- Raises seat height 5" above the toilet.
- Easy-to-use locking system requires no tools and fits most toilets.
- Ergonomically correct seating surface for extra comfort.
- Available with wraparound foam padded armrests for a sure grip, comfort and support.
- Lightweight polypropylene resin is easy to clean.
- Weight capacity 250 lbs (113 kg).

Model:

G30270A (with padded arms)

G30260 (without arms)

Width Overall: 20.5" (52 cm)

Length Overall: 18" (46 cm)

Opening: 8" wide x 9.5" oblong (20 x 24 cm)

Height from seat surface:

to seat 5" (13 cm)

to arms 11" (28 cm)

Weight: 6 lbs (3 kg)

Weight capacity: 250 lbs (113kg)

ECONOMY RAISED TOILET SEAT

- Universal design is made of easy-to-clean polypropylene resin.
- Orthopedically designed to provide user with stability and comfort.
- Shaped to fit most toilets.

Model: G30250

Width Overall: 14.5" (37 cm)

Length Overall: 15.5" (39 cm)

Height from Seat:

front 5" (13 cm)

back 6" (15 cm)

Weight: 4 lbs (1.6 kg)

Weight capacity: 250 lbs (113kg)

TOILET SAFETY FRAME

- Handles are adjustable and rotate back to allow a wide range of comfortable and secure positions.
- Closed cell foam armrests provide a comfortable, sure grip.
- Adjustable height 26" - 31" (66-74 cm) accommodates standard or elevated toilet seats.
- Easy-to-clean aluminum frame mounts securely onto bowl with adjustable bracket 18"-24" (46-61 cm).
- Available in attractive retail packaging.

Model:

G30300 (bulk pack 2 per carton)

G30300-1 (retail package)

Height Adjustment

armrest to floor: 26"-31" (66-74 cm)

Width Adjustment

between armrests: 8"-24" (46-61 cm)

Overall Depth:

arm: 16.75" (43 cm)

at leg 12" (30 cm)

Weight: 10 lbs (4.5 kg)

Weight capacity: 250 lbs (114 kg)

Guardian®
Toilet
Assist
Aids



Guardian Easy Care® Commodes



Model G30213-4F

3 IN 1 STEEL COMMODE

- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Roomy 19" between arms is easy to use.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Steel frame is rust resistant and chrome plated finish offers greater durability.
- Includes pail, lid and splash guard
- Height adjusts 16.5" - 22.5"
- 350 lb weight capacity.
- Lifetime limited warranty.

Models:

- G30213-4F (Folding Frame - Painted)
- G30213-4 (Painted)
- G30211-4 (Chrome Plated)



Model G30210-4F

3 IN 1 ALUMINUM COMMODE

- Aluminum frame is lightweight and easy to transport.
- Folding frame design assembles in seconds without tools and reduces storage space by 35% while making delivery easy. Allows user to refold for transport or storage.
- 3 in 1 design allows for use as a bedside commode, or raised toilet seat and safety frame over a standard toilet.
- Full 19" between arms is easy to use.
- Wide base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Includes pail, lid and splash guard
- Height adjusts 17.5" - 21.5"
- 250 lb. weight capacity.
- Lifetime limited warranty.

Model:

- G30210-4F (Folding Frame - Aluminum)



Model G30214-2

HEAVY DUTY 3 IN 1 COMMODE

- Extra wide frame is easy to use with a full 21.12" between arms.
- Large base depth provides more stability.
- Clip-on seat can be easily removed for cleaning.
- Height adjusts 18" - 22"
- Weight capacity of 400 lbs.
- Lifetime limited warranty.

Model: G30214-2

Guardian® Commodes



Model G30216B



Model G98202



Model G98204

GUARDIAN *Plus+* BARIATRIC COMMODE

- Extra wide frame with 23.5" between arms accommodates larger users.
- Reinforced steel frame offers great durability.
- Includes pail, lid and splash guard
- Height adjusts 16.5" - 23.5".
- 500 lb. weight capacity.
- Lifetime limited warranty.
- 2 per case.

Model: G30216B

NON-PADDED DROP-ARM COMMODE


- Designed to accommodate lateral transfers.
- Cantilevered arm design drops with just a push of a button, allowing unrestricted left-to-right leg movement for easy transfer.
- Hard, chrome plated steel frame is rust resistant and offers greater durability.
- A cost effective solution with a non-padded round seat.
- Foam handgrips won't slip, even if wet.
- Backrest can be removed for use as an over-the-toilet safety frame.
- Pail, lid and splash guard included.
- Height adjusts 18.25" - 23.25".
- Weight capacity 300 lbs.
- Lifetime limited warranty.

Model: G98202

GUARDIAN *Select* PADDED DROP-ARM COMMODE

- Same features as 98202 with thick, resilient foam padded seat and back for increased comfort.
- Weight capacity 350 lbs.
- Lifetime limited warranty.

Model: G98204



Specifications

SHOWER CHAIRS

Model	Color	Backrest	Padded Cushion	Commode Seat	Assembly Required	Seat Height Range	Seat Depth	Seat Width	Max. Wt. Capacity	Limited Warranty
G30400 Series	White/ Gray	G30400, G30402 Yes G30401, G30403 No	No	No	G30400, G30401 TF/No G30402, G30403 Yes	13"-20" 33-51 cm	12" 30 cm	19" 48 cm	250 lbs 113 kg	Lifetime
G98207	Beige	Yes	Yes	No	Tool-free	18"-22" (46-56 cm)	16.75" (43 cm)	17" (43 cm)	300 lbs (136 kg)	Lifetime
G30406B	White	Yes	No	No	Partial	15"-19" (38-48 cm)	15" (38 cm)	20" (51 cm)	500 lbs (227 kg)	Lifetime

TRANSFER BENCHES

Models	Color	Backrest	Padded Cushion	Commode Seat	Assembly Required	Seat Height Range	Seat Depth	Overall Depth	Seat Width	Overall Width	Max. Wt. Capacity	Limited Warranty
G98308A	Grey	Yes	No	No	Tool-free	17"-22" (43-56 cm)	16.25" (41 cm)	23"-24" (58-61 cm)	26.5" (67 cm)	30.5" (77 cm)	300 lbs (136 kg)	Lifetime
G98309	White	Yes	No	No	Tool-free	17"-21.5" (43-55 cm)	14.5" (37 cm)	21" (53 cm)	31.5" (80 cm)	31.5" (80 cm)	300 lbs (136 kg)	Lifetime
G98309G	White	Yes	No	No	Yes	17"-21.5" (43-55 cm)	14.5" (37 cm)	21" (53 cm)	31.5" (80 cm)	31.5" (80 cm)	300 lbs (136 kg)	Lifetime
G98338F	Beige	Yes	Yes	No	Tool-free	17"-22" (43-56 cm)	16.5" (42 cm)	21" (53 cm)	26.5" (67 cm)	30.5" (77 cm)	300 lbs (136 kg)	Lifetime
G98013A	Beige	Yes	Yes	Yes	Tool-free	17"-22" (43-56 cm)	16.5" (42 cm)	21.5" (55 cm)	25.5" (65 cm)	29" (74 cm)	300 lbs (136 kg)	Lifetime
G98310B	Grey	Yes	No	No	Tool-free	15.5"-20.5" (39-52 cm)	16.25" (41 cm)	23" (58 cm)	26.5" (67 cm)	30.5" (77 cm)	500 lbs (227 kg)	Lifetime

COMMODES

Models	Frame Material	Backrest	Assembly Required	Seat Height Range	Seat Depth	Width Between Arms	Overall Width	Unit Weight	Max. Wt. Capacity	Limited Warranty
G30210-4F	Aluminum	Yes	Tool-free	17.5"-21.5" (44-54 cm)	16" (41 cm)	19" (48 cm)	21.25" (54 cm)	9 lbs (4 kg)	250 lbs (113 kg)	Lifetime
G30211-4	Chromed Steel	Yes	No	16.5"-22.5" (42-57 cm)	16" (41 cm)	19" (48 cm)	21.25" (54 cm)	14 lbs (6 kg)	350 lbs (159 kg)	Lifetime
G30213-4F	Painted Steel	Yes	Tool-free	16.5"-22.5" (42-57 cm)	16" (41 cm)	19" (48 cm)	21.25" (54 cm)	14 lbs (6 kg)	350 lbs (159 kg)	Lifetime
G30214-2	Painted Steel	Yes	No	17.5"-21.5" (46-56 cm)	16" (41 cm)	21.125" (54 cm)	23.75" (60 cm)	16 lbs (7 kg)	400 lbs (182 kg)	Lifetime
G98202	Chromed Steel	Yes	Tool-free	18.25"-23.25" (46-59 cm)	18"-20" (46-51 cm)	18.5" (47 cm)	22" (56 cm)	24 lbs (11 kg)	300 lbs (136 kg)	Lifetime
G98204	Chromed Steel	Yes Padded	Tool-free	20.25"-25.25" (51-64 cm)	18"-20" (46-51 cm)	18.5" (47 cm)	22" (56 cm)	26 lbs (12 kg)	350 lbs (159 kg)	Lifetime
G30216B	Painted Steel	Yes	No	16.5"-23.5" (42-62 cm)	18.5" (47 cm)	23.5" (60 cm)	27" (69 cm)	22 lbs (10 kg)	500 lbs (227 kg)	Lifetime

BATH SAFETY ACCESSORIES

Model	Description	G30210-4F	G30211-4	G30213-4F	G30214-2	G30400-4	G30401-4	G30402	G30403	G98202	G98204	G98013F	G98308A	G98338F	G98309/G	G98310B
G30215	Pail with lid and handle	•	•	•	•					•	•					
G98018	Pail Hanger												•			•
G98340	Suction Cup Footpiece, grey												•	•	•	•
G98240	Splash Guard	•	•	•	•					•	•	•				



Model NON770100W

CERVICAL COLLARS

Designed to offer patient comfort and support.

- Firm density foam.
- Contoured for uniform support.
- Hook and loop closure.
- Includes an extra synthetic cover that can be washed.

Model: NON770100W 3.5"H x 22"L (Retail Packaging)



Model MSC192102

INVALID RING

The perineal cushion is designed for the new mother or home convalescent with patient discomfort. The resilient molded foam ring provides economical pressure relief and includes a removable, washable cover.

Model: MSC192102



Model MED192105

INFLATABLE RUBBER CUSHION

This inflatable PVC ring is ideal for use at home, office or travel. The durable, PVC conforms to the natural contour of your body for added comfort. Easily cleaned with a damp cloth.

Model: MED192105 (Retail Packaging)



Model MSC04155

CERVICAL PILLOWS

Therapeutic pillow helps support correct sleeping posture. Durable nylex cover can be easily wiped clean.

Model: MSC04155 15" L x 5" D, covered



Model MSC019022

WHEELCHAIR CUSHION


Designed with soft, even support for maximum comfort, the Medline® wheelchair cushion is made of resilient polyurethane foam and fits standard wheelchairs.

- Durable nylex cover can easily be wiped clean.
- Size: 18" x 18" x 2"
- Color: Blue

Model: MSC019022

Medline Home Comfort Products

Guardian® home comfort and positioning products help to support and align the body for more comfortable sitting or sleeping. Proper positioning and body alignment can relieve stress that causes muscle aches and pains.



Our most popular professional-quality walking aids, bath safety and personal care products are available in attractive retail packaging designed to make it easy for consumers to shop. Ask about our merchandising programs that can help turn your dealership into a retail showroom.

Guardian® products offer the peace-of-mind of a generous warranty. Refer to actual product warranty for specific terms or limitations.

Medline Industries is committed to continuous product improvement, therefore specifications are subject to change without notice.

U.S. and foreign patents pending.



Medline United States
1-800-MEDLINE (633-5463)

www.medline.com | info@medline.com

Medline Canada
1-800-396-6996

www.medline.ca | canada@medline.com

Medline México
01-800-831-0898

www.medlinemexico.com | mexico@medline.com