




# FALL PREVENTION

MEDLINE PATIENT SAFETY ALARMS





**Fall Prevention** is a great concern facing health care today. The Centers for Disease Control and Prevention (CDC) has concluded that falls are the leading cause of injury deaths and the most common cause of nonfatal injuries and hospital admissions for trauma among patients 65 years and older.<sup>1</sup> Each year in the United States, nearly one third of older adults experience a fall.<sup>2</sup> Recent compliance regulations, as outlined by the Joint Commission's Patient Safety Goals and the Centers for Medicare & Medicaid Services (CMS) guidelines, have created increased urgency for the prevention of patient/resident falls.<sup>3</sup>

Medline can help health care facilities improve fall management by providing a number of useful products to be used in conjunction with a comprehensive fall management program. One of these products, patient safety alarms, has been highlighted by JCAHO as a tool that can be used to promote improved monitoring of persons at risk for falling.<sup>4</sup> When used properly, patient alarms can alert caregivers when a fall-risk candidate is on the move.

*For a complete listing of Medline alarms and accessories, visit [www.medline.com](http://www.medline.com) or contact your Medline sales representative.*

<sup>1</sup> Falls Among Older Adults: An Overview:  
<http://www.cdc.gov/ncipc/factsheets/adultfalls.htm>

<sup>2</sup> Costs of Falls Among Older Adults:  
<http://www.cdc.gov/ncipc/factsheets/fallcost.htm>

<sup>3</sup> 2009 National Patient Safety Goals:  
[http://www.jointcommission.org/GeneralPublic/NPSG/09\\_npsgs.htm](http://www.jointcommission.org/GeneralPublic/NPSG/09_npsgs.htm),  
Hospital-Acquired Conditions:  
[http://www.cms.hhs.gov/HospitalAcqCond/06\\_Hospital-Acquired\\_Conditions.asp](http://www.cms.hhs.gov/HospitalAcqCond/06_Hospital-Acquired_Conditions.asp)

<sup>4</sup> FAQs for the 2007 National Patient Safety Goals:  
[http://www.jointcommission.org/NR/rdonlyres/D4844675-25D7-4B5B-A47D-C549D939F9E5/0/07\\_NPSG\\_FAQs\\_9.pdf](http://www.jointcommission.org/NR/rdonlyres/D4844675-25D7-4B5B-A47D-C549D939F9E5/0/07_NPSG_FAQs_9.pdf)



## PRESSURE-SENSING SAFETY ALARMS

**Quick Alert® Ultra & Ultra II** *Our top sellers.* Quick Alert alarms are self-activated, auto-sensing systems that are always working. Each monitor is constructed from heavy-duty plastics that securely house sophisticated electronics. Quick Alert alarms have been our best-selling pressure-sensing alarms for nearly 10 years.



**NEW! Quick Alert Ultra II** MDT8175

- Adjustable delayed alarm setting limits false alarms
- Volume control button
- Output for connection to a nurse call system
- Protective silicone cover reduces breakage
- Adjustable strap for use with beds or chairs
- Power indicator light
- Low battery indicator light



**Quick Alert Ultra** MDT8100

*Perfect for wheelchair monitoring.*

- Contains all the great standard features of the Quick Alert Ultra II except volume control, delayed alarm and nurse call



**Quick Alert Ultra** pressure mats & accessories

ITEM	DESCRIPTION	PKG
MDT81180C5	Chair sensor mat, 180-day warranty	5/cs
MDT81180B5	Bed sensor mat, 180-day warranty	5/cs
MDT8190C*	Chair sensor mat, 90-day warranty	1 ea 5/cs
MDT8190B*	Bed sensor mat, 90-day warranty	1 ea 5/cs
MDT8130C*	Chair sensor mat, 30-day warranty	1 ea 5/cs
MDT8130B*	Bed sensor mat, 30-day warranty	1 ea 5/cs
MDT8150NURSE	Nurse call cable for Ultra II	1 ea/cs

*\*Items can be ordered by adding 1 for a single mat or 5 for case of 5 mats*

## PRESSURE-SENSING SAFETY ALARMS

**MDT8400/8450** *All the features you need at a great price.* The MDT8400 series patient alarms come packed with some of the most sought-after technological features typically found only on more expensive alarms, including autosensing without the need for an on/off switch. The best part is that you can choose from a number of optional accessories as needed.



### Personal Safety Monitor MDT8400

- Tamper-resistant 2-button design limits patient/resident monitoring interference
- Adjustable delayed alarm setting limits false alarms
- Volume control
- Output for connection to a nurse call system
- Low battery indicator light
- Extra-durable wheelchair attachment clip included
- Optional: Protective cover MDT84CP  
Bed bracket MDT84BRKT



### Tamper-Resistant Deluxe Personal Safety Monitor MDT8450

*All the great features of the MDT8400 plus more.*

- Built-in protective shock-absorbing rubber bumpers surround monitor
- Wireless remote control pendant for easy resetting
- AC power adapter included — perfect for bedroom monitoring
- Optional: Bed bracket MDT84BRKT



### Personal Safety Monitor pressure mats & accessories

ITEM	DESCRIPTION	PKG
MDT8445TC*	Timed chair sensor mat, 45-day warranty	1 ea 6/cs
MDT8445TB*	Timed bed sensor mat, 45-day warranty	1 ea 6/cs
MDT84180TC*	Timed chair sensor mat, 180-day warranty	1 ea 6/cs
MDT84180TB*	Timed bed sensor mat, 180-day warranty	1 ea 6/cs
MDT84365TC*	Timed chair sensor mat, 365-day warranty	1 ea 6/cs
MDT84365TB*	Timed bed sensor mat, 365-day warranty	1 ea 6/cs
MDT84NURSE	Nurse call cable for MDT8400 & MDT8450	1 ea

*\*Items can be ordered by adding 1 for a single mat or 6 for case of 6 mats*



## PULL CORD SAFETY ALARMS

Numerous features for varying situations. Medline offers a variety of pull cord alarms - choose the exact feature you want in a proven safety alarm.



### **Elite Personal Safety Alarm** MDT829900-1 & MDT829900-6

*Pacemaker safe.*

- Magnet-free plastic tab on activation cord makes for safe use near pacemakers
- Volume adjustment
- Interlocking cradle and bracket for easy installation
- 5 alarm tone options
- Battery tester button
- Hook and loop fastener option



### **Advantage Magnetic Alarm** MDT5000 & MDT5000Z

*Value priced - lowest cost per unit in its class.*

- Works from any angle
- Cradle and adjustable strap included
- Battery tester button



### **Safety Sentry Alarm** MDT5100

*Tamper-resistant.*

- No on/off switch prevents tampering
- Volume control
- Tamper-resistant safety locking clip
- Optional: Protective cover (MDT5100CP)



### **Deluxe Personal Safety Alarm** MDT5200

- Unique magnetic alarm with nurse call compatibility
- No on/off switch prevents tampering
- Volume control
- Tamper-resistant safety locking clip
- Optional: Protective cover MDT5200CP  
Bed bracket MDT84BRKT

## SPECIALTY ALARMS & ACCESSORIES



### **Velcro® Seatbelt Alarm** MDT5500 & MDT5500L

*Provides added security.*

- Seatbelt alarm system notifies of attempts of unassisted wheelchair exits
- Non-restrictive resident-release design
- Shock-absorbing protective cover
- Adjustable wheelchair attachment strap
- Seatbelt mounting ties
- Low battery indicator light
- 2 sizes: 36"-46" (91 cm - 117 cm) MDT5500  
46"-58" (117 cm - 147 cm) MDT5500L



### **Infrared Alarm** MDTIRM1

*Zero contact monitoring.*

- Motion-detecting infrared beam
- Bedside or headboard positioning
- Nurse call compatible
- 2 alarm tone options
- Optional: Metal mounting bracket MDTIRM1BRKT  
Nurse call cable MDTIRM1NURSE



### **Floor Sensor Pad** UMPFSM1

*Use alone or as a backup.*

- Alarm activates upon stepping onto mat
- Works with Tamper-Resistant Deluxe Personal Safety Monitor MDT8450



### **Under Seat Wheelchair Alarm** STM009

*Invisible wheelchair monitoring.*

- Weight-sensing system does not require disposable sensor mats
- Discreetly mounts under the wheelchair seat

### **Safe Sense™ Toilet Monitoring System** MDT8500KIT

*Restroom monitoring. (not pictured)*

- Mountable toilet sensor indicates unassisted bathroom exits
- Kit includes: alarm unit, toilet sensor strip and mesh pouch
- Replacement toilet strip sensor MDT8500TS

