

Omega Green Supreme MRI Service Vacuum by Atrix™



PRODUCT HIGHLIGHTS

WARRANTY

3 Year Limited

SERVICE VACUUM WITH IMMEDIATE CONTAINMENT

15' Nonmetallic hose and nonmetallic accessories for servicing near MRI equipment

POWERFUL

1,000 Hour Rated Ametek Advantek II Motor

QUIET

69 Decibels

SAFE

HEPA Filtration along with ESD Safety draining damaging static voltages to ground and built-in electronic line filter to minimize EMI/RFI noise emissions

APPLICATIONS

Hospital and MRI equipment cleaning



Ultrafine Vacuums Built to Last. Quality Assured.

www.atrix.com | sales@atrix.com | 1-800-222-6154 | Burnsville, MN, USA

Omega Green Supreme MRI Service Vac by Atrix™

SPECIFICATIONS

VACOSMRI

DIMENSIONS

Length: 20" (508 mm)
 Width: 6.25" (159 mm)
 Height: 9.25" (235 mm)

WEIGHT 11.15 lbs. (5.05 kg)

HOSE DIAMETER

Inside - 1.00" (25.4 mm)
 Outside - 1.25" (31.75 mm)

POWER

120, 628W

CFM

99 CFM

INCLUDED CONTENTS

- A - Omega Green Supreme MRI Service Vac by Atrix
- B - 31651-120V Detachable Power Cord
- C - 31654- NONMETALLIC ESD Safe Crevice Tool
- D - 31657- NONMETALLIC ESD Safe Flexible Gooseneck Tool
- E - 31692- 15' NONMETALLIC ESD Safe Hose
- F - OF612HE- HEPA Filter Cartridge
- G - 410017- Filter Plug
- H - AVPA016- Adjustable Shoulder Strap

FILTRATION

HEPA

99.97% Efficient at .3 Micron
 1 Gallon Dry
 56 Pleat, 3 Layer
 One Cellulose Layer
 Two Glass Layers

CERTIFICATION



GENERAL INFORMATION

PACKAGING

Dimensions: L W H
 Inches: 20.25 7.75 10.75
 mm: 514 197 273
 Weight: 16 lbs.
 7.25 kg

HS Code

8508.11.0000

UPC

8 38067 00152 0

OPTIONAL ACCESSORIES

Omega ULPA Filter Cartridge OF712UL



20' Non Metallic ESD Safe Hose 31690



HEPA CERTIFIED Omega Series Filter Cartridge - OF612HE-IC

