



Foot Pillow

Long periods of pressure over bony prominences put a patient at an increased risk for the development of pressure ulcers. The heel is one of the most prone areas to this problem. Those who are neurologically impaired, the acutely hospitalized patients, and elderly individuals are generally of great concern for this condition.



The foot pillows provide comfort and general protection from ulcers for those patients prone to this condition. By reducing pressure, blood flow is increased to the heel region, preventing ulcers from developing.

- Extra-thick polyfoam reduces pressure and cushions sensitive heels.
- Covered with soft terry cloth reduces skin shear.
- Oversized heel protectors are designed to stay in place, even on restless residents.
- Secure Velcro® closure makes securing and removal quick and easy.
- Machine washable.
- Meets CA #117.

Anti-Abduction Knee Spreader

- Comfortably holds hips in abduction while distributing pressure over thighs.
- Available with either a polyester cover or in an easy to clean vinyl.
- Strap is 46" long (117cm).
- Positions thighs 4" (10cm) apart.
- Latex free.

