


Protekt O₂ Cushion

Medicare code E2622



The Protekt O₂ cushion is Proactive Medicals premier wheelchair cushion for the prevention and treatment of pressure sores. The soft interconnected air cells allow for deep immersion providing constant low pressure across the seating surface.

- Air adjustable 4" or 2" cells offer deep immersion to provide excellent pressure redistribution for tissue protection and treatment.
- Dual valve compartments address seating challenges such as pelvic obliquities, sacral seating and enhances posture control.
- Dual compartments can be adjusted independently to increase positioning and stability for side to side control.
- Nylon breathable cover with non-skid bottom and safety straps.
- Easy to use hand pump.
- Flame resistant, latex free neoprene.
- 2 year warranty.
- Patch kit included.

78001	16"x16"x2"	78004	16"x16"x4"
78002	18"x16"x2"	78005	18"x16"x4"
78003	20"x16"x2"	78006	20"x16"x4"