


# Protekt™ 100

Pressure Redistribution Foam Mattress



 **Proactive**  
medical products

Medicare  
Code  
E0184

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

[www.proactivemedical.com](http://www.proactivemedical.com) • [sales@proactivemedical.com](mailto:sales@proactivemedical.com)

The Protekt™ 100 pressure redistribution foam mattress is specifically designed for the value oriented facilities with limited budgets. The three (3) zones are comprised of 670 individual cells that are independent from one another. This design enables the mattress to conform to the patient's body and customize support as the patient moves. Excellent pressure redistribution is achieved by varying cell depths throughout the mattress. Special attention is given to the vulnerable heel section to decrease interface pressure

## FEATURES AND BENEFITS:

- **3 Unique Pressure Zones: Channel Cut High Density Polyurethane Foam** 670 individual cells offer superior pressure redistribution. Foam core is fire-retardant and anti-microbial/anti-bacterial.
- **Top Cover:** Anti-microbial and anti-fungal; low moisture vapor transmission(MVTR); fluid resistant; conforms to NFPA702 class 1 and TB 117E.
- 360 degree zipper allows for easy removal for cleaning or replacement.
- **Side Rails: PTSS Rail System** (Patient transfer support system) Provides firm support for safe and easy patient egress and ingress from bed.
- **Anti-Contamination Flap:** This flap covers the zipper reducing the possibility of contaminating the foam core.
- **Bottom Cover:** Triple laminate vinyl; fluid resistant; anti-microbial/anti-bacterial; fire retardant.
- **Fire Retardant:** Conforms to 16 CFR Part 1632 & 1633 and NFPA 702 Class 1 .
- **Weight Capacity:** 300 Lbs.
- **Warranty:** 4 Year, Non-Prorated Warranty.

Interface pressure readings performed by independent certified testing company are available upon request.

81011	Protekt™ 100	Pressure Redistribution Foam Mattress	36"x76"x6"
81012	Protekt™ 100	Pressure Redistribution Foam Mattress	36"x80"x6"
81013	Protekt™ 100	Pressure Redistribution Foam Mattress	36"x84"x6"
81014	Protekt™ 100 Plus	Pressure Redistribution Foam Mattress w/ 3" Raised Rails	36"x76"x9"
81015	Protekt™ 100 Plus	Pressure Redistribution Foam Mattress w/ 3" Raised Rails	36"x80"x9"
81016	Protekt™ 100 Plus	Pressure Redistribution Foam Mattress w/ 3" Raised Rails	36"x84"x9"



Optional raised side rails



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com