


Protekt™ 600

Bariatric Pressure Redistribution Foam Mattress



 **Proactive**
medical products

Medicare
Code
E0184

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ 600 is constructed to provide optimal pressure redistribution and support for the high risk bariatric patient. The mattress was designed utilizing a combination of the highest density foams to create the quintessential sleeping surface for the heavier patient.

FEATURES AND BENEFITS:

- Top-layer of high resilient high density die cut foam provides optimal support and comfort.
- Heel section is made of Lura-Quilt™ foam. It is an aerated foam which allows air to continuously pass thru keeping the patient cool and dry.
- Bottom layer of foam is constructed of ultra high density high resilient foam which provides unprecedented patient support and stability.
- The top of the cover is a four way stretch nylon material with a low moisture-vapor transmission rate.
- Water resistant anti-microbial/bacterial cover is perfect for infection control.
- Anti-Contamination Flap covers the zipper reducing the possibility of fluids contaminating the foam core.
- Weight capacity: 1000 lbs.
- 16 CFR part 1632 & 1633.
- Warranty: One year non-prorated.

81061	Protekt 600-42	Bariatric Pressure Redistribution Foam Mattress	42"x76"x6"
81062	Protekt 600-42	Bariatric Pressure Redistribution Foam Mattress	42"x80"x6"
81063	Protekt 600-42	Bariatric Pressure Redistribution Foam Mattress	42"x84"x6"
81064	Protekt 600-48	Bariatric Pressure Redistribution Foam Mattress	48"x76"x6"
81065	Protekt 600-48	Bariatric Pressure Redistribution Foam Mattress	48"x80"x6"
81066	Protekt 600-48	Bariatric Pressure Redistribution Foam Mattress	48"x84"x6"
81067	Protekt 600-54	Bariatric Pressure Redistribution Foam Mattress	54"x76"x6"
81068	Protekt 600-54	Bariatric Pressure Redistribution Foam Mattress	54"x80"x6"
81069	Protekt 600-54	Bariatric Pressure Redistribution Foam Mattress	54"x84"x6"



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com