

# Protekt™ Aire 5000

Low Air Loss/Alternating Pressure Mattress System



 **Proactive**  
medical products

Medicare  
Code  
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 5000 is the 1st defense for facilities to treat pressure ulcers stages I-IV. The high volume user friendly pump combines low air loss and alternating pressure in one unit. The air flow and pressure range can be customized according to each individual patients weight. The “static” mode provides firm support for easy patient transfer. The mattress is comprised of (20) 8” air bladders that allow for deep immersion and excellent pressure redistribution. The top cover is quilted for added patient comfort.

## FEATURES AND BENEFITS

### MATTRESS & COVER

- 18 laser air holes to keep patient cool and dry.
- 5” mattress with 3” foam base provides a safety support in case of power failure.
- Removable, water resistant, non-shear, anti-microbial nylon cover is quilted for extra comfort.
- Non-skid cover prevents the mattress from slipping.
- CPR valve provides rapid deflation of mattress.
- Dimensions: 35”x80”x8”.
- Weight: 9 lbs.

### PROTEKT™ AIRE 5000 PUMP

- High flow pump with 8 liters of air per minute.
- Adjustable comfort settings allows pressure to be customized according to each individual patients therapeutic requirements.
- Low pressure alarm indicates when the pumps airflow may be compromised.
- Static function for safe and easy patient transfer.
- Low pressure alarm mute button.
- UL/CSA approved power cord.
- Weight Capacity: 350 lbs.
- Dimension: 11”x8”x4”.
- Weight: 6.2 lbs.
- 18 month warranty.



Transport Handles




Easy Valve Exchange



Low Air Loss



Transport Connector



Non-skid Bottom



CPR Release



Quilted Cover

80050	Protekt™ Aire 5000	8” Alternating/Low Air Loss Mattress System
80051	Protekt™ Aire 5000 Pump	Pump only for Protekt™ Aire 5000
80052	Protekt™ Aire 5000 Mattress	Mattress only for Protekt™ Aire 5000
80053	Protekt™ Aire 5000 Cover	Cover only for Protekt™ Aire 5000
80054	Protekt™ Aire 5000 Bladder	Individual Bladder for Protekt™ Aire 5000



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com