

Protekt™ Aire 8000BA

Low Air Loss/Alternating Pressure Bariatric Mattress System



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 8000BA Bariatric Mattress Series is the first choice by clinicians to treat and prevent bariatric patients from pressure ulcers.

FEATURES

- Alternating pressure cycles of 10, 15, 20 & 25 allows clinicians to customize pressure redistribution and skin surface interface.
- Low air loss function controls and manages moisture to prevent skin maceration.
- Cell-on-cell mattress design prevents patient from bottoming out in case of power failure.
- 10" air cells are constructed of Nylon/PU for added comfort and support.
- Audible low pressure alarms alerts clinicians when preset pressure levels decrease.
- Static mode allows for easy transfer in and out of bed.
- Pillow top feature prevents head and neck from unnecessary movement during alternating mode.
- CPR pull tab quickly releases air from mattress.
- 2 year warranty.
- Weight: 9 lbs.

SPECIFICATIONS

- Air Volume: 15 liters per minute
- Cycle Time: 10, 15, 20 and 25 minutes, selectable.
- Pump Weight: 8.8 lbs.
- Pump Dimensions: 13.4"x4.9"x9.7".
- Mattress Dimensions: 42"x80"x10", 48"x80"x10".
- Cell Depth: 10".
- UL/CSA approved power cord.
- Power: 120 VAC.
- Input Frequency: 60Hz.
- Pressure Range: 20-70 mm Hg.
- Weight Capacity: 650 lbs.

| | | |
|-------|----------------------------------|---|
| 80080 | Protekt™ Aire 8000BA-42 | Low Air Loss/Alternating Pressure Mattress System-42" |
| 80081 | Protekt™ Aire 8000BA-42 Pump | Pump only for Protekt™ Aire 8000BA-42 |
| 80082 | Protekt™ Aire 8000BA-42 Mattress | Mattress only for Protekt™ Aire 8000BA-42 |
| 80083 | Protekt™ Aire 8000BA-42 Cover | Cover only for Protekt™ Aire 8000BA-42 |
| 80084 | Protekt™ Aire 8000BA-42 Bladder | Individual Bladder for Protekt™ Aire 8000BA-42 |
| 80085 | Protekt™ Aire 8000BA-48 | Low Air Loss/Alternating Pressure Mattress System-48" |
| 80086 | Protekt™ Aire 8000BA-48 Pump | Pump only for Protekt™ Aire 8000BA-48 |
| 80087 | Protekt™ Aire 8000BA-48 Mattress | Mattress only for Protekt™ Aire 8000BA-48 |
| 80088 | Protekt™ Aire 8000BA-48 Cover | Cover only for Protekt™ Aire 8000BA-48 |
| 80089 | Protekt™ Aire 8000BA-48 Bladder | Individual Bladder for Protekt™ Aire 8000BA-48 |



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com