

**PURSUIT**

## OUTSTANDING OUTDOOR PERFORMANCE WITH LUXURIOUS FEATURES

### PURSUIT

The **Pursuit** features a powerful drivetrain and large 13" wheels to handle rugged outdoor terrain. Luxury touches like a wraparound delta tiller, low profile, non-marking tires and a full lighting package with directional signals, combine seamlessly with the Pursuit's sleek, sporty look.



**13" low-profile  
wheels**



**Full suspension**



**Full lighting package  
with directional  
signals**



**Delta tiller**

**Pride**<sup>®</sup>  
Mobility Products Corp.  
— LIVE YOUR BEST —

# Pursuit

Model number:	SC713
Weight capacity:	400 lbs.
Maximum speed <sup>1</sup> :	Up to 9.32 mph - high speed/4.66 mph - low speed
Ground clearance <sup>2</sup> :	3.56"
Turning radius:	67.25"
Overall length <sup>2</sup> :	53.25"
Overall width <sup>2</sup> :	25.25"
Front tires:	4" x 13" solid
Rear tires:	4" x 13" solid
Seat to deck:	17"-18.25"
Suspension:	Front & Rear
Range per charge <sup>1</sup> :	Up to 14 miles
Seat-to-ground:	24.75" - 26.25"
Total weight without batteries:	218.5 lbs.
Heaviest piece when disassembled:	95 lbs. (front section)
Standard seating:	Type: CRS reclining high-back seat with headrest Weight: 45 lbs. Material: Black Vinyl Dimensions: width - 18" (overall width) depth - 18" (effective depth)
Drive system:	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual braking system:	Regenerative and electromechanical
Battery requirements <sup>3</sup> :	Two 12-volt, deep cycle Size: NF-22
Battery weight:	NF-22: 38 lbs. each
Battery charger:	5-amp, off-board
Colors:	



Candy Apple Red



Viper Blue



Silver

<sup>1</sup>Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tyre condition

<sup>2</sup>Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%)

<sup>3</sup>AGM or Gel-Cell type required

<sup>4</sup>Battery weight may vary based on manufacturer

<sup>5</sup>Includes standard seat

NOTE: All specifications subject to change without notice

## Features

- Infinitely adjustable tiller angle w/conveniently located handle
- Feather-touch disassembly
- Powerful drive train, features 120 amp microprocessor based controller
- 24V, 4-pole DC motor for increased power
- Deluxe reclining high-back seat with headrest and sliders
- Full lighting package including directional signals and angle-adjustable, lower headlight
- Delta tiller with wraparound handles
- Large 13" low profile solid tires for excellent outdoor performance
- Front and rear suspension
- Easily accessible tie down points (for transport of unoccupied scooter)
- Front basket standard
- Rear view mirror standard

## Options

- Rear basket
- Safety flag
- Durable weather cover



The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %.

800-800-8586 (US) • 888-570-1113 (Canada) • [www.pridemobility.com](http://www.pridemobility.com)

