

Restore Foam Dressings with Silicone

Performance. Comfort. Value.




- Fluid handling you need
- Comfort your patient demands
- Value you can count on




Restore Foam Dressings with Silicone

Exceptional Fluid Management,
Secure Yet Gentle Hold.




Large perforations designed to improve fluid absorption.

- Perforations are 2.5X larger than a leading foam dressing, which speeds absorption and helps to minimize the risk of maceration.*



Restore foam dressing with silicone



Mepilex® Border foam dressing


*Compared to Mepilex Border dressing.¹
Images shown at 30X magnification

Performance.

Large perforations and absorbent foam make the difference.


- Large perforations allow for rapid absorption of exudate
- Highly absorptive foam wicks exudate away from the wound bed

Speed of absorption¹



Error bars represent one standard deviation. (N=10)
*Brands with silicone interface over pad

Total fluid handling¹



Error bars represent one standard deviation. (N=10)

Comfort.

During application. During wear. During removal.

- **Silicone interface**
 - Provides secure, yet gentle hold
 - Provides gentle tack for easier application
 - Minimizes adhesion to the wound bed
 - Reduces risk of trauma and pain at dressing change
- **Smooth, film backing to ease sliding across bedding and clothing**
- **Soft, conformable dressing for difficult to dress areas**



Value.

Simply performs.

- Fluid handling you need
- Comfort your patient demands
- Priced lower than leading competitive brands**
- Quality backed by Hollister

Restore
Foam Dressings with Silicone

** Pricing compared to similar sizes of Mepilex, Mepilex Border, Alleyn Gentle, and Alleyn Gentle Border dressings. Based on 03.15 analysis of FSS and GHX 2013 data.

Indications for Use

Restore foam dressings with silicone are indicated for moderate to heavily exuding acute and chronic wounds such as pressure ulcers, diabetic ulcers, venous stasis ulcers, arterial ulcers, leg ulcers, malignant wounds, surgical wounds, donor sites, partial-thickness burns, skin tears, skin abrasions, and postoperative wounds.

Ordering Information:

	Stock No	Product Size	Quantity/Box	HCPCS
	Restore foam dressing with silicone – border			
	520003	3.5" x 3.5" (9.0cm x 9.0cm) <i>pad size 2.5" x 2.5" (6.4cm x 6.4cm)</i>	10	A6212
	520005	5" x 6.2" (12.7cm x 15.7cm) <i>pad size 3.4" x 4.6" (8.6cm x 11.7cm)</i>	10	A6212
	520000	7" x 7" (18.0cm x 18.0cm) <i>pad size 5.5" x 5.5" (14.0cm x 14.0cm)</i>	10	A6213
	Restore foam sacral dressing with silicone - border			
	520001	7.3" x 7.4" (18.6cm x 18.8cm) <i>pad size 6.1" x 5.7" (15.4cm x 14.6cm)</i>	10	A6213
	Restore foam dressing with silicone - non-border			
	520023	2.5" x 2.5" (6.4cm x 6.4cm)	10	A6209
	520025	4" x 5" (10.2cm x 12.7cm)	10	A6210
	520027	6.5" x 8" (16.5cm x 20.3cm)	10	A6211
	Restore foam heel dressing with silicone - non-border			
	520020	8" x 5" (20.3cm x 12.7cm)	10	A6210

1. Internal Data on File.



Hollister Incorporated

2000 Hollister Drive
Libertyville, Illinois 60048 USA
1.888.740.8999

www.hollister.com

The Hollister logo and Restore are trademarks of Hollister Incorporated. All other trademarks are the property of their respective owners.

©2015 Hollister Incorporated.
922817-0415

Restore
Foam Dressings with Silicone